

Aug 06, 2018

Endo H Denatured Protocol for Deglycosylating Glycoproteins

DOI


dx.doi.org/10.17504/protocols.io.sedeba6

Mike Gibson¹

¹QA-Bio Inc



Mike Gibson



OPEN ACCESS



DOI: dx.doi.org/10.17504/protocols.io.sedeba6

External link: <http://www.qa-bio.com/docs/E-EH02.QA-Bio.specsheet.pdf>

Protocol Citation: Mike Gibson 2018. Endo H Denatured Protocol for Deglycosylating Glycoproteins. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.sedeba6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 06, 2018

Last Modified: August 06, 2018

Protocol Integer ID: 14501

Keywords: Endo H, Endoglycosidase H



Abstract

Endo H cleaves Asparagine-linked hybrid or high mannose oligosaccharides, but not complex oligosaccharides. It cleaves between the two N-acetylglucosamine residues in the diacetylchitobiose core of the oligosaccharide, generating a truncated sugar molecule with one N-acetylglucosamine residue remaining on the asparagine. In contrast, PNGase F removes the oligosaccharide intact. Detergent and heat denaturation may increase the rate of cleavage for some glycoproteins.

Materials

MATERIALS

 Endo H **QA-Bio Inc Catalog #E-EH02**



- 1 Add up to 200 μg of glycoprotein to an Eppendorf tube. Adjust to 37.5 μl final volume with de-ionized water.
- 2 Add 10 μl 5x Reaction Buffer 5.5 and 2.5 μl of Denaturation Solution. Heat at 100°C for 5 minutes.(NOTE: It is not necessary to add Triton X-100. SDS will not inactivate Endo H.)
- 3 Add 2.0 μl of Endo H to the reaction. Incubate 3 hours at 37°C.
- 4 Monitor cleavage by SDS-PAGE.