



Sep 04, 2024

Endo H/ PNGase F Assay for JESS Automated Western Blot

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2qzyjl1y/v1>

Charis Ma¹, Yu Chen¹, sidranse¹

¹National Human Genome Research Institute



Charis Ma

National Institutes of Health

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2qzyjl1y/v1>

Protocol Citation: Charis Ma, Yu Chen, sidranse 2024. Endo H/ PNGase F Assay for JESS Automated Western Blot .
protocols.io <https://dx.doi.org/10.17504/protocols.io.4r3l2qzyjl1y/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 03, 2024



Last Modified: September 04, 2024

Protocol Integer ID: 106891

Keywords: jess automated western blot system, jess automated western blot, western blot, assay, endo, protocol detail

Funders Acknowledgements:

ASAP

Grant ID: ASAP-024297

Abstract

This protocol details an optimized Endo H and PNGase F digestion assay that produces samples that can be run with the JESS Automated Western Blot System.

Materials

10X Glycoprotein Denaturing Buffer (NEB, Catalog# B0701S)

10X GlycoBuffer 3 (NEB, Catalog# B0701S)

10X GlycoBuffer 2 (NEB, Catalog# B0701S)

10% NP-40 (NEB, Catalog# B0701S)

EndoH (NEB, Catalog# P0703S)

PNGase F (NEB, Catalog# P0704S)

10K Amicon Ultra Centrifugal Filter (Millipore, UFC5010BK)



Glycoprotein Denaturation

- 1 In an eppendorf tube, add  40 μg protein,  2 μL 10X Glycoprotein Denaturing Buffer, and water for a total reaction volume of  20 μL . Mix well.
- 2 Denature glycoprotein by heating reaction at  100 $^{\circ}\text{C}$ for  00:10:00 . 10m
- 3 Split reaction volume into 2 eppendorf tubes with 10 μL in each tube. One tube will be used for the EndoH reaction, while the other will be used for the PNGaseF reaction.

Endo H

- 4 Add the following to one reaction tube from step 3:  2 μL 10X GlycoBuffer 3,  2 μL Endo H,  6 μL water. Total volume will be  20 μL .
- 5 Mix gently.
- 6 Incubate reaction at  37 $^{\circ}\text{C}$ for  01:00:00 1h

PNGase F

- 7 Add the following to the second reaction tube from step 3:  2 μL 10X GlycoBuffer 2,  2 μL NP-40,  6 μL water,  1 μL PNGaseF. Total volume will be  21 μL .
- 8 Mix gently.
- 9 Incubate reaction at  37 $^{\circ}\text{C}$ for  01:00:00 . 1h

Desalting and Washing for JESS Automated Western Blot System (Bio-Techne)



preparation

- 10 Add  480 μL and  479 μL of 0.1X Sample Buffer (Bio-Techne) to the Endo H and PNGaseF reaction tubes respectively to make the final volume  500 μL .
- 11 Add all  500 μL of each reaction to a 10K Amicon Ultra Centrifugal Filter
- 12 Spin at  14000 x g for  00:20:00 . You will recover a final volume of  20 μL , and the protein concentration recovered will be ~ [M] 1 mg/mL .

20m

JESS Sample Preparation

5m

- 13 Add 4 parts of the washed Endo H/ PNGaseF sample from step 12 with 1 part of 5X Fluorescent Master Mix (Bio-Techne).
- 14 Vortex to mix
- 15 Incubate sample at  95 $^{\circ}\text{C}$ for  00:05:00 .
- 16 Load into JESS microplate (Bio-techne) to be detected by the Automated Western Blot System.

5m

Protocol references

<https://www.neb.com/en-us/protocols/2012/10/18/endo-hf-protocol>

<https://www.neb.com/en-us/protocols/2014/07/31/pngase-f-protocol>