

Jun 15, 2023

Endosome isolation

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn3bp6l5d/v1>

Chuyu Chen¹

¹Northwestern University, Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815



Chuyu Chen

Northwestern University, Aligning Science Across Parkinson's...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn3bp6l5d/v1>

Protocol Citation: Chuyu Chen 2023. Endosome isolation . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.14egn3bp6l5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 15, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 83503

Keywords: ASAPCRN, late endosome, subcellular fractionation, centrifugation step, series of centrifugation step

Funders Acknowledgements:

Aligning Science Across Parkinson's [ASAP-020600] through the Michael J. Fox Foundation for Parkinson's Research (MJFF)
Grant ID: ASAP-020600

Abstract

Subcellular fractionation to isolate early and late endosomes (EEs and LEs) by performing a series of centrifugation steps.



Continuous sucrose gradient tube preparation

- 1 prepare 0.5M and 2M sucrose buffer
- 2 prepare buffer S2 by mixing 0.5M and 2M sucrose at 1:2 and buffer S3 at 2:1 ratio
- 3 place Beckman ultracentrifuge tube (No.344059) on dry ice
- 4 add 2.5ml 2M sucrose in the tube and wait for  00:10:00 or until it frozen. 10m
- 5 add 2.5 ml buffer S2 on top of frozen 2M sucrose and wait for  00:10:00 or until it frozen 10m
- 6 add 2.5 ml buffer S3 on top of frozen buffer S2 and wait for  00:10:00 or until it frozen 10m
- 7 add 2.5 ml 0.45M sucrose on top of frozen buffer S3 and wait for  00:10:00 or until it frozen 10m
- 8 keep the frozen tubes at -80 until experiment day
- 9 place the frozen tubes at 4C for  05:00:00 or until defrosted before use 5h

Fractionation

- 10 Dissect striatum from mice, fresh freeze the tissue with dry ice
- 11 Prepare homogenization buffer (HB): 150mM NaCl, 10mM HEPES, 1mM EGTA, 0.1M MgCl



12 Homogenize striatum tissue with 1ml HB containing 1x Halt protease and phosphatase inhibitor (EDTA-free) by Dounce homogenization on ice (30 strokes)

13 Centrifugation ⌚ 00:10:00 , 3,000g, 4C

10m

14 load postnuclear supernatants onto the top of continuous sucrose gradient

15 Centrifugation ⌚ 18:00:00 , 200,000g, 4C (Accel and Decel → slow)

18h

16 Collect fractions from top to bottom of the gradient

TCA precipitation

50m 3s

17 add equal volume of ice-cold 55% TCA and incubate on ice for ⌚ 00:30:00

30m

18 centrifugation ⌚ 00:10:00 , 20,000g, 4C

10m

19 wash pellet with 1ml ice-cold acetone centrifugation ⌚ 00:10:00 , 20,000g, 4C

10m

20 ➡ go to step #19 for 2 times. Total 3 washes

21 resuspend pellet with 70ul of 1% SDS buffer (1% SDS, 10mMTris, 2mMEDTA) containing 1x Halt protease and phosphatase inhibitor (EDTA-free), sonication for ⌚ 00:00:03 and followed by incubation on orbital shaker for 1hr at room temperature.

3s

22 Store samples at -80 for western blot assay