

Oct 17, 2019

## Energy solution preparation for OnePot PURE cell-free system

DOI

[dx.doi.org/10.17504/protocols.io.8ayhsfw](https://dx.doi.org/10.17504/protocols.io.8ayhsfw)



Konstantinos Ragios<sup>1</sup>

<sup>1</sup>EPFL - EPF Lausanne

iGEM EPFL



Konstantinos Ragios

EPFL - EPF Lausanne

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8ayhsfw](https://dx.doi.org/10.17504/protocols.io.8ayhsfw)

External link: <https://pubs.acs.org/doi/10.1021/acssynbio.8b00427>

**Protocol Citation:** Konstantinos Ragios 2019. Energy solution preparation for OnePot PURE cell-free system. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8ayhsfw>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 16, 2019

**Last Modified:** October 17, 2019

**Protocol Integer ID:** 28728

**Keywords:** energy, PURE, cell-free, iGEM, expression, tRNA, NTPs

## Abstract

In this protocol we explain the procedure to create the Energy solution used for protein expression in OnePot PURE cell-free system.

## Materials

Material/Consumables:

- Amino Acids (Sigma-Aldrich: LAA21-1KT)
- Magnesium acetate (Sigma-Aldrich: M0631)
- Potassium glutamate (Sigma-Aldrich: 49601)
- DTT (SantaCruz Biotech: sc-29089B)
- ATP (ThermoFisher: R0481)
- GTP (ThermoFisher: R0481)
- CTP (ThermoFisher: R0481)
- UTP (ThermoFisher: R0481)
- tRNA (Roche:10109541001)
- Creatine phosphate (Sigma-Aldrich: 27920)
- Folinic acid (Sigma-Aldrich: PHR1541)
- Spermidine (Sigma-Aldrich: S2626)
- HEPES (Sigma-Aldrich: H0887-100ML)
- biDistilled water

Equipment:

- Vortex
- Nanodrop Spectrophotometer

## Before start

Keep all the components in ice while preparing the Energy solution.



- 1 For a 2.5x Energy Solution add the materials needed to a tube. The final concentration of the components is presented in Table 1.

**Note**

Before adding each component make sure it is totally melted and then vortex for a few seconds

2

Compound	Concentration	Units
Amino Acids	0.75	<i>mM</i>
Magnesium acetate	29.5	<i>mM</i>
Potassium glutamate	250	<i>mM</i>
DTT	2.5	<i>mM</i>
ATP	5	<i>mM</i>
GTP	5	<i>mM</i>
CTP	2.5	<i>mM</i>
UTP	2.5	<i>mM</i>
tRNA	130	UA260/mL
Creatine phosphate	50	<i>mM</i>
Folinic acid	0.05	<i>mM</i>
Spermidine	5	<i>mM</i>
HEPES	125	<i>mM</i>



#### Note

Because tRNA degrades over time before adding it to the solution you need to check the  $A_{260}$  units.

- 3 To complete the solution add biDistilled water until you reach the desirable volume, and then vortex for a few seconds
- 4 The Energy solution is stored at  $-80^{\circ}\text{C}$ .