



Sep 01, 2025

Enrichment of phosphopeptides from lysate using TiO_2

DOI

<https://dx.doi.org/10.17504/protocols.io.rm7vzqdp2vx1/v1>

Jasmin Jansen^{1,2,3}

¹University of Konstanz, Department of Biology; ²Aligning Science Across Parkinson's;

³Konstanz Research School of Chemical Biology (KoRS-CB)



Jasmin Jansen

University of Konstanz

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzqdp2vx1/v1>

Protocol Citation: Jasmin Jansen 2025. Enrichment of phosphopeptides from lysate using TiO_2 . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzqdp2vx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it has worked in previous experiments.



Created: May 22, 2025

Last Modified: September 01, 2025

Protocol Integer ID: 218798

Keywords: ASAPCRN, phosphoenrichment, TiO₂-MOAC, phosphopeptide, sample preparation for phosphopeptide enrichment, phosphopeptide enrichment, enrichment of phosphopeptide, enriched phosphopeptide, rich repeat kinase, phosphopeptide, phosphorylation, phosphorylation activity, cell lysate, wide interactomes of leucine, lysis of these cell, mass spectrometry, based mass spectrometry, proteome, lysi

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP))

Grant ID: ASAP-000519

Abstract

This protocol describes the sample preparation for phosphopeptide enrichment from lysed HEK293T cells that were transfected with TurboID, TurboID-LRRK1 or TurboID-LRRK2 and treated with MLI-2, IN04 or DMSO. It is meant to accompany the methods section and manuscript of "Probing the proteome-wide interactomes of Leucine-rich Repeat Kinases 1 and 2 and alterations in their phosphorylation activity".

For details on culturing and treatment of cells as well as lysis of these cells, please refer to Protocol "TurboID-phospho protocol" (doi: 10.17504/protocols.io.ewov1drpyvr2/v1) where these steps are explained in detail. This protocol is only meant to describe the sample processing from a cell lysate to enriched phosphopeptides ready to be measured by DIA-based mass spectrometry.

This protocol starts with the cell lysate obtained in the above mentioned protocol at step 14 and ends with enriched phosphopeptide that can be analysed by mass spectrometry according to steps 45ff of the same protocol.

All steps are described per sample.

Materials

-  cOmplete™, EDTA-free protease inhibitor cocktail **Roche Catalog #11873580001**
-  Tris(2-carboxyethyl)phosphin -hydrochlorid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C4706**
-  Iodoacetamide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I1149-5G**
-  Trypsin/Lys-C Mix, Mass Spec Grade, 5 × 20ug **Promega Catalog #V5073**
-  Sequencing Grade Modified Trypsin **Promega Catalog #V5113**
-  HiPPR®; Detergent Removal Spin Column Kit **Thermo Fisher Catalog #88305**
-  Sep-Pak tC18 1 cc Vac Cartridge 50 mg Sorbent per Cartridge **Waters Catalog #WAT054960**
-  Titansphere TiO₂, 5 µm, 500 mg **GL Science Catalog #020-75000**
-  Glycolic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #124737**
-  3M™ Empore™ C18 47 mm Extraction Disc Model 2215 20 pack 3 packs per case **3M corporation Catalog #2215**
-  Lactic acid solution, 85 % **Merck MilliporeSigma (Sigma-Aldrich) Catalog #252476**

Protocol materials

-  Tris(2-carboxyethyl)phosphin -hydrochlorid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C4706**
-  Iodoacetamide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I1149-5G**
-  Sequencing Grade Modified Trypsin **Promega Catalog #V5113**
-  Sep-Pak tC18 1 cc Vac Cartridge 50 mg Sorbent per Cartridge **Waters Catalog #WAT054960**
-  Lactic acid solution, 85 % **Merck MilliporeSigma (Sigma-Aldrich) Catalog #252476**
-  Glycolic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #124737**
-  cOmplete™, EDTA-free protease inhibitor cocktail **Roche Catalog #11873580001**
-  Trypsin/Lys-C Mix, Mass Spec Grade, 5 × 20ug **Promega Catalog #V5073**
-  3M™ Empore™ C18 47 mm Extraction Disc Model 2215 20 pack 3 packs per case **3M corporation Catalog #2215**
-  HiPPR®; Detergent Removal Spin Column Kit **Thermo Fisher Catalog #88305**
-  Titansphere TiO₂, 5 µm, 500 mg **GL Science Catalog #020-75000**



Safety warnings

⚠ After adding acetone to your proteins, always work on ice until you have dried the pellet. If your sample warms up during the precipitation process, you might lose proteins.

After performing the phosphopeptide enrichment, make sure that the exposition to basic pH is as short as possible because this can cause a loss of phosphorylations

Perform the enrichment shortly before measurement to keep native phosphorylation intact.

Before start

You will need 6 volumes of ice-cold acetone for 1 volume of your sample to perform the precipitation. Please calculate the total sample volume before you start and choose the right tube size for the precipitation.

Acetone-Protein precipitation

- 1 Thaw  200 μL cell lysate in NP-40 buffer ( 25 millimolar (mM) Tris-HCl  7.4 ,  150 millimolar (mM) NaCl,  1 % (v/v) NP-40,  5 millimolar (mM) MgCl_2 ,  1 millimolar (mM) Na_3VO_4 ,  5 millimolar (mM) NaF,  5 millimolar (mM) β -glycero phosphate, 1x  cComplete™, EDTA-free protease inhibitor cocktail **Roche Catalog #11873580001** ,  1 millimolar (mM) DTT) or other lysis buffers at  0 °C .
- 2 Add 1 volumes (here  200 μL) of ice-cold acetone and mix sample.
 - 2.1 Add 2 volumes (here  400 μL) of ice-cold acetone and mix sample.
 - 2.2 Add 3 volumes (here  600 μL) of ice-cold acetone and mix sample.
- 3 Incubate at  -80 °C for  00:15:00 , then at  -20 °C for  01:30:00 .
- 4 Centrifuge precipitated samples for  00:10:00 at  16000 x g, 4°C .
 - 4.1 Afterwards, you should be able to see a whitish pellet of precipitated proteins (depending on the protein amount).
 - 4.2
- 5 Carefully remove supernatant without disturbing the pellet.
 - 5.1 Optionally, wash your pellet with ice-cold acetone one or two more times. This depends on the downstream processing. Washing of the pellet can for example reduce the amount of detergents in the sample.
- 6 Dry the pellet at  Room temperature with an open lid.



Protein reduction, alkylation and digestion

7 Resuspend protein pellet in  150 μ L  8 Mass Percent urea.

7.1 Incubate for  00:15:00 at  37 °C and  800 rpm to resuspend proteins.

8 Add TCEP

 Tris(2-carboxyethyl)phosphin -hydrochlorid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #C4706**

for  3 millimolar (mM) final concentration to reduce cysteines.

8.1 Incubate for  00:30:00 at  37 °C and  650 rpm .

9 Add IAA (

 Iodoacetamide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #I1149-5G**) to

 6 millimolar (mM) final concentration to alkylate cysteines.

9.1 Incubate for  00:30:00 at  Room temperature and  650 rpm **in the dark.**

10 Dilute urea to  4 Mass Percent using  50 millimolar (mM) ammonium bicarbonate.

11 Add  1 μ g

 Trypsin/Lys-C Mix, Mass Spec Grade, 5 $\times$ 20ug **Promega Catalog #V5073** to sample.

11.1 Incubate for  03:30:00 at  37 °C and  650 rpm .

12 Dilute urea concentration to  1 Mass Percent using  50 millimolar (mM) ammonium bicarbonate.

13 Add  2 μ g  Sequencing Grade Modified Trypsin **Promega Catalog #V5113** .





13.1 Incubate for 18:00:00 at 37 °C and 650 rpm .



14 Remove residual detergents with the

HiPPR[®]; Detergent Removal Spin Column Kit **Thermo Fisher Catalog #88305**

(200 µL per sample) according to the manufacturer's instructions.

15 Acidify flow-through by adding formic acid to a final concentration of 2 % (v/v) .

15.1 After mixing, check if 2 or lower. This is important since acidic pH stops trypsin digestion and is crucial for a successful C18 desalting.



16 Add 1 % (v/v) acetonitrile (ACN) and desalt your sample using

Sep-Pak tC18 1 cc Vac Cartridge 50 mg Sorbent per Cartridge **Waters Catalog #WAT054960**

and a vacuum manifold.

16.1 Wet cartridges with 1 mL pure ACN.

16.2 Equilibrate cartridges with 2 mL 1 % (v/v) ACN, 0.1 % (v/v) formic acid.

16.3 Add sample and apply it slowly to the C18 resin. Do not open valve fully to load sample slowly.

16.4 Wash peptides with 2 mL 1 % (v/v) ACN, 0.1 % (v/v) formic acid.

16.5 Elute peptides into low-binding 1.5 mL tube with 1 mL 50 % (v/v) ACN, 0.1 % (v/v) formic acid.

17 Dry desalted peptides in a vacuum concentrator.



Phosphopeptide enrichment



- 18 Weigh  2 mg beads per sample (if ~  200 μg peptide input)
 Titansphere TiO, 5 μm , 500 mg **GL Science Catalog #020-75000** and resuspend
in  10 μL per 1 mg beads **Loading buffer 1** ( 0.1 Molarity (M)
 Glycolic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #124737** ,
 70 % (v/v) ACN,  5 % (v/v) trifluoroacetic acid (TFA)).
- 18.1 To minimize unspecific binding, incubate TiO_2 beads in **Loading buffer 1** for
 00:20:00 at  Room temperature and  1200 rpm prior to peptide
incubation.
- 19 Resuspend approx.  200 μg of dried peptides in  200 μL **Loading buffer 1**.
- 19.1 Incubate for  00:10:00 at  37 $^{\circ}\text{C}$ and  650 rpm for proper resuspension.
- 20 Add equilibrated beads to peptides at a **peptide:bead ratio of 1:10**.
- 21 Incubate for  00:20:00 at  Room temperature on rolling incubator.
- 22 Prepare self-made, single-layered C8 Stage tips. See

Citation

Rappsilber J, Mann M, Ishihama Y (2007)
. Protocol for micro-purification, enrichment, pre-fractionation and storage of peptides
for proteomics using StageTips.

<https://doi.org/>

LINK

for details.

Use



3M™ Empore™ C18 47 mm Extraction Disc Model 2215 20 pack 3 packs per
case **3M corporation Catalog #2215**

to make the tips.



- 23 Equilibrate C8 Stage tips with **Loading buffer 1**.
- 24 Settle incubated TiO₂ beads by centrifugation  10000 x g for  00:02:00 .
- 24.1 Transfer  150 µL of supernatant to fresh low-binding tube and keep for second enrichment step, see  [go to step #31](#) .
- 25 Use residual supernatant to transfer beads to equilibrated C8 Stage tips.
- 25.1 Centrifuge for  00:02:00 at  500 x g and add flow-through to supernatant fraction from step  .
- 26 Wash beads that were retained by C8 Stage tips with  50 µL **Wash buffer 1** ( 80 % (v/v) ACN,  1 % (v/v) TFA) and centrifuge at  500 x g for  00:02:00 .
- 27 Wash beads with  50 µL **Wash buffer 2** ( 10 % (v/v) ACN,  0.2 % (v/v) TFA) and centrifuge at  500 x g for  00:02:00 . Transfer tips to new vial.
- 28 Prepare low-binding tubes with  60 µL  10 % (v/v) formic acid for phosphopeptide elutions to directly acidify them and minimize loss of phosphorylation by ammonium hydroxide.
- 29 Elute bound phosphopeptides with  30 µL **Elution buffer 1** ( 1 % (v/v) NH₄OH) by centrifugation at  500 x g for  00:02:00 and directly transfer elution to prepared tubes ().
- 30 Elute residual bound phosphopeptides with  30 µL **Elution buffer 2** ( 5 % (v/v) NH₄OH,  25 % (v/v) ACN) by centrifugation at  500 x g for  00:02:00 and directly transfer elution to prepared tubes ().
- 31 Repeat phosphopeptide protocol with retained supernatant from  and repeat protocol with **Loading buffer 2** ( 20 % (v/v)).



Lactic acid solution, 85 % **Merck MilliporeSigma (Sigma-Aldrich) Catalog #252476**

, [M] 70 % (v/v) ACN, [M] 5 % (v/v) TFA) instead of **Loading buffer 1**.

This second enrichment step maximizes yield of phosphopeptide by repeating the enrichment and changing the specificity by using lactic acid instead of glycolic acid.

32 Combine phosphopeptide elutions and dry them by vacuum evaporation.

33 Samples are now ready for mass spectrometric analysis. Please refer to protocol "TurboID-phospho protocol" ([dx.doi.org/10.17504/protocols.io.ewov1drpyvr2/v1](https://doi.org/10.17504/protocols.io.ewov1drpyvr2/v1)) for details.

Acknowledgements

This research was funded by Aligning Science

Across Parkinson's (grant ASAP-000519) through the Michael J. Fox

Foundation for Parkinson's Research (MJFF). This work was also supported by

funding of the German Research Foundation (DFG, 496470458 and 516836828).