

Feb 11, 2025

Environmental DNA (e-DNA) extraction from soil

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp9295vzp/v1>

Nuttapon Pombubpa^{1,2}, Ariya Chindamporn^{3,2}, Nattapol Kraisitudomsook^{4,2}, Neil Andrew Robert Gow^{5,2}, Syahriar Nur Maulana Malik Ibrahim^{1,2}, Fikran Aranda Fahrudin^{6,2}

¹Department of Microbiology, Faculty of Science, Chulalongkorn University, Bangkok 10330, Thailand;

²Southeast Asia Antifungal Resistance Monitoring Initiative (SEA-ARMi), Department of Microbiology, Faculty of Science, Chulalongkorn University, Bangkok 10330, Thailand;

³Department of Microbiology, Faculty of Medicine, Chulalongkorn University, Bangkok 10330, Thailand;

⁴Department of Biology, Faculty of Science and Technology, Muban Chombueng Rajabhat University, Ratchaburi 70150, Thailand;

⁵Medical Research Council Centre for Medical Mycology, University of Exeter, Exeter, United Kingdom;

⁶Program of Biotechnology, Faculty of Science, Chulalongkorn University, Bangkok 10330, Thailand



Syahriar Nur Maulana Malik Ibrahim

Chulalongkorn University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp9295vzp/v1>

Protocol Citation: Nuttapon Pombubpa, Ariya Chindamporn, Nattapol Kraisitudomsook, Neil Andrew Robert Gow, Syahriar Nur Maulana Malik Ibrahim, Fikran Aranda Fahrudin 2025. Environmental DNA (e-DNA) extraction from soil. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp9295vzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 22, 2025

Last Modified: February 11, 2025

Protocol Integer ID: 118861

Keywords: WHO FPPL, e-DNA, Environmental DNA, Filamentous fungi, Metagenomic, Soil, extracting fungal genomic dna, resistant fungi, accurate responses to fungal threat, fungal genomic dna, antifungal resistance, fungal pathogen, antifungal resistance to major drug class, several pathogenic fungi, fungal threat, agricultural biosecurity, strategies against fungal infection, including fungi, fungal infection, including fusarium group, fusarium group, leveraging genomic technology, specific resistance gene, comprehensive genome database, genomic technology, environmental dna, diverse microorganism, agricultural soil, dna extraction, collecting soil sample, genomic dna, prevalent in thailand, eumycetoma agent, soil soil, extraction from soil soil, soil sample, soil, sustainable agriculture, significant human infection, network of mycologist, dna, fungi, agriculture

Funders Acknowledgements:

FAILSAFE Awarded projects - Seed Grant

Grant ID: FAILSAFE FR1-23

Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

Soil harbors diverse microorganisms, including fungi that play crucial environmental roles but can also harm plants, animals, and humans through infections or spore exposure. In 2022, the WHO released several pathogenic fungi, including Fusarium group, Mucorales and Eumycetoma agents, which are prevalent in Thailand's agricultural soils and have caused significant human infections. Despite this, the genetic mechanisms behind their infections and antifungal resistance remain poorly understood, particularly in Southeast Asia. Antifungal resistance to major drug classes, such as azoles, echinocandins, and polyenes, is increasing, driven by specific resistance genes. To address these challenges, we propose establishing a Southeast Asian (SEA) network of mycologists to study fungal pathogens, focusing on collecting soil samples, extracting fungal genomic DNA, and creating a comprehensive genome database. This initiative will improve detection and management of drug-resistant fungi, enhance clinical treatments, and strengthen agricultural biosecurity. By leveraging genomic technologies and fostering international collaboration, the project aims to develop targeted strategies against fungal infections and resistance, benefiting public health, food security, and sustainable agriculture across SEA. Establishing a fungal biobank and expert network will support long-term efforts to combat antifungal resistance, providing critical resources for swift and accurate responses to fungal threats.

Attachments



[Hand book DNeasy](#)

[Pow...](#)

470KB

Materials

1. Soil Sample: 0.5 g
2. Micropipette and Tips: 1000-10 μ L
3. DNeasy PowerSoil Pro Kit (QIAGEN)
4. Eppendorf Tubes
5. Agarose and Electrophoresis Set
6. TBE Buffer
7. DNA Ladder



Soil Sample Preparation:

- 1 Weigh 0.5 g of soil (remove debris like stones or large particles).

DNA Extraction Using the DNeasy PowerSoil Kit:

- 2 Add the soil sample into the PowerBead Pro Tube.
- 3 Add Solution CD1 to the tube as per kit instructions.
- 4 Homogenize the sample using one of the following Vortex Adapter method.

Cell Lysis

- 5 Ensure proper mixing and cell lysis through homogenization.

Inhibitor Removal

- 6 Add Solution CD2 to remove impurities.
- 7 Add Solution CD3, then load the sample into the MB Spin Column.

Washing

- 8 Wash the column with Solution EA followed by Solution C5 to ensure purity.

Elution

- 9 Elute the DNA using Solution C6 to collect the extracted DNA.



DNA Qualification and Quantification:

- 10 Perform electrophoresis on an agarose gel to confirm DNA quality.
- 11 Use a Nanodrop spectrophotometer to measure DNA concentration and purity.

Submission of Environmental DNA:

- 12 Send the extracted e-DNA to the FAILSAFE Project Team at Chulalongkorn University for further analysis.