

Dec 27, 2023 Version 2

Environmental DNA (eDNA) 12S Metabarcoding PCR Protocol (with Platinum SuperFi II Taq) V.2

 Version 1 is forked from [PCR Protocol Template](#)

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwdm39lmk/v2>

Kathleen Pitz¹, Jacoby Baker¹

¹MBARI

Better Biomolecular Oce...



Kathleen Pitz

Monterey Bay Aquarium Research Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwdm39lmk/v2>

Protocol Citation: Kathleen Pitz, Jacoby Baker 2023. Environmental DNA (eDNA) 12S Metabarcoding PCR Protocol (with Platinum SuperFi II Taq). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.e6nvwdm39lmk/v2> Version created by [Kathleen Pitz](#)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 27, 2023

Last Modified: February 06, 2024

Protocol Integer ID: 92774

Keywords: jbaker@mbari.org, mifish primer pcr protocol, version of the mifish primer pcr protocol, metabarcoding environmental dna, metabarcoding pcr protocol, dna 12s rrna gene, mitochondrial dna, environmental dna, rrna gene in eukaryote, rrna gene, hypervariable region of the mitochondrial dna, pcr inhibitor, pcr, platinum superfi ii taq, dna, mifish, fish, platinum superfi ii, resistance to pcr inhibitor

Funders Acknowledgements:

David and Lucile Packard Foundation

NASA

Grant ID: 80NSSC21M0032

NOAA

Grant ID: NA22NOS0120184



Abstract

The 12S protocol is aimed at amplifying the hypervariable region of the mitochondrial DNA 12S rRNA gene in eukaryotes. The primers (MiFish-U-F & MiFish-U-R) used in this protocol were developed by Miya et al., 2015 for metabarcoding environmental DNA (eDNA) from fishes.

This protocol follows an updated version of the MiFish primer PCR protocol. The Platinum SuperFi II was designed to have high fidelity and have increased resistance to PCR inhibitors.

MIOP: Minimum Information about an Omics Protocol

1

MIOP Term	Value
methodology category	omics analysis
project	Marine Biodiversity Observation Network (MBON)
purpose	PCR [OBI:0000415]
analyses	PCR [OBI:0000415]
geographic location	Monterey Bay [GAZ:00002509]
broad-scale environmental context	marine biome ENVO_00000447
local environmental context	oceanic epipelagic zone biome [ENVO:01000033]
environmental medium	sea water [ENVO:00002149] DNA extraction [OBI:0000257]
target	DNA extraction [OBI:0000257]
creator	Jacoby Baker, https://orcid.org/0000-0002-0673-7535
materials required	agarose gel electrophoresis system [OBI:0001134] PCR instrument [OBI:0000989]
skills required	sterile technique pipetting skills
time required	420
personnel required	1
language	en
issued	2023-11-14
audience	scientists
publisher	Monterey Bay Aquarium Research Institute, Chavez Lab
hasVersion	V.3
license	CC BY 4.0
maturity level	Mature

See <https://github.com/BeBOP-OBON/miop/blob/main/model/schema/terms.yaml> for list and definitions.

AUTHORS

2

PREPARED BY All authors known to have contributed to the preparation of this protocol, including those who filled in the template.	AFFILIATION	ORCID (visit https://orcid.org/)
Jacoby Baker	MBARI	0000-0002-0673-7535
N. Kobun Truelove	MBARI	0000-0002-2236-1000
Kathleen J. Pitz	MBARI	0000-0002-4931-1000
Francisco Chavez	MBARI	

MBARI : Monterey Bay Aquarium Research Institute, Moss Landing, CA

RELATED PROTOCOLS

3

PROTOCOL NAME AND LINK	ISSUER / AUTHOR	RELEASE DATE This is the date corresponding to the version listed to the left
https://mbari-bog.github.io/MBON-Protocols/eDNA_12S_SupFi2_PCR_V3.html	Jacoby Baker	2023-11-07
https://mbari-bog.github.io/MBON-Protocols/Bead_cleanup.html	Jacoby Baker	2023-11-07
Environmental DNA (eDNA) 12S Metabarcoding Illumina MiSeq NGS Protocol with size selection V.2 ; dx.doi.org/10.17504/protocols.io.eq2lyj1kelx9/v2	Jacoby Baker	2023-12-27

This is a list of other protocols which should be known to users of this protocol. Please include the link to each related protocol.

ACRONYMS AND ABBREVIATIONS

4

ACRONYM / ABBREVIATION	DEFINITION
eDNA	environmental DNA
NTC	No Template Control

GLOSSARY

5

SPECIALISED TERM	DEFINITION
Content Cell	Content Cell
Content Cell	Content Cell

BACKGROUND

6 **Summary**

The 12S protocol is aimed at amplifying the hypervariable region of the mitochondrial DNA 12S rRNA gene in eukaryotes. The primers (MiFish-U-F & MiFish-U-R) used in this protocol were developed by Miya et al., 2015 for metabarcoding environmental DNA (eDNA) from fishes.

This work was supported by the David and Lucile Packard Foundation, and NASA award 80NSSC21M0032 and NOAA award NA22NOS0120184 in support of the CeNCOOS MBON.

7 **Method description and rationale**

This protocol follows an updated version of the MiFish primer PCR protocol. The Platinum SuperFi II was designed to have high fidelity and have increased resistance to PCR inhibitors.

8 Spatial coverage and environment(s) of relevance

This protocol has been used to amplify extracted DNA from filtered sea water samples taken from marine coastal stations off the western coast of North America (primarily off of California).

sea water [ENVO:00002149]

http://purl.obolibrary.org/obo/ENVO_00002149

9 Personnel Required

1 technician.

10 Safety

Identify hazards associated with the procedure and specify protective equipment and safety training required to safely execute the procedure

11 Training requirements

Sterile technique, pipetting skills.

12 Time needed to execute the procedure

Total time is 7 hours.

PCR preparation and running the PCR protocol takes 3 hours. Running the following gel is 1 hour, bead cleanup setup preparation and process takes 2 hours, and then another gel is run for 1 hour.

EQUIPMENT

13

DESCRIPTION e.g. filter	PRODUCT NAME AND MODEL Provide the official name of the product	MANUFACTURER Provide the name of the manufacturer of the product
Durable equipment		
Agarose gel electrophoresis system		
PCR Thermal Cycler		
Consumable equipment		
PCR plates	SuperPlate PCR Plate, 96-well, semi-skirted	ThermoFisher Scientific
Plate seals	PCR Plate Seals	Bio Rad
Chemicals		
2X Platinum SuperFi II PCR Master Mix	2X Platinum SuperFi II PCR Master Mix	ThermoFisher Scientific

STANDARD OPERATING PROCEDURE

14 In the following SOP, please use the exact names of equipment as noted in the table above.

Provide a step-by-step description of the protocol. The identification of difficult steps in the protocol and the provision of recommendations for the execution of those steps are encouraged.

PREPARATION

15 BEFORE STARTING

Disinfect work surfaces with 10% bleach or RNase Away followed by a MilliQ / DI water rinse and 70% ethanol wipe. Clean pipet surfaces with RNase Away and ethanol wipe.

UV pipets, molecular grade water, and tube racks for 30 minutes prior to starting protocol.

PCR

- 16 1. eDNA template & PCR processing were performed at the Monterey Bay Aquarium Research Institute (MBARI). PCR reactions for the 12S locus were performed with a two-step amplification protocol for each sample using the MiFish_U primers (Miya et al. 2015) with Fluidigm adapters CS1 & CS2. All primers listed in the 5' to 3' direction. MiFish primers are in bold.

Fluidigm CS1+**MiFish_U (forward):**

ACACTGACGACATGGTTCTACA **GTCGGTAAAACTCGTGCCAGC**

Fluidigm CS2+**Mifish_U (reverse):**

TACGGTAGCAGAGACTTGGTCT **CATAGTGGGGTATCTAATCCCAGTTTG**

PCR Primer Name	Direction	Sequence (5' → 3')
Fluidigm CS1 + 12S MiFish_U	forward	ACACTGACGACATGGTTCTA CAGTCGGTAAAACTCGTGCC AGC
Fluidigm CS2 + 12S MiFish_U	reverse	TACGGTAGCAGAGACTTGGT CTCATAGTGGGGTATCTAATC CCAGTTTG

- 17 2. The primary PCR amplifications were carried out in singleton 50 µl reactions using:
- 25 µL 2X Platinum SuperFi II PCR MM
 - 3 µL forward primer (10 µM)
 - 3 µL reverse primer (10 µM)
 - 3 µL eDNA extract template
 - 16 µL molecular-grade, nuclease-free water
- 18 3. PCR reactions were performed in 96-well plates with a no-template control (NTC) for each PCR plate, for a total of 3 PCR negative controls. An artificial community was used as a positive control.
- 19 4. Primary 12S cycling parameters, using the manufacturer's recommendation with a slight modification. The Platinum SuperFi II is designed to be used with a universal 60°C annealing temperature, however, we found we had better results with an annealing temperature of 62°C to reduce bacterial co-amplification:

PCR step	Temperature	Duration	Repetition
denaturation	98°C	30 seconds	1
denaturation	98°C	10 seconds	38 cycles
annealing	62°C	10 seconds	38 cycles
extension	72°C	30 seconds	38 cycles
final extension	72°C	5 minutes	1
HOLD	4°C	HOLD	1

Quality Control and Product Clean-up

- 20 1. After primary PCR amplification of the marker region the PCR product was run through a 2% agarose gel to confirm the presence of target bands (~270 bp) and absence of non-specific amplification (bacterial band ~370 bp) across environmental samples.

Primary PCR clean-up (**Bead Clean-up Protocol**)



- 21 2. Primary PCR products were purified and size selected using the Agencourt AMPure XP bead system (Beckman Coulter, USA) at 0.9x volume beads to product.
- 22 3. A second agarose gel was run to confirm primer removal and retention of target amplicons after purification. NTCs were also tested using a Qubit dsDNA 1x high sensitivity kit to ensure no amplification.

Secondary Amplification

23 The following steps are performed by MSU's RTSF Genomics Core

1. Secondary amplification and NGS were performed at Michigan State University's Research Technology Support Facility (RTSF). An aliquot of 20 µL from each purified primary PCR product was sent to RTSF Genomics Core at MSU for secondary PCR amplification with primers which targeted the CS1/CS2 ends of the primary PCR products and added dual indexed, Illumina compatible adapters with barcodes.

- PE1-BC-CS1 (forward):

```
AATGATACGGCGACCACCGAGATCT-[i5-BC(index 2)]-ACACTGACGACATGTTCTACA
```

- PE2-BC-CS2 (reverse):

```
CAAGCAGAAGACGGCATACGAGAT-[i7-BC(index 1)]-TACGGTAGCAGAGACTTGGTCT
```

PCR Primer Name	Direction	Sequence (5' → 3')
PE1-BC-CS1	forward	AATGATACGGCGACCACCGA GATCT-[i5-BC(index 2)]- ACACTGACGACATGTTCTA CA
PE2-BC-CS2	reverse	CAAGCAGAAGACGGCATACG AGAT-[i7-BC(index 1)]- TACGGTAGCAGAGACTTGGT CT

- 24 2. The secondary PCR amplifications were carried out in 15 µL reactions, using 1 µL of primary PCR product.
 - 6 µl 2.5X HotMaster Mix
 - 7 µl DI water
 - 1 µl Primer Mix (6uM)
 - 1 µl primary eDNA PCR product
- 25 3. Secondary 12S cycling parameters:

PCR step	Temperature	Duration	Repetition
denaturatio n	95°C	3 minutes	1
denaturatio n	95°C	15 seconds	15 cycles
annealing	60°C	30 seconds	15 cycles
extension	72°C	1 minute	15 cycles
final extension	72°C	3 minutes	1
HOLD	25°C	HOLD	1

QUALITY CONTROL



- 26 An agarose gel was run after secondary PCR to confirm the presence of target bands and absence of non-specific amplification across environmental samples as well as the absence of amplification in NTCs.

Products from the protocol are then used to create a pooled library and sequenced following a separate sequencing protocol (such as "Environmental DNA (eDNA) 12S Metabarcoding Illumina MiSeq NGS Protocol with size selection V.2", dx.doi.org/10.17504/protocols.io.eq2lvj1kelx9/v2)

BASIC TROUBLESHOOTING GUIDE

- 27 Identify known issues associated with the procedure, if any.

Provide troubleshooting guidelines when available.

REFERENCES

- 28
1. Miya M et al. 2015 MiFish, a set of universal PCR primers for metabarcoding environmental DNA from fishes: detection of more than 230 subtropical marine species. R.Soc.opensci. 2: 150088. <http://dx.doi.org/10.1098/rsos.150088>
 2. Kawato, M., Yoshida, T., Miya, M., Tsuchida, S., Nagano, Y., Nomura, M., Yabuki, A., Fujiwara, Y. and Fujikura, K., 2021. Optimization of environmental DNA extraction and amplification methods for metabarcoding of deep-sea fish. MethodsX, 8, p.101238. <https://doi.org/10.1016/j.mex.2021.101238>
 3. **Platinum SuperFi II PCR Master Mix User Guide**

APPENDIX A: DATASHEETS

- 29 Link templates (e.g. preformatted spreadsheets) used to record measurements and report on the quality of the data as well as any documents such as manufacturer specifications, images, etc that support this protocol. Please include a short note describing the document's relevance.