



Apr 15, 2019

Enzo's CGH Labeling Kit for Oligo Arrays

 [PLOS One](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.zjuf4nw>

Hendrik F van Essen¹

¹Amsterdam UMC



Hendrik F HFvE Van Essen

Amsterdam UMC

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.zjuf4nw>

External link: <https://doi.org/10.1371/journal.pone.0223827>

Protocol Citation: Hendrik F van Essen 2019. Enzo's CGH Labeling Kit for Oligo Arrays. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.zjuf4nw>

Manuscript citation:

Vincenten JPL, Essen HFv, Lissenberg-Witte BI, Bulkman NWJ, Krijgsman O, Sie D, Eijk PP, Smit EF, Ylstra B, Thunnissen E (2019) Clonality analysis of pulmonary tumors by genome-wide copy number profiling. PLoS ONE 14(10): e0223827. doi: [10.1371/journal.pone.0223827](https://doi.org/10.1371/journal.pone.0223827)



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 27, 2019

Last Modified: April 15, 2019

Protocol Integer ID: 21844

Keywords: cgh labeling kit

Materials

STEP MATERIALS

 PB binding buffer **Qiagen Catalog #19066**

Protocol materials

 PB binding buffer **Qiagen Catalog #19066**



Denature DNA and anneal random primers

1 pipette the following in a nuclease free thermocycler strip:

- a. DNA X μ l
- b. water 19 – X μ l (Vial W)
- b. Primers/Reaction buffer 20 μ l (Vial 1)

Total volume: 39 μ l

2 Add cap, flick to mix contents and briefly spin down

3 Heat at  99 °C for  00:10:00

4 Place in a cooled 96 well rack on ice  4 °C for  00:05:00

5 Centrifuge briefly

Extend Primers with Klenow Exo-DNA Polymerase

6 Add, while samples are on ice  4. °C :

- a. To test sample 10 μ l Cy3 (Vial 2)
- b. To reference sample 10 μ l Cy5 (Vial 3)
- c. Klenow Exo-DNA Pol 1 μ l (Vial 4)
- Total 50 μ l

7 Add cap, flick to mix contents and briefly spin down

8 Incubate at  37 °C  04:00:00 for

9 Add 5 μ l of Stop Buffer (Vial 5), mix, and briefly centrifuge

Removal of uncoupled nucleotides

- 10 Transfer labelled DNA's to 1.5 ml reaction tubes and keep Cy3 and Cy5 labelled DNA's separated
- 11 Add 275 µl of Binding Buffer (PB)
 PB binding buffer **Qiagen Catalog #19066**
- 12 Load column and centrifuge for 60 seconds at 16,000g and discard flow through
- 13 Add 650 µl of Wash Buffer (PE), centrifuge 60 seconds at 16,000g and discard flow through
- 14 Centrifuge for 60 seconds at 16,000g, discard centrifuge tube and place column in fresh 1.5 ml centrifuge tube
- 15 Add 10.5 µl of Elution Buffer (EB) (middle of the filter) and incubate for 1 minute.
- 16 Centrifuge for 60 seconds at 16,000g
- 17 Add 10.5 µl of Elution Buffer (EB) (middle of the filter) and incubate for 1 minute. .
- 18 Measure yield and specific dye incorporation with NanoDrop, Program MicroArray. Labelled, cleaned up gDNA can be stored at -20°C in the dark for 7 days.