


Jul 17, 2020

Version 2

Enzymatic Assay of Protease Using Azocasein as Substrate V.2

DOI

dx.doi.org/10.17504/protocols.io.bhqnj5ve



Neilier Junior¹

¹University of Manitoba



Neilier Junior

University of Manitoba, Universidade Federal de Viçosa

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bhqnj5ve>

Protocol Citation: Neilier Junior 2020. Enzymatic Assay of Protease Using Azocasein as Substrate. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bhqnj5ve>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 19, 2020


Last Modified: July 17, 2020

Protocol Integer ID: 38382

Keywords: Enzyme, substrate, kinetic, enzymology, protein, enzymatic assay of protease, azocasein as substrate, enzymatic assay, protease, using azocasein, assay

Materials

MATERIALS

 Calcium chloride **Merck Millipore (EMD Millipore) Catalog #1.02378.0500**

 Trichloroacetic acid (TCA) **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T6399**

 Sodium hydroxide **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S8045**

 Trizma[®] base **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T4661**

 Azocasein **Catalog #A2765**

Troubleshooting

Safety warnings

 Wear personal protective equipment: gloves, lab coat and mask.

Before start

Organize your workspace

Make sure all solutions and equipment are available.



Reagent Preparation

- 1
 - 100 mM Tris-HCl buffer, pH 8.0, 20 mM CaCl₂, at 37 °C.
 - 2.0% (w/v) Azocasein Solution
Heat gently (do not boil) to 50 - 60 °C for 10 min with stirring.
Adjust the pH to 8.0 at 37 °C, if necessary, with either 1.0 M NaOH or 1.0 M HCl.
 - 110 mM Trichloroacetic Acid Reagent (TCA). Dilute with deionized water.
 - 500 mM Sodium Hydroxide (NaOH) Solution. Prepare in deionized water.

Check how many samples will be analyzed to calculate the required volume of each solution to be prepared.

Procedure

- 2
Pipette (in microliters) the following reagents into 2.0 mL microtubes.

	Blank	Test
Tris-HCl buffer	750 μL	450 μL
Azocasein	750 μL	750 μL
<i>Mix and equilibrate to the at desired temperature. Then add:</i>	*	
Sample (enzyme source)	-	300 μL
<i>Mix and incubate at desired temperature for exactly 30 min.</i>	*	
<i>Remove a 1 mL aliquot from both (test and blank) solutions and place into 2.0 mL microtubes. Then add:</i>	*	
TCA	1000 μL	1000 μL
<i>Centrifuge at 20,000 g for 10 min. Remove a 1 mL aliquot from supernatant (test and blank) and place into 2.0 mL microtubes. Then add:</i>	*	

	NaOH	1000 μL	1000 μL
	<i>Mix and transfer the Test and Blank solutions to suitable cuvettes. Measure the A_{440nm} for Test and Blank using a spectrophotometer.</i>	*	

Calculation

3 $\Delta A_{440nm} = A_{440nm}^{Test} - A_{440nm}^{Blank}$