


Jul 16, 2020 Version 2

# Enzymatic Assay of Protease Using Azocasein as Substrate V.2

DOI

[dx.doi.org/10.17504/protocols.io.bhqnj5ve](https://dx.doi.org/10.17504/protocols.io.bhqnj5ve)



Neillier Junior<sup>1</sup>

<sup>1</sup>University of Manitoba



**Neillier Junior**

University of Manitoba, Universidade Federal de Viçosa

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bhqnj5ve](https://dx.doi.org/10.17504/protocols.io.bhqnj5ve)

**Protocol Citation:** Neillier Junior 2020. Enzymatic Assay of Protease Using Azocasein as Substrate. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bhqnj5ve>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** June 19, 2020

**Last Modified:** July 16, 2020

**Protocol Integer ID:** 38382

**Keywords:** Enzyme, substrate, kinetic, enzymology, protein



## Materials

### MATERIALS

⊠ Calcium chloride **Merck Millipore Catalog #1.02378.0500**

⊠ Trichloroacetic acid (TCA) **Sigma – Aldrich Catalog #T6399**

⊠ Sodium hydroxide **Sigma – Aldrich Catalog #S8045**

⊠ Trizma® base **Sigma Aldrich Catalog #T4661**

⊠ Azocasein **Catalog #A2765**

## Safety warnings

! Wear personal protective equipment: gloves, lab coat and mask.

## Before start

Organize your workspace

Make sure all solutions and equipment are available.

## Reagent Preparation

- 1
  - 100 mM Tris-HCl buffer, pH 8.0, 20 mM CaCl<sub>2</sub>, at 37 °C.
  - 2.0% (w/v) Azocasein Solution  
Heat gently (do not boil) to 50 - 60 °C for 10 min with stirring.  
Adjust the pH to 8.0 at 37 °C, if necessary, with either 1.0 M NaOH or 1.0 M HCl.
  - 110 mM Trichloroacetic Acid Reagent (TCA). Dilute with deionized water.
  - 500 mM Sodium Hydroxide (NaOH) Solution. Prepare in deionized water.

**Check how many samples will be analyzed to calculate the required volume of each solution to be prepared.**

## Procedure

- 2  
Pipette (in microliters) the following reagents into 2.0 mL microtubes.

	Blank	Test
<b>Tris-HCl buffer</b>	750 μL	450 μL
<b>Azocasein</b>	750 μL	750 μL
<i>Mix and equilibrate to the at desired temperature. Then add:</i>	*	
<b>Sample (enzyme source)</b>	-	300 μL
<i>Mix and incubate at desired temperature for exactly 30 min.</i>	*	
<i>Remove a 1 mL aliquot from both (test and blank) solutions and place into 2.0 mL microtubes. Then add:</i>		
<b>TCA</b>	1000 μL	1000 μL
<i>Centrifuge at 20,000 g for 10 min. Remove a 1 mL aliquot from supernatant (test and blank) and place into 2.0 mL microtubes. Then add:</i>	*	
<b>NaOH</b>	1000 μL	1000 μL
<i>Mix and transfer the Test and Blank solutions to suitable cuvettes. Measure the A440nm for Test and Blank using a spectrophotometer.</i>	*	



## Calculation

$$3 \quad \Delta A_{440\text{nm}} = A_{440\text{nm}}^{\text{Test}} - A_{440\text{nm}}^{\text{Blank}}$$