


Nov 14, 2019

Enzymatic Assay of Trypsin Inhibition

DOI

dx.doi.org/10.17504/protocols.io.9a4h2gw



Neillier Junior¹

¹University of Manitoba



Neillier Junior

University of Manitoba, Universidade Federal de Viçosa

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.9a4h2gw>

Protocol Citation: Neillier Junior 2019. Enzymatic Assay of Trypsin Inhibition. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.9a4h2gw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 14, 2019

Last Modified: November 14, 2019

Protocol Integer ID: 29756

Keywords: Enzymatic assay, Trypsin inhibitor, enzymatic assay of trypsin inhibition, trypsin inhibition, analysis of protease inhibitor, protease inhibitor, soybean leaf extract, enzymatic assay, other biological sample


Abstract

This protocol has been standardized for analysis of protease inhibitors in soybean leaf extract, but can be easily adjusted for other biological samples.

Materials


MATERIALS

 Tris

 DMSO Merck MilliporeSigma (Sigma-Aldrich) Catalog #D1435

 Calcium Chloride Merck MilliporeSigma (Sigma-Aldrich) Catalog #C4904

 Trypsin from bovine pancreas Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8003

 N α -Benzoyl-L-arginine 4-nitroanilide hydrochloride (L-BApNA) Merck MilliporeSigma (Sigma-Aldrich) Catalog #B3279

Troubleshooting

Before start

Check that all reagents and equipment are available. Plan the experiment!

Prepare the solutions and the workspace

1 Preparation of solutions

Trypsin solution:

Dilute 1.25 mg of bovine trypsin in 1 mL of water pH 3.0 (adjust with HCl).

L-BApNA stock solution:

Dilute 130.47 mg L-BApNA in 5 mL DMSO (concentration: 6.0×10^{-2} M). Store at -20°C and protected from light.

L-BApNA usage solution (freshly prepared):

Dilute 200 μL of stock solution in 10 mL of 100 mM Tris-HCl buffer, pH 8.2 and 20 mM CaCl_2 . Protect from light.

2 Separate three microtubes and name them "**blank**", "**control** (uninhibited test)" and "**test** (inhibited test)".

Pipette the following reagents.

Blank: 500 μL BA ρ NA usage solution and 500 μL buffer.

Control: 100 μL trypsin solution and 400 μL buffer.

Test: 100 μL enzyme, 100 μL leaf extract (source of inhibitors) and 300 μL buffer.

3 Mix the three microtubes by inversion and equilibrate to 25°C for 5 min

Zero spectrophotometer with **blank** content at 410 nm

4 To the **control** microtube, add 500 μL of the BA ρ NA usage solution

Immediately mix by inversion and mark the time and pour the contents into a cuvette.

After 30 s of reaction onset, monitor readings at 410 nm for 120 s

5 Add 500 μL of the BA ρ NA usage solution to the **test** tube

Immediately mix by inversion and mark the time and pour the contents into a cuvette.

After 30 s of reaction onset, monitor readings at 410 nm for 120 s

6 Calculations

$$\% \text{ Inhibition} = (\Delta A_{410 \text{ nm}}^{\text{control}} - \Delta A_{410 \text{ nm}}^{\text{test}}) * 100 / (\Delta A_{410 \text{ nm}}^{\text{control}} - A_{410 \text{ nm}}^{\text{blank}})$$

or

Trypsin Inhibitor Units / mL = $(\Delta A_{410 \text{ nm}}^{\text{control}} - \Delta A_{410 \text{ nm}}^{\text{test}}) / (8800 * \text{time} * \text{leaf extract volume})$

$A_{410 \text{ nm}}^{\text{blank}}$ = Absorbance in the **blank** at 410 nm

$\Delta A_{410 \text{ nm}}^{\text{control}}$ = Absorbance variation in the **control** sample at 410 nm within 120 seconds


$\Delta A_{410 \text{ nm}}^{\text{test}}$ = Absorbance variation in the **test** sample at 410 nm within 120 seconds

8800 = extinction coefficient of *p*-nitroanilide at 410 nm

time = 120 seconds

leaf extract volume = Volume of inhibitor source used (in milliliters)

The presence of trypsin inhibitors in the leaf extract decreases the enzymatic activity and this inhibition can be represented in the Lineweaver-Burk graph, where, as the inhibitor concentration increases, the slope of the line also increases.



Lineweaver-Burk plot analysis of the inhibitory activity of soybean leaf extract toward trypsin. In blue, the kinetics in the absence of inhibitors (**control**). The kinetics in the presence of inhibitors (**test**) at 0.5 K_i and K_i are represented in orange and gray, respectively.