

Feb 04, 2020

# Enzymatic disaggregation of human myometrium for 10X Single Cell RNA-seq

DOI

<https://dx.doi.org/10.17504/protocols.io.bb5miq46>

Aymara Mas<sup>1</sup>, Alba Machado-Lopez<sup>1</sup>, Patricia Escorcía<sup>1</sup>, Carlos Simón<sup>1,2,3</sup>

<sup>1</sup>Igenomix Foundation, INCLIVA, Valencia, Spain; <sup>2</sup>Valencia University, Valencia, Spain;

<sup>3</sup>Harvard University, Boston, USA



Aymara Mas

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.bb5miq46>

**Protocol Citation:** Aymara Mas, Alba Machado-Lopez, Patricia Escorcía, Carlos Simón 2020. Enzymatic disaggregation of human myometrium for 10X Single Cell RNA-seq. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bb5miq46>

### Manuscript citation:

Mas A, et al. Identification and characterization of the human leiomyoma side population as putative tumor-initiating cells. *Fertil Steril*. 2012 Sep;98(3):741-751.

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Working

**Working.** We and other researchers are using and published the results using this protocol (Mas et al., 2015; Prusinski et al., 2018; Corachán et al., 2019).

**Created:** February 04, 2020

**Last Modified:** February 04, 2020

**Protocol Integer ID:** 32653

**Keywords:** enzymatic disaggregation of human myometrium, dissociating myometrial sample, enzymatic disaggregation, myometrial sample, human myometrium, 10x single cell rna, using collagenase iv, collagenase iv, fertil steril, seq this protocol, seq, single cell

## Abstract

This protocol describes the procedure for dissociating myometrial samples, which is based in an enzymatic disaggregation using Collagenase IV and DNase I. This protocol is adapted from Mas et al. 2012 *Fertil Steril* with minor modifications.

## Guidelines

### Storage conditions of reagents

|  | Reagent             | Storage condition                                     |
|--|---------------------|-------------------------------------------------------|
|  | Hypothermosol       | 4°C                                                   |
|  | HBS S               | 4°C                                                   |
|  | HEPES 1M            | 4°C                                                   |
|  | DNAse type I        | Resuspended in PBS and store 100 µl aliquots at -20°C |
|  | Collagenase type IV | 4°C                                                   |
|  | DMEM low glucose    | 4°C                                                   |
|  | FBS                 | Store 30 ml aliquots at -20°C                         |
|  | Penicillin - Stre   | Store 5 ml aliqu                                      |



|  |                                             |                        |
|--|---------------------------------------------|------------------------|
|  | pto<br>myci<br>n<br>(100<br>00<br>U/ml<br>) | ots<br>at<br>-20°<br>C |
|--|---------------------------------------------|------------------------|

**Required equipment**

|  | <b>Equi<br/>pme<br/>nt</b>                                                      | <b>Sup<br/>plier</b>         | <b>Cata<br/>log<br/>n°</b> |
|--|---------------------------------------------------------------------------------|------------------------------|----------------------------|
|  | Hera<br>cell<br>150i<br>CO2<br>incu<br>bato<br>r                                | Ther<br>mo<br>Scie<br>ntific | 5102<br>628<br>0           |
|  | New<br>Brun<br>swic<br>k<br>Incu<br>bato<br>r<br>shak<br>er<br>Inno<br>va<br>40 | Epp<br>end<br>orf            | M12<br>99-<br>009<br>2     |
|  | Cent<br>rifug<br>e<br>Meg<br>a<br>Star<br>600                                   | VWR<br>colle<br>ction        | 521-<br>1893               |
|  | Micr<br>osco<br>pe<br>Niko<br>n<br>eclip<br>se<br>TE2<br>00                     | Niko<br>n                    | -                          |
|  | Clas<br>s II<br>micr<br>obiol<br>ogic<br>al<br>safe                             | Ther<br>mo<br>Scie<br>ntific | 5102<br>5411               |



|  |                                               |  |  |
|--|-----------------------------------------------|--|--|
|  | ty<br>cabi<br>net<br>MSC<br>adva<br>ntag<br>e |  |  |
|--|-----------------------------------------------|--|--|

**The protocol workflow is as follows:**

1. Tissue collection and cold transportation until processing
2. Dissociation: mechanic and enzymatic
3. Remove red blood cells if necessary
4. Prepare cells for Chromium

**Materials****MATERIALS**

- ✕ BRAND® counting chamber BLAUBRAND® Neubauer improved New without clips, double ruled **Merck MilliporeSigma (Sigma-Aldrich) Catalog #BR717805**
- ✕ Penicillin-Streptomycin **Gibco - Thermo Fisher Scientific Catalog #15140122**
- ✕ Collagenase Type 4 **Worthington Biochemical Corporation Catalog #LS004188**
- ✕ 100 µm Cell Strainer **Falcon Catalog #352360**
- ✕ HypoThermosol® FRS Preservation Solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H4416**
- ✕ TrypLE® Select Enzyme (1X), no phenol red **Thermo Fisher Catalog #12563029**
- ✕ HEPES (1 M) **Thermo Fisher Catalog #15630080**
- ✕ Nalgene® Rapid-Flow® Sterile Disposable Filter Units with PES Membrane, 250mL, 0.2µm pore, 50mm membrane **Thermo Fisher Catalog #568-0020**
- ✕ ACK Lysing Buffer **Thermo Fisher Catalog #A1049201**
- ✕ HBSS (1x) **Gibco - Thermo Fisher Scientific Catalog #14175-095**
- ✕ Deoxyribonuclease I from bovine pancreas **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D4513**
- ✕ Dulbecco's Modified Eagle's Medium - low glucose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D5921**
- ✕ Fetal bovine serum (FBS) **BioWest Catalog #S181B-500**
- ✕ Falcon 40 µm Cell Strainer **Corning Catalog #352340**
- ✕ Sterile surgical blades **Aesculap Catalog #7806607**



## Safety warnings

! Human samples including tissue, blood and bodily fluids have the potential to harbour HG2 and Hazard Group 3 (HG3) organisms, specifically Blood Borne Viruses (BBVs). As security warning, please follow the procedures related to Biological Safety at Containment Level 2.

## Before start

### Prepare Wash buffer:

Make stock of wash buffer and store at 4 °C. Aliquot 49.5 mL of HBSS in a 50 ml conical Falcon tube and add 500 µl of Pen/Strep (100 U/ml).

### Prepare Enzyme buffer:

To make 100 ml of enzyme buffer, combine 95.6 ml of HBSS with 1 ml of Pen/Strep (100 U/ml); 1 ml of HEPES (1M); 4.6 ml of wash buffer; 95 µl of DNase and 135 mg of collagenase IV and filter in a 0.2 µm bottle-top sterile filter unit.

### Prepare resuspension buffer:

Add 12 ml of inactivated FBS and 1 ml of Pen/Strep (10000 U/ml) to 87 ml of DMEM medium (low glucose) to prepare 100 ml of resuspension buffer.

## Tissue collection and transportation

- 1 Uterine tissues should be transported in preservation solution (HypoThermosol® FRS) at  4 °C from the surgery room to the lab.

## Dissociation

- 2 Rinse the sample with wash buffer, removing blood or mucus.
- 3 Place wet tissue under a petri dish. Isolate with a sterile scalpel the myometrial layer by removing the endometrium, serosa and laser-burnt zones.
- 4 Mechanical disaggregation by chopping and thoroughly mincing up the tissue into small pieces (<1 mm<sup>3</sup>).
- 5 Transfer contents to 50 ml falcon tubes containing  30 mL of enzyme buffer. Tighten lid, seal with parafilm and incubate at  37 °C overnight horizontally.
- 6 Afterwards, cell suspensions should be filtered through 100 to 50-µm polyethylene filters to remove cellular clumps and undigested tissue.
- 7 Centrifuge the filtered medium  400 x g, 00:05:00 . Remove supernatant and add at least  1 mL of resuspension buffer (depending of the pellet content).

## Remove red blood cells

- 8 [Optional] If after centrifugation there is a ring of red blood cells, add ACKL lysis buffer, consisting in hypotonic shock, and incubate at  37 °C for  00:05:00 . Then, follow step 7 again.

## Dissociation

- 9 Add  400 µL Tryple Select Enzyme to resuspended media containing cells and mix thoroughly by pipetting. Incubate for  00:10:00 -  00:15:00 at  37 °C and



centrifuge at  300 rpm .

10 Add  100  $\mu$ L DNase I to digest the extracellular genomic DNA and mix thoroughly by pipetting. Incubate for  00:05:00 at  Room temperature .

11 Add at least  1 mL of resuspension buffer. Filter through a 40  $\mu$ m cell strainer and centrifuge at  400 x g, 00:05:00 .

12 Remove supernatant and resuspend the content in at least  1 mL of resuspension buffer (depending of the pellet content).

## Prepare cells for chromium

13 Count the total number of cells (at least twice on two different cell counters).

14 Proceed to 10x experiment.

15 [Optional] Resuspend cells in freezing medium to a concentration of  $1 \times 10^6$  cells for storage in cryogenic vials.