




Jun 16, 2023

Version 1

🌐 Enzymatic fragmentation of plant chromatin for Hi-C libraries V.1

DOI

dx.doi.org/10.17504/protocols.io.j8nlk4pnxg5r/v1



ignacio.carvajal Carvajal¹, Elena Hilario¹

¹Plant and Food Research

Protist Research to Opti...

High molecular weight D...

2 more workspaces



Elena Hilario

Plant and Food Research

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk4pnxg5r/v1>

Protocol Citation: ignacio.carvajal Carvajal, Elena Hilario 2023. Enzymatic fragmentation of plant chromatin for Hi-C libraries. protocols.io <https://dx.doi.org/10.17504/protocols.io.j8nlk4pnxg5r/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 06, 2021

Last Modified: June 16, 2023

Protocol Integer ID: 52149

Keywords: chromatin conformation capture, Hi-C, DNase I, nuclei, plant, restriction enzyme, enzymatic fragmentation of plant chromatin, plant chromatin sample, plant chromatin, chromatin sample, enzymatic fragmentation, sample ready for ngs library preparation, ngs library preparation, several quality checks along the process

Abstract

The aim of this protocol is to learn how to prepare, evaluate and optimize a plant chromatin sample to produce a proximity ligated sample ready for NGS library preparation. We introduce several quality checks along the process to avoid wasting precious material or time working with a suboptimal sample (quality and quantity-wise).

Image Attribution

New Plymouth Coastal Walkway, Taranaki, New Zealand (Elena Hilario, 2021)



Guidelines

Every sample brings a different challenge and no chromatin protocol is universal. It is recommended to practice with a less precious sample before attempting the preparation of the final Hi-C library. There are a few commercial kits available to prepare Hi-C libraries and perform well under most circumstances, especially since often they have good technical support to help you along the way. This protocol describes an in vitro enzymatic fragmentation of plant chromatin and aims to help understand the importance of QC steps along the process.

Ideally, the chromatin fragmentation should be random, using DNase I, but in practice, this approach can take some time to optimize and if your sample size is limited you may consider using a 4-base **restriction enzyme**, selected based on the genome's GC content. Or you could also use a mix of two or more 4-base restriction enzymes to target as many sites as possible, either by having a different restriction site or a different methylation sensitivity. Although it is easier to control a restriction enzyme digestion, beware that the genome coverage is reduced, compared to a randomly fragmented chromatin sample. We strongly recommend to perform digestion tests on nuclear gDNA preps with several restriction enzymes until you find a combination that will produce an even digestion profile.

Using a biotinylated bridge to mark the contact ends adds extra steps and some assurance that the contacts detected are real, but you could also omit this step and fill in the 5'-end overhangs with biotinylated dNTPs. Beware the length of the carbon chains that links the biotin to the nucleotide (the standard is 6 carbon atoms, used in this protocol) has an impact on how well it will bind to the streptavidin molecule. The longer the carbon chain, the more efficient the binding.

If you have trouble isolating a clean nuclei pellet, you could purify it with a Percoll gradient, as described **here**. The wash volumes of the chromatin bound to magnetic beads can be increased if the sample is hard to clean, or extra washes can be implemented too. The lysis buffer volume can also be increased to reduce the amount of contaminants even before starting the chromatin washes.



Materials

~ 3 g of leaves **ground with liquid nitrogen** used immediately or previously ground and stored in a 50 mL Falcon tube at -80°C

Sorbitol Wash Buffer: 100 mM Tris-HCl pH 8.0, 0.35 M Sorbitol, 5 mM EDTA pH 8.0, 1 % (w/v) Polyvinylpyrrolidone Molecular wt. 40,000 (PVPK40). Autoclave and store at 4°C, it will last for at least 6 months. Add β -mercaptoethanol (1 % v/v) before the extraction

NEB complete buffer: 0.5 M Mannitol, 10 mM PIPES-KOH pH 6, 10 mM $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$, 2% PVP K40, 200 mM L-lysine monohydrochloride, 6 mM EGTA. Add Triton X-100 (final concentration 0.5% v/v), sodium bisulfite (final concentration 18 mM) and β -mercaptoethanol (final concentration 0.04% v/v) just before use

NEB- β ME: 0.5 M Mannitol, 10 mM PIPES-KOH pH 6, 10 mM $\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$, 2% PVP K40, 200 mM L-lysine monohydrochloride, 6 mM EGTA. Add sodium bisulfite (final concentration 18 mM) just before use

1X PBS buffer, sterile. 50 mL

10X Wash Buffer: 100 mM Tris-HCl pH 8, 1 M NaCl, sterile. Add Tween 20 to a final concentration of 0.5% after autoclaving and store at 4°C

1X Wash Buffer: Prepare at least 50 mL by diluting the 10X stock with sterile deionized water and store at 4°C

Quenching Solution: 1X Wash Buffer

FA buffer 50 mM HEPES-KOH pH 7.5, 140 mM NaCl, 1 mM EDTA pH 8, 1% Triton X-100, 0.1% sodium deoxycholate, 0.1% SDS. Prepare the buffer without the detergents, autoclave, and when cooled to room temperature, add the detergents. Store at 4°C, in the dark

1X CutSmart 1% SDS: Prepare 1 mL with sterile deionized water

100X TE pH 7.5: 1 M Tris pH 7.5, 100 mM EDTA, autoclave and keep at room temperature

80% ethanol, freshly prepared, 10 mL

Deionized sterile water (ddH_2O)

TE Buffer: Dilute 100X TE pH 7.5 1:100 with sterile deionized sterile water in a 50 mL sterile Falcon tube

10 mM Tris-HCl pH 8.0, sterile

0.5 M EDTA pH 8.0, sterile



5 M NaCl, sterile

DBBB: 10 mM Tris-HCl pH 7.5, 1 mM EDTA, 20% PEG 8000, 2.5 M NaCl, 0.025% Tween 20. Prepare this solution from sterile stocks and store at 4°C

1% agarose gel and 1XTAE buffer

Bridge adaptor,: Biotinylated bridge adaptor (+) (5' /5Phos/GCTGAGGGA/iBiotin-dT/C) and Bridge adaptor (-) (5' /5Phos/CCTCAGCT). The biotin hapten is connected to the oligonucleotide by a 6-carbon atom chain.

☒ HEPES Sodium salt **Merck MilliporeSigma (Sigma-Aldrich) Catalog #H7006**

☒ Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T8787-50ML**

☒ Sodium deoxycholate **Catalog #D6750**

☒ Tween 20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P9416-50ML**

☒ Agencourt AMPure XP **Beckman Coulter Catalog #A63880**

☒ DNase I, RNase free **Thermo Fisher Scientific Catalog #EN0525** supplied with **MnCl₂**

☒ NEB 10X CutSmart Buffer **New England Biolabs Catalog #B7204S**

☒ SDS, 10% Solution **Life Technologies Catalog #AM9822**

☒ Proteinase K (2 ml) **Qiagen Catalog #19131** 20 mg/mL

☒ RNase A **Qiagen Catalog #19101** 100 mg/mL

☒ Qubit® dsDNA HS Assay Kit **Thermo Fisher Scientific Catalog #Q32854**

☒ HS Genomic DNA Assay 75 - 20000 bp **Agilent Technologies Catalog #DNF-488-0500**

☒ HS NGS Fragment Assay 1-6000 bp 500 reactions **Agilent Technologies Catalog #DNF-474-0500**

☒ 1 Kb Plus DNA Ladder **Invitrogen - Thermo Fisher Catalog #10787018**

☒ Lambda DNA **Thermo Fisher Catalog #SD0011**

☒ SYBR SAFE DNA stain **Invitrogen - Thermo Fisher Catalog #S33102**

☒ NEBNext Ultra II End Repair/dA-Tailing Module - 24 rxns **New England Biolabs Catalog #E7546S**

☒ NEBNext FFPE DNA Repair Mix - 24 rxns **New England Biolabs Catalog #M6630S**

☒ T4 DNA Ligase, LC (1 U/μL) **Thermo Fisher Catalog #EL0016**

☒ NEBNext® Multiplex Oligos for Illumina® (Index Primers Set 3) **New England Biolabs Catalog #E7710L**

☒ NEBNext Ultra II Q5 Master Mix - 250 rxns **New England Biolabs Catalog #M0544L**

☒ Dynabeads™ M-270 Streptavidin **Thermo Fisher Scientific Catalog #65305**

☒ Polyethyleneglycol 8000 50% w/v **Jena Bioscience Catalog #CSS-256**

Equipment

- Bench top centrifuge with swing bucket rotor, refrigerated
- Fume hood
- Disposable 1 µL inoculation loops, sterile
- Ice bucket
- Orbital shaker
- Miracloth square (10 × 10 cm) with small funnel
- 50 mL Falcon tubes, sterile with rack
- 100 and 40 µm cell strainers
- Dounce homogenizer** with B pestle, 7 mL size
- 1.5 mL Eppendorf tubes, screw capped tubes and microcentrifuge
- 0.2 mL PCR tubes
- Magnet for 1.5 mL tubes
- 0.5 mL microcentrifuge tubes compatible with Qubit fluorimeter
- FA-WASTE container, β-ME-WASTE container and MnCl₂-WASTE container
- 1 mL Wide Bore tips
- 200 µL Wide Bore tips
- Mini gel box and powerpack
- Water bath
- Timer

Equipment	
ThermoMixer® C	NAME
Eppendorf	BRAND
Catalog No. 2231000680	SKU
https://online-shop.eppendorf.us/US-en/Temperature-Control-and-Mixing-44518/Instruments-44519/Eppendorf-ThermoMixerC-PF-19703.html	LIN K
	

Equipment

Gel Doc XR+ Gel Documentation System

NAME

Gel Documentation System

TYPE

Bio-rad Laboratories

BRAND

1708195

SKU

<https://www.bio-rad.com/en-us/product/gel-doc-xr-gel-documentation-system?ID=O494WJE8Z>^{LINK}

Equipment

Fragment Analyzer

NAME

capillary based nucleic acid fragment size separation

TYPE

Agilent

BRAND

M5311AA

SKU

<https://www.agilent.com/en/product/automated-electrophoresis/fragment-analyzer-systems/fragment-analyzer-systems/5300-fragment-analyzer-system-365721>^{LINK}

Troubleshooting

Safety warnings

- ! Prepare waste containers for the formaldehyde and manganese solutions labelled with the proper warning signs:

FA-WASTE container, β -ME-WASTE container and $MnCl_2$ -WASTE container

Before start

Characterize your sample

The protocol is developed for plant leaf tissues that will yield at least **5 µg nuclear gDNA per gram of leaf sample**. It is crucial that you characterize your biological sample before starting. This will give you a good understanding of the expected yields along the process and make decisions if you fall below the minimum amount of DNA to safely continue to the next stage.

Adjust the amount of starting leaf tissue to obtain the desired amount of ngDNA and make sure you have at least 6 tubes of nuclei before starting because it will take at least two tries before you get the optimal digestion conditions. Each tube should contain at least 2 µg of ngDNA.

How to anneal the bridge adaptor

- Dissolve the oligos in TE pH 7.5, at 1 mM. Store at -20°C after use
- Heat up a water bath to 65°C, containing at least 3 L of water
- Prepare the following mix in a 1.5 mL screw capped tube:

A	Volume, µL
TE pH 7.5	82
25 mM MgCl ₂	8
1 mM Biotinylated bridge adaptor (+)	5
1 mM Bridge adaptor (-)	5
Total volume	100


- Vortex briefly and do a quick spin
- Incubate the adaptor mix at 65°C for 5 min
- Transfer ~1.5 L from the water bath to a shallow plastic container (2 L capacity) on the lab bench
- Transfer the adaptor mix tube to the container, and place a thermometer
- Let it cool down to ~ 25°C, it should take about 2.5 h, but you can leave it overnight
- When ready, do a quick spin and store at -20°C
- This solution is enough for five bridge ligation reactions

Crosslinking Solution

0.5 mL 1X PBS + 27 µL

⊗ Formaldehyde solution Merck MilliporeSigma (Sigma-Aldrich) Catalog #F8775-25ML , freshly prepared

Streptavidin magnetic bead binding buffers

The following buffers can be prepared ahead of time and stored at  4 °C for at least 6 months

	A	1X B/W + T20	1X B/W	2X PEGBB
	ddH ₂ O	39.25 mL	39.75 mL	3 mL
	100X TE pH 7.5	250 µL	250µL	200 µL
	5 M NaCl	10 mL	10 mL	8 mL
	50% w/v PEG 8000	--	--	8 mL
	2.5% v/v Tween 20	500 µL	--	800 µL
	Total volume	50 mL	50 mL	20 mL

Abbreviations:

B/W + T20 = Bind and Wash plus Tween 20 buffer (0.5X TE pH 7.5, 1 M NaCl, 0.025% v/v Tween 20)

B/W = Bind and Wash Buffer (0.5X TE pH 7.5, 1 M NaCl)

PEGBB = Polyethylene glycol 8000 Binding Buffer (1X TE pH 7.5, 2 M NaCl, 20% v/v PEG 8000, 0.1% v/v Tween 20)





Nuclei isolation and integrity check


1




Safety information

Work in the fume hood

Add  20 mL sorbitol wash buffer to the ground tissue by gently dislodging it with an inoculation loop. Stir gently and drag the lump up against the tube wall until it is resuspended

2 Centrifuge  3500 rpm, 10°C, 00:05:00 . Pour off the supernatant into designated β -ME-WASTE container. The supernatant will have the consistency of light syrup and it might be light green with fine sediment

3  go to step #1 at least another 2 times. It can take up to 4 rounds of washing to obtain a supernatant that pours off like water and is clear. Place the tube on ice

4 Add  20 mL NEB complete buffer and gently resuspend the ground tissue by tapping the bottom of the tube and mixing by inversion. Add another  20 mL NEB complete buffer and mix by inversion. Place the tube horizontally over an ice bucket, and on the orbital  undetermined, Room temperature , 00:10:00 Slow speed

5 Assemble the funnel over a 50 mL Falcon tube and place the Miracloth. Pre-wet the miracloth with some NEB complete buffer


6 Filtrate the extract through one sheet of Miracloth

7 Filtrate the extract through a **100 μ m cell strainer** assembled over a 50 mL Falcon tube


8 Filtrate the extract through a **40 μ m cell strainer** assembled over a 50 mL Falcon tube. Lift the cell strainer slightly and press the mesh against the tube inner wall to speed up the filtration

**Note**


The flow through the 40 μ m cell strainer can be very slow. If it clogs, continue with a new 40 μ m cell strainer or stir the liquid very gently with a sterile 1 μ L inoculation loop

9 Collect the nuclei by centrifugation at  1800 x g, 10°C, 00:15:00

9.1 Discard the supernatant in the designated β -ME-WASTE container


10 Add  25 mL NEB β -ME and resuspend the pellet gently by tapping the bottom of the tube and mixing by inversion


11  go to step #9  go to step #9.1

12 Add  1 mL NEB β -ME and detach the pellet gently by tapping the bottom of the tube. Resuspend the nuclei using a 1 mL pipette set at 500 μ L and wide bore tips by drawing the liquid in and out very slowly. Avoid introducing air bubbles. It is very important the nuclei suspension is **homogeneous**

**Note**

A Dounce homogenizer with a B pestle can be used for this task

12.1 Transfer  50 μ L resuspended nuclei to a 1.5 mL tube to check the DNA integrity (see below)

















13 Split the nuclei suspension in **eight** 1.5 mL screw capped tubes, spin down 1 min at maximum speed, discard the buffer and the nuclei pellet tubes at  -80 °C



Note



- *To defrost a nuclei sample, place it on ice and flick the bottom of the tube gently until the pellet is resuspended*
- *If the cell density of the leaf sample is low, split the nuclei suspension in 5 or 6 tubes instead*


13.1 It is recommended to check the **DNA integrity** before preparing the library


1. Add  50 μ L 1X CutSmart Buffer + 1% SDS (prewarmed @ 55°C) to the nuclei
  go to step #12.1 vortex briefly
2. Add  1 μ L Proteinase K and vortex 2 seconds
3. Place the tube in the thermomixer programmed as follows:
55°C 15 min 1250 rpm \rightarrow 68°C 45 min 1250 rpm
1. Spin down the tube briefly and add  100 μ L AMPure XP beads mix and let stand
 00:05:00 room temperature
2. Place  00:02:00 magnet or until the solution is clear. Discard the supernatant
3. While the tube is on the magnet, add  200 μ L 80% ethanol over the beads and remove it  00:00:30 exactly . Repeat this step one more time
4. Quick spin the tube, place it on the magnet and remove all ethanol
5. Remove tube from magnet and add  50 μ L TE buffer . Tap the bottom of the tube to resuspend the beads and let the DNA elute at least
 00:10:00 room temperature
6. Place the tube on the magnet and transfer the eluted DNA to a new tube
7. Quantify  1 μ L eluted DNA with the HS dsDNA Qubit kit
8. Calculate the amount of nuclear genomic DNA in 50 μ L and then extrapolate to
  go to step #12 1 mL of resuspended nuclei
9. Calculate how much ngDNA is per tube (~118 μ L if split in 8 tubes, or ~158 μ L if split in 6 tubes). Ideally you will have ≥ 3 μ g of ngDNA per tube
10. Analyze  200 ng nuclear genomic DNA in a 1% agarose gel against the 1 kb+ ladder and  200 ng Lambda DNA (48.5 kbp) . The extracted DNA should be >20 kbp. If the DNA is not intact, **do not proceed** with crosslinking and DNase I test. Prepare a new nuclei prep and QC the sample again

Chromatin Crosslinking



14 Add  0.5 mL 1X PBS to one of the defrosted nuclei tubes  [go to step #13](#) . Tap gently the bottom of the tube, or use an inoculation loop to dislodge the nuclei from the tube walls. The nuclei should be fully resuspended


15 Add  0.5 mL crosslinking solution to the nuclei pellet, and mix by inversion

16 Incubate at  00:10:00 room temperature with gentle mixing using the orbital shaker by placing the tube horizontally

10m


**Note**

If you can't find an orbital shaker, a rotisserie oven set at room temperature or a hula mixer are good alternatives


17 Centrifuge at  3500 rpm, 00:05:00 , at room temperature and discard the supernatant in designated FA-WASTE container

Note

Some plant nuclei might not form a solid pellet at this buffer/centrifugal speed/time combination. You can increase the speed to 6000 or 13000 rpm instead.

18 Quench the crosslinking reaction by adding  1 mL 1X Wash buffer and resuspend the pellet gently by tapping the bottom of the tube, or by using an inoculation loop, and mix by inversion





19 Centrifuge at  3500 rpm, 00:05:00 , (6000 or 13000 rpm) and discard the supernatant in designated FA-WASTE container



Safety information**Move back to your lab bench**



Lysate characterization



- 20 Resuspend nuclei in  500 μ L FA buffer and incubate on the thermal mixer set at  1250 rpm, 37°C, 00:15:00


Note

Optional: RNA can be removed after lysis by adding  5 μ L RNase A (100 mg/uL) and incubating  Room temperature 5 min



- 21 Take  50 μ L lysate and  go to step #13.1 DNA QC to extract the lysate DNA and run 100 or 200 ng in a 1% agarose gel. The DNA should be > 20 kbp. **Do not proceed** if the DNA is degraded.

Keep the rest of the extracted lysate DNA at 4°C for quality check in the Fragment Analyzer at a later stage, together with the DNase I titration experiment (see below).

Note



After the incubation, the lysate can be stored at  -80 °C if you can't continue with the next section

Chromatin Capture

- 22 Calculate the volume of lysate to have **2 μ g** of DNA  go to step #21 and add 2 volumes of AMPure XP beads, mix and let it stand  00:05:00 at room temperature


Note


The chromatin lysate in the presence of magnetic beads becomes messy. The beads might stick to the inside and outside of the pipette tip and won't come off easily by pipetting in and out. Instead, dispense the beads above the lysate and mix the beads/lysate by tapping the bottom of the tube to avoid losing captured chromatin every time you add any liquid to the tube.

- 23
- Place the tube on the magnet for at least  00:02:00 or until the solution is clear.
 - Discard the supernatant, and wash the beads with  1 mL 1X Wash buffer .





4m


- Place the tube on the magnet again for at least  00:02:00 or until the solution becomes clear and discard it

24  [go to step #23](#) at least one more time, but it might take 4 or more rounds of washing until the solution becomes clear when the beads are on the magnet. The total number of washes depend on the plant species.

Note


If you can't proceed with the rest of the steps, add  200 µL 1X Wash buffer to the beads and store at  4 °C overnight

DNase I Digestion Dilution Series

25 **This section explains how to determine the optimal condition to digested chromatin. You will use one of the the nuclei tubes  [go to step #13](#) .**
Once you have found the optimal condition, you will apply to them to a new nuclei tube (section DNase I Digestion and onwards)

Prepare the following solutions:

	A	DNase I solution	Beads solution
	Sterile ddH2O	450 µL	400 µL
	10X DNase I reaction buffer	50 µL	50 µL
	100 mM MnCl2	--	50 µL
	Total volume	500 µL	500 µL

25.1 Set the thermomixer at  1250 rpm, 37°C ∞ min

- 26
- Mix the **DNase I 1 U/µL stock** by inversion 20 times, quick spin
 - Mix by inversion again 20 times, and quick spin

26.1

- Transfer 20 µL of DNase I solution to a new tube and add 5 µL **DNase I 1 U/µL stock**



- Mix by inversion 20 times and quick spin
- The **DNase I** concentration of this dilution is **0.2 U/μL**
- Place tube on ice

26.2 ▪ Prepare the dilution series as follows:

A	DNase I solution	DNase I 0.2 U/μL	U/μL
Dilution A	94 μL	6 μL	0.012
Dilution B	95 μL	5 μL	0.01
Dilution C	96 μL	4 μL	0.008
Dilution D	97 μL	3 μL	0.006

- Mix the dilution by tapping the bottom of the tube 20 times and quick spin
- Place on ice

27 After removing the last 1X Wash buffer solution from

➡ go to step #24 captured chromatin add 🧴 500 μL Beads solution

➡ go to step #25 and resuspend by vortexing briefly

28 Prepare one tube labelled **T = 0 h** and four tubes as follows:

A	Dilution A	Dilution B	Dilution C	Dilution D
T = 1 h	A-1 h	B-1 h	C-1 h	D-1 h

- 29
- Transfer 🧴 50 μL captured chromatin from ➡ go to step #27 to the tube labelled **T = 0 h**
 - Transfer 🧴 100 μL captured chromatin to each tube **A-1 h to D-1 h**
 - Place the 5 tubes on ice

- 30
- Start the digestion by adding **10 μL** of its corresponding dilution to tubes **A-1 h to D-1 h** and flick the bottom of the tube to mix the beads
 - Transfer the 5 tubes to the thermomixer
 - 🕒 01:00:00 TIMER

- 31
- Stop the reaction by adding 🧴 25 μL 0.5 M EDTA pH 8.0 to tubes **A-1 h to D-1 h** and vortex 3 s
 - Add 🧴 12.5 μL 0.5 M EDTA pH 8.0 to tube **T = 0 h** and vortex 3 s

- 32 Add 🧴 200 μL 1X Wash buffer to each tube, vortex briefly and place them on the magnet until the solution is clear

1h





- 32.1 Discard the supernatant in the MnCl_2 -WASTE container and go to step #32 to repeat the wash step
- 33 After removing the supernatant, add 200 μL 1X Wash buffer and 2 μL 0.5 M EDTA pH 8 to each tube, vortex briefly and store at 4 °C overnight

STEP CASE

Time series option

From 83 to 84 steps

1. Prepare **100 μL DNase I solution**
2. Mix the DNase I stock solution (1 U/ μL) by inversion and do a quick spin to collect the solution at the bottom of the tube
3. Add **1 μL DNase I 1 U/ μL** to 100 μL DNase I solution and mix by tapping the bottom of the tube and do a quick spin. **The DNase I is now diluted at 0.01 U/ μL . Place the tube on ice**
4. After removing the last 1X Wash buffer solution from the captured chromatin add **500 μL Beads solution** and resuspend by vortexing briefly
5. Prepare the tubes for the digestion time course. Label 9 tubes at 10 min intervals from T= 10 min to T = 1 h 30 min and add 12.5 μL 0.5 M EDTA pH 8.0 to each tube. The EDTA will stop the DNase I digestion.
6. Remove **50 μL captured chromatin** and transfer to a tube labelled ZERO and place it in the thermomixer
7. Add **2 μL DNase I 0.01 U/ μL** to the **450 μL remaining captured chromatin**, vortex gently and place the tube in the thermomixer.
8. START THE TIMER
9. Sample **50 μL captured chromatin digestion** every 10 min. Transfer the bead sample to the tube labelled with the time point. Vortex briefly and leave at room temperature while you finish sampling the rest of the beads. **Note:** *The bead lysate mixture is messy. Avoid dipping the pipette tip too far into the solution so you don't carry too much mixture outside the tube.*
10. Once the time course sampling is finished, add **200 μL 1X Wash buffer** to each tube, mix and place them on the magnet until the solution is clear
11. Discard the supernatant in the MnCl_2 -WASTE container and repeat the wash step one more time
12. After removing the supernatant, add **200 μL 1X Wash buffer plus 2 μL 0.5 M EDTA pH 8.0** mix and store at 4 °C overnight

Note: *If the digestion is incomplete, repeat the experiment but add **6 μL of DNase I at 0.01 U/ μL** in step 7 and continue with the rest of the time course.*

Chromatin Quality Check

- 34
- Incubate 1 mL 1X CutSmart Buffer NEB + 1% SDS for 55 °C 10 min before use
 - Vortex and quick spin before opening the tube
- 34.1 Set the thermomixer program as follows:
55°C 15 min 1250 rpm → 68°C 45 min 1250 rpm



- 35
- Equilibrate the tubes [go to step #33](#) at room temperature ~ 5 minutes
 - Place the tubes on the magnet for at least 00:02:00 until the solution becomes clear
 - Discard the supernatant and place the tubes on a rack

36 Add 100 μ L prewarmed 1X CutSmart Buffer NEB + 1% SDS to each tube

36.1 Add 1 μ L Proteinase K (20 mg/mL) to each tube, vortex briefly and quick spin

Note

Most samples will be fully digested with 20 μ g of Proteinase K under these conditions, but you could increase it to 40 or 50 μ g if you have experienced problems digesting your sample with this enzyme in other protocols (total genomic DNA extraction, for example).

37 Place the tubes in the thermomixer [go to step #34.1](#) and start the program

38 When the Proteinase K digestion is finished, quickly spin down the tubes and place them on the magnet for at least 00:02:00

39 Transfer the **supernatant containing the DIGESTED CHROMATIN to a new tube** and discard the beads







- 40
- Add 180 μ L AMPure XP beads, mix by flicking the bottom of the tube until the solution is homogeneous
 - Let it stand 00:05:00 at room temperature

40.1 Capture the beads with the magnet for at least 00:02:00 and discard the supernatant in the MnCl_2 -WASTE container


- 40.2
- Without removing the tube from the magnet add 200 μ L 80% ethanol over the beads and wait for **exactly** 00:00:30 before removing the ethanol
 - Repeat the ethanol wash one more time


40.3 Spin down the tubes 00:00:02 and place them on the magnet again and remove **all** the liquid

- 41
- Add  50 µL TE buffer directly onto the beads, close the tube and remove from the magnet. Mix the beads by gently flicking the bottom of the tube
 - Do a quick spin (~ 1 second) to collect the beads at the bottom
 - Let the DNA elute at least  00:10:00 at room temperature



- 41.1 Place the tube on the magnet for  00:02:00 until all the beads are against the magnet and transfer the eluted DNA to a new tube. Keep at  Room temperature

Note

If you can't proceed with the DNA quantification step, store the tubes at  4 °C overnight

- 42 Quantify  1 µL eluted DNA with the High Sensitivity dsDNA Qubit kit and estimate the total amount of DNA for each time point (see table below)

- 43 Prepare the samples for the Fragment Analyzer run. Use the HS NGS kit 1-6000 bp:

- If the concentration is **> 5 ng/µL**, dilute it to 2.5 ng/µL with TE buffer and use  2 µL of 2.5 ng/uL dilution
- If the concentration is **< 2.5 ng/µL** use  2 µL undiluted .
- Open the Fragment Analyzer outfile with ProSize 3.0 and perform the Smear Analysis by selecting the following fragment ranges, and fill out the table below:

44

	A	B	C	D	E	F	G	H	I	J
		% Sme ar (bp)	% Sme ar (bp)	% Sme ar (bp)	% Smea r (bp)	Qu bit	Qu bit	Qu bit	Sample (1 nuclei tube)	Sample (1 nuclei tube)
		100- 250 0	300- 200 0	100- 300	2000 -250 0	ng/ µL	Vo l µL	Yie ld ng	Vol µL	Yield ng
		Opti mal	Middl e	Low end	High end					
	Nucl ei samp le						50		equals 1 mL by number of nuclei tubes (step 13)	Aim for 2000

	A	B	C	D	E	F	G	H	I	J
	Lysate						50		500 (step 20)	Expect 1900-2000
	T = 0 h						50			
	Dilution A						50			
	Dilution B						50			
	Dilution C						50			
	Dilution D						50			

Note

The nuclei, lysate and T = 0 h samples will run over the upper marker (6000 bp), this is fine. It is more important that the resolution is higher between 1 and ~ 3000 bp. But if you prefer, use the HS Genomic DNA kit instead. However, the ProSize 3.0 software can't combine data from two different methods in the "Project" option. The amplified Hi-C library (see below) is best analyzed with the HS NGS method because the resolution between 1 and 500 bp is higher. If you want to compare this data with the DNase I digest and later the intramolecular ligated chromatin sample, use the HS NGS method.

45

- The amount of DNA in the nuclei sample allows you to estimate the DNA yield per tube and per gram of tissue
- The lysate should produce intact DNA, with same profile as your nuclei QC sample described in [go to step #13 Nuclei prep](#). The lysate QC step allows you calculate how much chromatin to expect per lysis event. The lysate and T = 0 h are internal controls to monitor if DNA degradation occurred during sample manipulation
- The goal of the DNase I dilution series is to determine the optimal amount of DNase I where most of the DNA is within the 100-2500 bp range

≥ 50% of the input DNA is within 100-2500 bp

- If you have at 0.9-1 µg of optimally digested chromatin you can now take another nuclei tube and process it in the same way (Section DNase I onwards) to produce either a Nanopore or an illumina library




- If you have ~ 150 ng of optimally digested chromatin you can only prepare an illumina library under these condition
- If you have 80-100 ng of optimally digested chromatin you are at the lower limit for preparing an illumina library, but it is till possible to produce one

< 50% of the input DNA is within 100-2500 bp

- If most of the smear is ≤ 300 bp for all dilutions (or time points - see step case option), the sample was **over digested**. Repeat the DNase I dilution series but incubate the digestion for 30 min only. If you did a time series
[⇒ go to step #33 step case](#) use a 0.005 U/ μ L dilution of DNase I
- If most of the smear is ≥ 2500 bp, the sample was **under digested**. Before repeating or modifying the chromatin prep, check that the DNase I stock is working correctly. Perform a control digestion of 1 μ g of Lambda DNA, 20 min only, in 1X DNase I buffer + 10 mM $MnCl_2$, 0.02 U DNase I, in a total volume of 20 μ L, 37°C. Stop the reaction with 1 μ L 0.5 M EDTA pH, add 4 μ L 6X gel loading buffer and load it **all** in a 1% agarose gel and include an aliquot of undigested Lambda DNA as control (1X TAE, 140 V, 30 min, stain the gel with SybrSafe after the run or add it directly on the gel). If the test shows no digestion, buy a new DNase I stock (with $MnCl_2$)
- If the under digested result is accompanied by a high DNA yield (**>5 μ g**), QC the lysate again and make sure you only use **2 μ g** for the chromatin capture step

DNase I Digestion

- 46 This section describes how to process one nuclei tube from crosslinking to the DNase I digestion step once the optimal reaction conditions have been established. 

You will collect samples of lysate, T = 0 h and optimal digestion time for QC later. This is a precaution to ensure the same results observed earlier are reproduced

- 46.1
- Select the time point that contains **$\geq 50\%$ of the DNA within 100-2500 bp** to prepare your library
 - Take one tube [⇒ go to step #13 Nuclei prep](#) from the -80°C freezer and perform the steps from [⇒ go to step #14 Chromatin crosslinking](#) up to [⇒ go to step #24 Chromatin capture](#)
- 46.2 Carry out the **DNase I digestion** with two time points only: **T = 0 h** and **Optimal QC** as follows:



- Transfer 50 µL captured chromatin resuspended in 500 µL Beads solution to a tube labelled **T = 0 h**
- Start the reaction by adding the [go to step #26.2](#) Optimal DNase I dilution to the remaining 450 µL captured chromatin and place both tubes in the thermomixer

For example: if the optimal digestion was observed with dilution B

	A	Captured chromatin	Dilution B @ 0.01 U/µL	U/captured chromatin (µL)
	Dilution series	100 µL	10 µL	0.001
	Library prep	450 µL	45 µL	0.001

- When the time is up, transfer 50 µL DNase I digested captured chromatin to a tube labelled **Optimal QC** and add 12.5 µL 0.5 M EDTA pH 8.0 , mix and leave it on a rack at room temperature
- Add 12.5 µL 0.5 M EDTA pH 8.0 to the **T = 0 h**, mix and leave it on the rack
- Add 100 µL 0.5 mM EDTA pH 8.0 to the remaining **400 µL** of digested chromatin, mix and leave it on the rack

46.3




- Add 800 µL 1X Wash Buffer to each tube, mix and place on the magnet for at least 00:02:00 until the solution is clear
- Discard all the solution in the MnCl₂-WASTE container
- Remove the tubes from the magnet and add 1 mL 1X Wash buffer to each tube, mix and place on the magnet again
- When the solution is clear transfer it to the MnCl₂-WASTE container

2m

46.4

- Extract the DNA from **T = 0 h** and **Optimal QC** tubes as described before
 [go to step #34](#) Chromatin Quality Check

These samples, together with the leftover lysate ([go to step #21](#)) and the intra molecular ligated chromatin will be quantified and loaded in the Fragment Analyzer later (see Chromatin Reverse Crosslinking Section and Quality Check section below)



- 46.5
- Add  200 µL 1X Wash Buffer and  2 µL 0.5 M EDTA pH 8.0 to the **Digested chromatin** tube and store at  4 °C overnight
 - This is the sample to be used for preparing your library (see section **End Repair**)

End Repair

47

- From the Optimal QC sample determine how much DNA is digested and captured in the AMPure XP beads
- The NEBNext Ultra II End repair can only repair up to **1 µg** of DNA. If there is more than 1 µg of digested chromatin needs to be repaired, up-scale the End Repair according to the table shown below
- The **NEBNext FFPE DNA repair reaction** is only needed for Nanopore library preparation and can only repair up to 1 µg of DNA. If you are preparing an Illumina library, omit this reagent
- Make sure the buffers are fully defrosted and dissolved. Vortex and quick spin all the reagents

48


Place the **digested captured chromatin** from  [go to step #46.5](#) on the magnet and discard the clear supernatant. Wash the beads twice with  200 µL 1X Wash buffer and discard all the liquid


49

Add the following reagents of the desired library to be prepared to the washed beads in the order shown below. Mix the contents before and after adding the enzymes:

Reagent	Nanopore library	Illumina library
	Volume µL	Volume µL
Total volume	60	60
Deionized sterile water	48	50
NEBNext End Repair buffer	3.5	7
NEBNext FFPE DNA Repair buffer	3.5	
End repair enzyme mix	3	3
FFPE enzyme mix	2	



49.1

Place the tube in the thermal mixer set at  1250 rpm, 20°C, 01:00:00


49.2 Quick spin the tube and place it on the magnet for ~  00:02:00 magnet and discard all the liquid

2m

49.3

- Remove the tube from the magnet and add  1 mL 1X Wash buffer and gently resuspend the beads by tapping with a pen
- Quick spin the tube and place it for at least  00:02:00 magnet and discard all the liquid

2m

49.4 Repeat  [go to step #49.3](#) one more time and proceed to the **Bridge Ligation** step immediately


Bridge Ligation

50 Defrost the ligase buffers and vortex until no solid precipitate (DTT) and no syrupy pellet (polyethylene glycol) are visible. Keep the buffers on ice




51 Add the following reagents to the beads in the **order specified**:

Reagent	Volume μL
Deionized sterile water	50
Annealed Biotinylated Bridge 50 μM	20
Mix beads gently	
10X T4 DNA ligase buffer	10
50% PEG 4000	10
Mix beads gently	
T4 DNA ligase 1 U/ μL	10
Mix beads gently	
Final volume	100

51.1 Place the tube in the thermomixer set at  1250 rpm, 22°C, 01:00:00


51.2 Quick spin the tube and place it for at least  00:02:00 magnet and discard all the liquid

2m


- 51.3
- Remove the tube from the magnet and add  500 µL 1X Wash buffer and gently resuspend the beads by tapping with a pen
 - Quick spin the tube and place it for at least  00:02:00 magnet and discard all the liquid
- 51.4
- Repeat  [go to step #51.3](#) one more time making sure all the liquid is discarded
 - Place the tube on a rack and continue with **Proximity Ligation** section

2m

Proximity Ligation

- 52 While performing the Bridge ligation, defrost the ligase buffers and vortex it until no solid precipitate (DTT) and no syrupy pellets are visible. Keep the buffers on ice
- 53 Add the following reagents to the beads  [go to step #51.4](#) in the **order specified**:

A	B
Reagent	Volume µL
Deionized sterile water	390
10X T4 DNA ligase buffer	50
50% PEG 4000	50
Mix beads gently	
T4 DNA ligase 1 U/µL	10
Mix beads gently	
Final volume	500

- 53.1 Place the Proximity Ligation reaction tube in the thermomixer set at  1250 rpm, 22°C, 16:00:00



Chromatin Reverse Crosslinking and Quality Check

- 54
- Program the thermomixer as follows: **55°C 15 min 1250 rpm** → **68°C 45 min 1250 rpm**
 - Prewarm the thermomixer at **55°C** for ~ 15 min without shaking




- 55
- Incubate 1 mL 1X CutSmart (NEB) + 1% SDS for 00:10:00 at 55 °C before use
 - Vortex and quick spin before opening the tube
- 55.1 Quick spin the [go to step #53.1](#) Proximity Ligation reaction tube , place it for at least 00:02:00 magnet and discard all the liquid
- 55.2 Add 100 µL 1X CutSmart buffer (NEB) + 1% SDS and vortex briefly
- 55.3 Add 1 µL Proteinase K , vortex briefly and quick spin the tube (see note [go to step #36.1](#))
- 55.4 Place the tube in the thermomixer and start the program [go to step #54](#)
- 56 Quick spin the tube and **TRANSFER THE LIQUID TO A NEW TUBE**



THIS IS YOUR INTRA-MOLECULAR LIGATED CHROMATIN (IMLC)

Structure of an intra molecular ligated chromatin contact



- 56.1 Add 50 μ L TE buffer pH 7.5 to the beads, vortex gently, place for 00:02:00 magnet and transfer the liquid to the rest of the IMLC from the previous step

The total volume is **150 μ L IMLC**. Discard the beads

Note

If you can't continue with the DNA purification, store the IMLC at -20 $^{\circ}$ C .

Defrost the tube on the bench and warm it up at 30 $^{\circ}$ C ~ 5 min until the solution is clear

- 57 Add 105 μ L AMPureXP beads to the **150 μ L IMLC** and vortex gently until the sample is homogeneous




Note

The bead to sample ratio at 0.7 will select fragments above 200 bp. See Left Side Selection graph [here](#)

- 57.1 Incubate 00:05:00 at room temperature
- 57.2 Place the tube on magnet and let it stand for at least 00:02:00 magnet or until all the beads are against the magnet and discard the liquid 2m
- 57.3
- While the tube is still on the magnet, add 200 µL 80% ethanol 30s
 - Remove it **exactly** after 00:00:30 ethanol wash
- 57.4 Repeat step [go to step #57.3](#) one more time
- 57.5
- Place the tube on the centrifuge and do a quick spin (~ 00:00:02) and place on magnet again Remove **all** traces of ethanol 2s
- 57.6
- Add 54 µL 10 mM Tris-HCl pH 8.0 directly onto the beads, vortex gently and do a quick spin 10m
 - Let the DNA elute at room temperature for **at least** 00:10:00 DNA elution
- 57.7
- Place the tube on the magnet for ~ 00:02:00 magnet until all the beads are collected 2m
 - Transfer the eluted DNA to a new tube
- 58
- Quantify 1 µL IMLC and also 1 µL of **T = 0 h** and **Optimal QC** samples
 - [go to step #46.4](#) with the HS dsDNA Qubit kit
 - Run 2 µL at 2.5 ng/uL each sample on the Fragment Analyzer (HS NGS kit 1-6000 bp)
- 58.1 Store the IMLC at -20 °C until ready to prepare the sequencing library

Nanopore library preparation

- 59 To prepare a Nanopore library you will need a least  1 µg IMLC . Check the [Nanopore Community page](#) to select the kit and best conditions for sequencing

Illumina library preparation


60 IMLC End Prep and Illumina adaptor ligation

The universal illumina adaptor and USER enzyme can be found in the NEBNext Multiplex Oligos for Illumina (Index Primers Set 3) kit

Set up the following reaction in a 200 µL PCR tube:

A	Volume µL
IMLC	51
NEBNext Ultra II End prep Buffer	7
NEBNext Ultra II End prep Enzyme	3
Total volume	61

Place the tube on a thermocycler programmed as follows: **20°C 30 min** → **65°C 30 min** → **12°C ∞**

- 61 Add  2.5 µL NEB Universal Illumina adaptor (15 µM) , vortex and quick spin

- 62 Add the following reagents to the End prep IMLC + illumina adaptor tube in the **order specified**:

Reagent	Volume µL
Deionized sterile water	11.5
10X T4 DNA ligase buffer	10
50% PEG 4000	10
Mix beads gently	
T4 DNA ligase 1 U/µL	5
Mix beads gently	
Final volume	100



- 62.1 Incubate at 20 °C 30 min . Or leave it in the refrigerator overnight
- 62.2 Add 3 µL USER enzyme mix , vortex, quick spin and incubate at 37 °C 15 min
- 62.3 Spin down the tube and transfer all the sample to a 1.5 mL tube
- 62.4
- Add 50 µL TE pH 7.5 to the PCR tube, vortex, quick spin and transfer all to the 1.5 mL tube
 - Total volume: **153 µL**
- 63 Add 120 µL AMPure XP beads and vortex gently until completely homogeneous
- 63.1
- Incubate Room temperature at least 5 min
 - Place the tube on magnet 00:02:00 or until the solution becomes clear and discard the solution
- 63.2
- Add 180 µL 80% ethanol to the tube while on the magnet
 - Discard all liquid after 00:00:30 exactly
 - Repeat this wash one more time
 - Quick spin the tube, place it on the magnet and remove **all** liquid
- 63.3
- Add 100 µL TE pH 7.5 to the tube, remove from the magnet, and resuspend the beads by gently vortexing the tube
 - Incubate Room temperature for at least 00:10:00
 - Return the tube to the magnet and when the solution becomes clear transfer it to a new 1.5 mL tube
- This is now your **IMLC-illumina library**
- 64 Quantify 1 µL IMLC-illumina library with the Qubit HS dsDNA kit

Note

If you can't continue with the Streptavidin magnetic bead capture or for long term storage, keep the IMLC-illumina library at -20 °C





- 65 The following steps describe how to prepare the streptavidin magnetic beads for binding an IMLC-illumina library at a concentration of 2.5 ng/ μ L or higher. 4m

If your IMLC-illumina library is less concentrated, resuspend the SAMB in a different volume. For example: The IMLC-illumina library is at 1.5 ng/ μ L. To have 100 ng for SAMB capture you need 66 μ L IMLC-illumina . Resuspend the prepared SAMB beads in

66 μ L 2X PEGBB instead of 40 μ L as described here:

Streptavidin magnetic bead (SAMB) preparation

1. Coat one 1.5 mL screw-capped tube by adding 1 mL 1X B/W + T20 , mix and remove all liquid
2. Add 1 mL 1X B/W to the tube
3. Vortex the SAMB until fully resuspended and transfer 20 μ L SAMB to the tube, mix by pipetting. This volume contains **200 μ g** of beads
4. Place the tube on the magnet for at least 00:02:00 and discard **half** of the liquid
5. Quick spin the tube and place it on the magnet
6. When the solution is clear, discard it **all**
7. Remove the tube from the magnet and add 20 μ L 1X B/W and mix by pipetting
8. Place the tube on the magnet for 00:02:00 and discard the liquid
9. Repeat this wash step two more times
10. Add 40 μ L 2X PEGBB and mix by pipetting




66 SAMB capture

36m

1. Calculate the volume required of **IMLC-illumina library** [go to step #64](#) to have **100 ng** and adjust it to a final volume of 40 μ L 1X TE pH 7.5 . Mix and quick spin the tube
2. Transfer the IMLC-illumina solution to the prepared **SAMB** [go to step #65](#) , vortex briefly and place in the thermomixer 1250 rpm, 25°C, 00:30:00
3. Quick spin the tube and place it on the magnet for at least 00:02:00 and discard the solution*
4. Remove the tube from the magnet and add 120 μ L 1X B/W + T20 and mix by pipetting
5. Place the tube on the magnet for at least 00:02:00 and discard the solution*

Note

** The clear solution from the unbound and the washes contain IMLC-illumina library, about 50% of the input material. You could discard it, but it is recommended to keep it in case you need to prepare another SAMB capture sample. See instructions on how to recover the unbound IMLC-illumina from these fractions in the last step-case of the Expected Results section*

6. Repeat steps 4 and 5 one more time
7. Remove the tube from the magnet and add  120 µL 1X B/W and mix by pipetting
8. Place the tube on the magnet for at least  00:02:00 and discard the solution*
9. Repeat steps 8 and 9 one more time
10. Resuspend the beads in  40 µL 10 mM Tris-HCl pH 8.0

The bead concentration is ~ **5 µg beads/µL**. Store at  4 °C


The sample is stable at this temperature for at least 2 weeks

67

Library amplification and size selection

Set up the following PCR master mix:

A	5 reactions, µL
Sterile ddH2O	100
NEBNext Ultra II Q5 master mix	125
Universal PCR primer 10 µM	12.5
NEBNext index primer 10 µM	12.5
IMLC-illumina-SAMB ~5 µg/µL	10
Total volume	250

- Vortex until the beads are mixed evenly
- Aliquot  50 µL PCR master mix per PCR tube and place them in the thermocycler



PCR profile:

98°C 30 sec → (98°C 10 sec - 65°C 75 sec) x 12 → 65°C 5 min → 12°C ∞









67.1 Spin down the tubes briefly and transfer all the solution to one 1.5 mL tube
This tube contains the amplified Hi-C library plus IMLC-illumina-SAMB


67.2 Wash the PCR tubes:

- Add  100 μ L 1X TE pH 7.5 to one of the PCR tubes and set the pipette volume to **150 μ L**
- Pipette the solution several times and transfer to the next tube. Repeat this step on all the PCR tubes
- Transfer the solution to the 1.5 mL tube  [go to step #67.1](#)

67.3 Recover the SAMB-illumina bound library:



4m

1. Place the 1.5 mL tube  [go to step #67.2](#) on the magnet for at least  00:02:00
2. Transfer the clear solution to a **new 1.5 mL tube** and set aside. **This is your amplified Hi-C library to be sequenced**
3. Add  50 μ L 1X TE pH 7.5 to the SAMB beads, remove it from the magnet, vortex gently and place it in the magnet again for at least  00:02:00
4. Transfer the solution to the 1.5 mL tube set aside earlier (step 67.3.2)
5. Add  10 μ L 10 mM Tris-HCl pH 8 to the SAMB-illumina bound library and store at  4 $^{\circ}$ C This sample can be used again to prepare more amplified library if needed, but use it within 1 week
6. Estimate the total volume of amplified library using a 1000 μ L pipette . Expect ~ 400 μ L (250 μ L PCRs + 100 μ L PCR tubes wash + 50 μ L beads wash)



- 67.4
- Add  0.8 - volumes of AMPure XP beads to the PCR reaction (step 67.3.2). For example: if 400 μ L were recovered, add 320 μ L of AMPure XP beads
 - Vortex briefly until the solution is homogeneous



Note

The amount of AMPure XP beads will remove unused primers only. Once the amplified library has been analyzed, you will need to do a double size selection (0.5X/0.3X) before sequencing to select for fragments with an average size of 670 bp. See instructions described in the last step case in Expected Results section below.


- 67.5
- Incubate at  Room temperature 5 min
 - Place the tube on the magnet for at least  00:02:00 or until the solution is clear
 - Discard the solution



- 67.6
- While still on the magnet, add  500 μ L 80% ethanol to the beads and remove it exactly after  00:00:30
 - Repeat this step one more time

- 67.7
- Quick spin the tube and place it on the magnet again
 - Remove **all** traces of ethanol
 - Remove the tube from the magnet and place it on a rack
 - Add  100 μ L 1X TE pH 7.5 and gently vortex the tube until the beads are resuspended
 - Incubate  Room temperature at least 10 min
 - Place the tube on the magnet and when the solution is clear, transfer it to a new 1.5 mL Eppendorf tube

This is the amplified Hi-C library cleaned with 0.8X AMPure XP beads

- 67.8
- Qubit  1 μ L amplified Hi-C library with Qubit HS dsDNA kit
 - Prepare 10 μ L of a 2.5 ng/ μ L dilution of the amplified library and run it in the Fragment Analyzer (HS NGS 1-6000 bp method)

- 68
- We have tested this protocol on several species and although not all of them showed ideal digestions and IMLC profiles (see step case below) the MiSeq QC run produced usable libraries.




Fig. 1 Fragment size distribution of the lysate, **DNase I digestion**, and **Intramolecular ligated chromatin (IMLC)**.

The samples were characterized with either the HS NGS method or the HS Genomic DNA method on the Fragment Analyzer. The upper marker of the NGS method is denoted with a yellow rhomboid. The sharp peak at ~ 600 bp shown in some samples could be RNA or aggregated DNA. The peak disappears after the reverse crosslinking of the IMLC.

The sea urchin and bat samples included in this section were grounded in liquid nitrogen, washed and resuspended in 1X PBS before crosslinking. After quenching the formaldehyde with 1X Wash buffer, the cells were filtered through a 200 μ m cell strainer and collected by centrifugation. The myna bird blood cells were defrosted in 1X TE pH



7.5, 100 mM EDTA as soon as they were taken out of the -80°C freezer, spun down, washed with 1X PBS twice before crosslinking and processed as the other animal samples.

The table below shows the expected yields along the protocol across 5 different plants.


		Diploid kiwifruit (<i>Actinidia chinensis</i>)	<i>Eucalyptus</i> sp.	Karaka (<i>Corynocarpus laevigatus</i>)	Grape (<i>Vitis vinifera</i>)	<i>Wahlenbergia albomarginata</i>
Nuclei prep stage	Unit					
Leaf sample	g	3	3	4	4	10
# Nuclei tubes	tube	8	8	8	8	6
Hi-C Library prep stage		Total per nuclei tube				
Lysate	ng DNA	2200	1960	2380	2285	5000
Intramolecular ligated chromatin (IMLC)	ng DNA	825	83	906	142	398
IMLC-illumina library	ng DNA	581	61	644	100	290
# Streptavidin-magnetic bead (SAMB) capture (100 ng IMLC-illumina each)	captures	5	1	6	1	3
# Library amplification reactions		100	20	120	20	60
Library amplification stage (2 µL SAMB-IMLC-illumina per reaction)		10 PCR				
Amount amplicon after 0.8X AMPure XP clean up	ng DNA	2060	206	756	970	3640
After 0.5X/0.3X size selection	ng DNA	664	108	308	432	1408
% amplicon to be sequenced		32%	52%	41%	45%	39%

It is crucial to use no more than 100 ng of IMLC-illumina sample for the streptavidin magnetic bead capture because not all of it will be captured. Under the conditions described in this method, about 50 ng of the initial 100 ng of IMLC-illumina library is bound to the streptavidin beads. It is possible that by using a longer linker connecting the biotin hapten to the Bridge (+) strand, the steric hindrance might be avoided.

The Hi-C library amplification double size selection might look a bit lumpy:




And the next table shows the metrics obtained with the MiSeq QC run:

 [Hi-C MiSeq comparisons.xlsx](#)

STEP CASE

Ideal DNase I and IMLC profiles 1 step

DNase I digestion and Intra molecular ligated chromatin



The sharp peak at ~650 bp observed on the digested chromatin might be an artifact produced when the crosslinked chromatin binds to the AMPure XP beads. We have observed this peak in two plants and on a bird blood library prep. The peak disappears after the intramolecular ligation and reverse crosslinking steps.

Acknowledgements

- 69 We would like to thank Ashley Jones (Australian National University, Canberra) for joining us in our first Hi-C marathon in 2022, and our colleagues for providing samples to test: Susan Thomson, Wendy Hall, Ed Morgan (Plant and Food Research, New Zealand), Anna Santure (University of Auckland, New Zealand), Annabel Whibley and Darrell Lizamore (Bragato Institute), Jessie Prebble (Landcare Research, New Zealand), Niel Gimmel and Joanne Gillum (University of Otago, New Zealand).

Protocol references

Literature

- Duan, Z. (2021). Targeted DNase Hi-C. Capturing Chromosome Conformation. B. Bodega and C. Lanzuolo. New York, NY, Humana.
- Gollosi, R., J. T. Sanders and R. P. McCord (2018). "Iteratively improving Hi-C experiments one step at a time." *Methods* 142: 47-58.
- Golov, A. K., S. V. Ulianov, A. V. Luzhin, E. P. Kalabusheva, O. L. Kantidze, I. M. Flyamer, S. V. Razin and A. A. Gavrillov (2020). "C-TALE, a new cost-effective method for targeted enrichment of Hi-C/3C-seq libraries." *Methods* 170: 48-60.
- Gridina, M., E. Mozheiko, E. Valeev, L. P. Nazarenko, M. E. Lopatkina, Z. G. Markova, M. I. Yablonskaya, V. Y. Voinova, N. V. Shilova, I. N. Lebedev and V. Fishman (2021). "A cookbook for DNase Hi-C." *Epigenetics & Chromatin* 14(1): 15.
- Hoffman, E. A., B. L. Frey, L. M. Smith and D. T. Auble (2015). "Formaldehyde crosslinking: a tool for the study of chromatin complexes." *The Journal of biological chemistry* 290(44): 26404-26411.
- Kadota, M., O. Nishimura, H. Miura, K. Tanaka, I. Hiratani and S. Kuraku (2020). "Multifaceted Hi-C benchmarking: what makes a difference in chromosome-scale genome scaffolding?" *GigaScience* 9(1).
- Kasem, S., N. Rice and R. Henry (2008). DNA extraction from plant tissue. *Plant Genotyping II: SNP Technology*. 2: 219-271.
- Lieberman-Aiden, E., N. L. van Berkum, L. Williams, M. Imakaev, T. Ragoczy, A. Telling, I. Amit, B. R. Lajoie, P. J. Sabo, M. O. Dorschner, R. Sandstrom, B. Bernstein, M. A. Bender, M. Groudine, A. Gnirke, J. Stamatoyannopoulos, L. A. Mirny, E. S. Lander and J. Dekker (2009). "Comprehensive Mapping of Long-Range Interactions Reveals Folding Principles of the Human Genome." *Science* 326(5950): 289.
- Ma, W., F. Ay, C. Lee, G. Gulsoy, X. Deng, S. Cook, J. Hesson, C. Cavanaugh, C. B. Ware, A. Krumm, J. Shendure, C. A. Blau, C. M. Disteche, W. S. Noble and Z. Duan (2014). "Fine-scale chromatin interaction maps reveal the cis-regulatory landscape of human lincRNA genes." *Nature Methods* 12: 71.
- Ma, W., F. Ay, C. Lee, G. Gulsoy, X. Deng, S. Cook, J. Hesson, C. Cavanaugh, C. B. Ware, A. Krumm, J. Shendure, C. A. Blau, C. M. Disteche, W. S. Noble and Z. Duan (2018). "Using DNase Hi-C techniques to map global and local three-dimensional genome architecture at high resolution." *Methods* 142: 59-73.
- Nagano, T., C. Várnai, S. Schoenfelder, B. M. Javierre, S. W. Wingett and P. Fraser (2015). "Comparison of Hi-C results using in-solution versus in-nucleus ligation." *Genome Biol* 16(1): 175.
- Niu, L., W. Shen, Y. Huang, N. He, Y. Zhang, J. Sun, J. Wan, D. Jiang, M. Yang, Y. C. Tse, L. Li and C. Hou (2019). "Amplification-free library preparation with SAFE Hi-C uses ligation products for deep sequencing to improve traditional Hi-C analysis." *Communications Biology* 2(1): 267.
- Padmarasu, S., A. Himmelbach, M. Mascher and N. Stein (2019). In Situ Hi-C for Plants: An Improved Method to Detect Long-Range Chromatin Interactions. *Plant Long Non-Coding RNAs: Methods and Protocols*. J. A. Chekanova and H.-L. V. Wang. New York, NY, Springer New York: 441-472.
- Ramani, V., D. A. Cusanovich, R. J. Hause, W. Ma, R. Qiu, X. Deng, C. A. Blau, C. M. Disteche, W. S. Noble, J. Shendure and Z. Duan (2016). "Mapping 3D genome architecture through in situ DNase Hi-C." *Nat Protoc* 11(11): 2104-2121.