



Mar 31, 2025

Enzyme-linked immunosorbent assay (ELISA)

DOI

dx.doi.org/10.17504/protocols.io.kxygx7eydl8j/v1

Karyna Tarasova¹

¹Veterinary Medicine University



Karyna Tarasova

Veterinary Medicine University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx7eydl8j/v1>

Protocol Citation: Karyna Tarasova 2025. Enzyme-linked immunosorbent assay (ELISA). **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygx7eydl8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 31, 2025



Last Modified: March 31, 2025

Protocol Integer ID: 125756

Keywords: immunosorbent assay, linked immunosorbent assay, using elisa quantitative sandwich kit, elisa quantitative sandwich kit, total protein concentration, supernatant, enzyme, assay, protein

Abstract

To analyze the expression of TNF- α , supernatants were collected at the end of each designated culture period, centrifuged and stored until subsequent enzyme-linked immunosorbent assay (ELISA). Total protein concentrations were determined by Bradford assay. The concentrations of TNF- α were measured using ELISA Quantitative Sandwich kits according to the manufacturer's instructions.

Materials

Bradford assay (5000006, Bio-Rad Laboratories, California, USA)

ELISA Quantitative Sandwich kits (MBS2701341, MyBiosource, San Diego, USA)

Troubleshooting



Safety warnings

! The ELISA Quantitative Sandwich kit (MBS2701341) has high sensitivity and excellent specificity for detection of TNFa.







Before start

The ELISA Quantitative Sandwich kit (MBS2701341) was acquired from MyBiosource (San Diego, USA).

Sample preparation

- 1 Sample preparation and total protein concentration determination.
- 1.1 Cell culture supernatant was collected.
- 1.2 The supernatant was centrifuged 2000g x 10min x 4°C. The samples were kept on ice. 
- 1.3 Total protein concentrations were determined by Bradford assay (5000006, Bio-Rad Laboratories, California, USA) by measuring the absorbance (595 nm) with a plate reader (Varioskan LUX, Thermo Scientific, USA). 

Enzyme-linked immunosorbent assay (ELISA)

- 2
- 2.1 The standard curve was prepared by serial diluting the stock standard solution in the diluent according to the manufacturer's instructions (MBS2701341, MyBiosource, San Diego, USA). 
- 2.2 The standards, blank and samples were added into the pre-coated plate and incubate for 1h at 37°C. 
- 2.3 The liquid of each well was removed without wash. 100µL of Detection Reagent A working solution was added to each well and incubated for 1 hour at 37°C. 
- 2.4 The wells were washed with 1× Wash Solution three times. The remaining liquid from all wells was completely removed by snapping the plate onto absorbent paper. 
- 2.5 100µL of Detection Reagent B working solution was added to each well and incubated for 30 minutes at 37°C. 
- 2.6 The wells were washed with 1× Wash Solution three times. The remaining liquid from all wells was completely removed by snapping the plate onto absorbent paper. 



- 2.7 90µL of Substrate Solution was added to each well and incubated for 20 minutes at 37°C.
- 2.8 Without removing or washing, 50µL of Stop Solution was added to each well.
- 2.9 The absorbance (450 nm) was measured immediately with a plate reader (Varioskan LUX, Thermo Scientific, USA). TNF-α levels were normalized to total protein concentrations.

