



Aug 19, 2020

Enzyme linked immunosorbent assay (ELISA) for determining the serum concentration of IL-17 in humans.

DOI

dx.doi.org/10.17504/protocols.io.bjyykpxw

Angel A Justiz-Vaillant¹

¹University of the West Indies St. Augustine

University of the West In...

angel.vaillant@sta.uwi.e...



Angel A Justiz-Vaillant

University of the West Indies St. Augustine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bjyykpxw>

Protocol Citation: Angel A Justiz-Vaillant 2020. Enzyme linked immunosorbent assay (ELISA) for determining the serum concentration of IL-17 in humans.. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjyykpxw>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: August 19, 2020

Last Modified: August 19, 2020

Protocol Integer ID: 40696

Keywords: interleukin17a, immunosorbent assay, cytokine, linked immunosorbent assay, inflammatory, immunology, journal of immunology, production of chemokine, il

Abstract

Interleukin17A is a pro-inflammatory cytokine. This protein is produced by a group of T helper cell known as T-helper17 lymphocytes in response to their stimulation with IL-23. After interaction with its receptor, IL-17 activates signalling pathways that stimulate the production of chemokines [1].

Reference

1. Cho ML, Kang JW, Moon YM, Nam HJ, Jhun JY, Heo SB, Jin HT, Min SY, Ju JH, Park KS, Cho YG, Yoon CH, Park SH, Sung YC, Kim HY (May 2006). Journal of Immunology. 176 (9): 5652–61.

Materials

MATERIALS

⊗ IL-17, Interleukin-17 IL-17A, human **Bio Basic Inc. Catalog #RC212-28.SIZE.1mg**

⊗ Set of one 96-well filter plate with 2 plate sealers **Merck Millipore (EMD Millipore) Catalog #MX-PLATE**

Troubleshooting



- 1 Ninety-six well ELISA plates are coated with monoclonal anti-human antibodies to interleukin-17(IL-17).
- 2 Patient serum samples are added to the plates.
- 3 The plate is incubate x 1.30 hour at RT.
- 4 The plate is washed 4 times with PBS-tween 20 buffer.
- 5 The wells are incubated with a biotin conjugated anti-human IL-17 for 1.30 hour at RT.
- 6 The plates are washed again as above.
- 7 To the plate a peroxidase-labeled streptavidin conjugate is added and incubated for 1 hour at RT.
- 8 After a further washing procedure a substrate solution reactive is added and allowed to produced a colored reaction in positive controls.
- 9 The level of IL-17 in the sample is proportional to the colored product developed.
- 10 The addition of acid stopped the reaction.
- 11 The absorbance is measured at 450 nm.
- 12 The IL-17 concentration can be calculated by generating an standard curve.