

Aug 13, 2020

## Enzyme linked immunosorbent assay for investigating the binding of Protein-L to diverse immunoglobulins

DOI

[dx.doi.org/10.17504/protocols.io.bjpfkmjn](https://dx.doi.org/10.17504/protocols.io.bjpfkmjn)

Angel A Justiz-Vaillant<sup>1</sup>

<sup>1</sup>University of the West Indies St. Augustine

University of the West In...

angel.vaillant@sta.uwi.e...



Angel A Justiz-Vaillant

University of the West Indies St. Augustine

---

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.bjpfkmjn](https://dx.doi.org/10.17504/protocols.io.bjpfkmjn)

**Protocol Citation:** Angel A Justiz-Vaillant 2020. Enzyme linked immunosorbent assay for investigating the binding of Protein-L to diverse immunoglobulins . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bjpfkmjn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 13, 2020

**Last Modified:** August 13, 2020

**Protocol Integer ID:** 40391



## Abstract

This SpL ELISA can be used to detect specific antibodies in various animal species including human, mouse, rat, dog, rabbit, chicken, monkey, pig and hamster [1].

1. De Chateau M, Nilson BH, Erntell M, Myhre E, Magnusson CG, Akerstrom B et al. On the interaction between protein L and immunoglobulins of various mammalian species. *Scand J Immunol* 1993; 37: 399–405

- 1 This ELISA is used to study the interaction of protein L with different immunoglobulin preparations.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 1 µl/mg per well of unlabelled protein-L (SpL) from *P. magnus* in carbonate-bicarbonate buffer pH 9.6.
- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 µl of animal serum (1 mg/ml) is added and incubated for 1h at room temperature and washed.
- 5 Then 50 µl of peroxidase-labeled SpL conjugate diluted 1:5000 in PBS-non-fat milk is added to each well and incubated for 1h at RT. The plate is washed 4X with PBS-Tween.
- 6 50 µl of 4 mg/ml o-phenylenediamine solution (OPD) is added and the plate is incubated 15 minutes at RT in the dark.
- 7 The reaction is stopped with 50 µl of 3M H<sub>2</sub>SO<sub>4</sub> solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 492 nm.
- 9 A cut-off point should be calculated as the mean of the optical density of negative controls x 3.