Aug 13, 2020

Enzyme linked immunosorbent assay for investigating the binding of protein-LA (SpLA) to immunoglobulins.

DOI

dx.doi.org/10.17504/protocols.io.bjpjkmkn

Angel A Justiz-Vaillant¹, Norma McFarlane-Anderson²

¹University of the West Indies St. Augustine; ²University of West Indies. Mona Campus

University of the West In... angel.vaillant@sta.uwi.e...



Angel A Justiz-Vaillant

University of the West Indies St. Augustine





DOI: dx.doi.org/10.17504/protocols.io.bjpjkmkn

Protocol Citation: Angel A Justiz-Vaillant, Norma McFarlane-Anderson 2020. Enzyme linked immunosorbent assay for investigating the binding of protein-LA (SpLA) to immunoglobulins. . **protocols.io**

https://dx.doi.org/10.17504/protocols.io.bjpjkmkn

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: August 13, 2020

Last Modified: August 13, 2020

Protocol Integer ID: 40395

Abstract

This SpL ELISA can be used to detect specific antibodies in various animal species including human, mouse, rat, dog, rabbit, chicken, monkey, pig and hamster [1].

1. De Chateau M, Nilson BH, Erntell M, Myhre E, Magnusson CG, Akerstrom B et al. On the interaction between protein L and immunoglobulins of various mammalian species. Scand J Immunol 1993; 37: 399–405

- 1 This ELISA is used to study the interaction of recombinant protein LA (SpLA) with different immunoglobulin preparations.
- 2 The 96 well microtitre plate is coated overnight at 4°C with 2 μg/μl per well of SpLA in carbonate-bicarbonate buffer pH 9.6.
- 3 Then plate is treated with bovine serum albumin solution and washed 4X with PBS-Tween.
- 4 50 μl of animal serum (1 mg/ml) is added and incubated for 1h at room temperature and the microplate is rewashed 4X with PBS-Tween.
- 5 Then 50 μl of peroxidase-labeled SpLA conjugate diluted 1:5000 in PBS-non-fat milk is added to each well and incubated for 1h at RT. The plate is washed 4X with PBS-Tween.
- 6 50 μl of 4 mg/ml o-phenylenediamine solution (OPD) is added and the plate is incubated 15 minutes at RT in the dark.
- 7 The reaction is stopped with 50 μ l of 3M H2SO4 solution.
- 8 The plate is visually assessed for the development of colour and read in a microplate reader at 492 nm.
- A cut-off point should be calculated as the mean of the optical density of negative controls x 3. The higher the OD value the higher the affinity of SpLA to immunoglobulin G.