

Nov 14, 2025

Epi Multiome ATAC + Gene Expression – Stanford Bone Marrow TMC

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkyjw5g5r/v1>

Patricia Favaro¹, Jumana Afaghani¹, Sean Bendall¹

¹Stanford University

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu



Patricia Favaro

Stanford University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkyjw5g5r/v1>

Protocol Citation: Patricia Favaro, Jumana Afaghani, Sean Bendall 2025. Epi Multiome ATAC + Gene Expression – Stanford Bone Marrow TMC. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkyjw5g5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 12, 2025

Last Modified: November 16, 2025

Protocol Integer ID: 232249

Keywords: ATAC, HuBMAP, Tissue Map Center (TMC), Multiome, Gene Expression, bone Marrow, 10x genomics single cell multiome nuclei isolation protocol, 10x genomics single cell multiome, epi multiome atac, quality nuclei suitable for downstream atac, stanford bone marrow tmc fresh, human whole bone marrow aspirate, gene expression multimodal profiling, gene expression, mononuclear cell, downstream atac, nuclei isolation protocol

Funders Acknowledgements:

HuBMAP NIH

Grant ID: U54HL165445

Abstract

Fresh human whole bone marrow aspirates (25 mL; AllCells) were collected and processed on the same day. Mononuclear cells were isolated by Ficoll density gradient centrifugation followed by red blood cell lysis. CD34⁺ hematopoietic stem and progenitor cells were enriched by magnetic selection (STEMCELL Technologies, catalog # 17856). Equal numbers of CD34⁺ and CD34⁻ cells (5×10^5 each) were then pooled at a 1:1 ratio. Nuclei were extracted following the 10x Genomics Single Cell Multiome nuclei isolation protocol (see links below), generating high-quality nuclei suitable for downstream ATAC + gene expression multimodal profiling.



Steps

- 1 Nuclei isolation from fresh CD34⁺ and CD34⁻ cells: [link](#)
- 2 Chromium Next GEM Single Cell Multiome ATAC + Gene Expression library preparation: [link](#)

Acknowledgements

This work was supported by HuBMAP NIH grant U54HL165445.