




Jun 06, 2018

Version 5

Epifluorescent Microscopy of Virus Particles Using SYBR Green V.5

DOI

dx.doi.org/10.17504/protocols.io.qsgdwbw



Steven W Wilhelm¹, Samantha R Coy¹, Alyssa Alsante¹

¹The University of Tennessee, Knoxville

The Aquatic Microbial E...

Samantha R Coy's Proto...



Ashley A Humphrey

University of Tennessee, Knoxville

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.qsgdwbw>

External link: <https://doi.org/10.1371/journal.pone.0184371>

Protocol Citation: Steven W Wilhelm, Samantha R Coy, Alyssa Alsante 2018. Epifluorescent Microscopy of Virus Particles Using SYBR Green. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.qsgdwbw>



Manuscript citation:

Gainer PJ, Pound HL, Larkin AA, LeClerc GR, DeBruyn JM, Zinser ER, Johnson ZI, Wilhelm SW (2017) Contrasting seasonal drivers of virus abundance and production in the North Pacific Ocean. PLoS ONE 12(9): e0184371. doi: [10.1371/journal.pone.0184371](https://doi.org/10.1371/journal.pone.0184371)

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 06, 2018

Last Modified: June 06, 2018

Protocol Integer ID: 12840

Keywords: epifluorescent microscopy of virus particle, determination of virus abundance, virus particle, epifluorescence microscopy, epifluorescent microscopy, virus abundance, microscopy, virus, using sybr green adapted, sybr green adapted, molecular biology

Abstract

Adapted from: Ortmann and Suttle (2009) *Determination of Virus Abundance by Epifluorescence microscopy*. Ch. 10 Methods of Molecular Biology.

Contact Dr. Steven Wilhelm (wilhelm@utk.edu) or Samantha Coy (srose16@vols.utk.edu) for additional information regarding this protocol.

Attachments



Epifluorescent Micro...

127KB

Troubleshooting



Solutions

- 1 Make a 1% p-phenylenediamine (antifade) solution -- 0.01g + 1mL Milli-Q water.
- 2 Vortex, keep in the dark.
- 3 Retrieve an 8 uL stock of 1,000x SYBR Green, DI water, and 225 uL 50:50 PBS/Glycerol solutions
- 4 Add 6 uL SYBR Green to 235 uL filtered DI. Keep in the dark.
- 5 Add 25 uL 1% p-phenylenediamine to the 225 uL aliquot of 50:50 PBS/Glycerol. Keep in the dark.
- 6 Place the SYBR Green and antifade solutions into a drawer while working.

Note

Only the red light can be on when the solutions are out of the drawer.

Sample Preparation

- 7 Fix sample with 0.5x glutaraldehyde.

Note

If you are going to freeze your sample between slide preparation and visualization, then you need to cool the fixed sample in the dark at 4°C for 15-30 min, and flash freeze in liquid nitrogen before storing at -80°C.

Vacuum Setup

- 8 Turn on the vacuum (no more than ~25 mm Hg).
- 9 Rinse wells with Milli-Q water.



- 10 While the knobs are turned to open (vacuum on), use the tweezers to add the 0.45 μ m nitrocellulose backing filter.

Note

This filter is used to protect the more fragile 0.02 μ m Anodisc filter that the viruses are collected on.

- 11 Once the filter is on, you should add a few drops of Milli-Q and then turn the knobs to closed (vacuum off).
- 12 Use tweezers to apply the 0.02 μ m Anodisc filter to the top of the backing filter.
- 13 Turn the vacuum back on by turning the knobs to open and add the sample in two 500 μ L aliquots.


SYBR Green

- 14 Pipette SYBR Green as dots onto a glass petri dish. Make enough dots for each filter.

 30 μ L

- 15 With the vacuum on, use a needle and tweezers to remove the Anodisc filter and place onto a dot. Turn the vacuum off.

- 16 Put the lid over the plate and back into the drawer to stain for 20 min.

 00:20:00

Slide Preparation

- 17 During the wait, label slides with sample ID, dilution, stain, sample volume, and date.

Note

You may turn the white light on during this period switching it off anytime the stain is exposed to light.

- 18 Re-apply Milli-Q to the backing filter and place the stained Anodisc on top. Add another drop of Milli-Q.
- 19 Turn the vacuum on again. After the stain has been washed away, use the needles and tweezers to remove the Anodisc--again, with the vacuum still on.
- 20 Place the filter onto a Kim Wipe to dry.
- 21 On the labeled side, add 18 uL antifade solution, apply the dried Anodisc filter and apply another 18 uL of antifade to the top.
- 22 Apply a coverslip and lightly press down so there are no bubbles.

Storage

- 23 Store at -20°C in the dark for up to four months (3-4 weeks recommended). .

Note

Thawing the sample repeatedly will result in fading of the stain, but you can attempt to re-stain the filter if it fades

Microscope

- 24 Turn on all three computers and the shutter.
- 25 Open the LAS AF program.
- 26 On the second menu, use the FLUO--Filter cubes L5 and the FLUO incident settings.
- 27 On the fourth menu, click set/clear focus position to desired bottom using the up and down arrows. Once positioned, click "set". You can then move the microscope slide to the desired position to visualize the slide and click "set". These settings will then be saved within the microscope.

- *Z: up/down (turn the back knob)
- *X/Y: Side to side (turn the top or bottom half of the front knob)

- 28 Once the options are set and the microscope slide is set into place, you can then turn the shutter on/off on the first menu--'IL--Shutter ON/OFF'.
- 29 Rod all of the way in on the microscope -- image only on the microscope itself
Rod half way in -- image on the microscope and the computer screen
Rod all of the way out -- image only on the computer screen

Computer Imaging Software

- 30 Under the acquire → acquisition tab:
*To look at a preview of your image, click "live" (delay may occur due to exposure)
*Set desired exposure, gain and intensity.
*To save a desired image, click "capture image". The image will appear on the right screen, where you can adjust the background fluorescence by moving the bar on the black to green scale at the bottom lefthand corner.
- 31 Under the acquire → experiment tab:
*The captured images will automatically save to the experiment .lif file. This can only be viewed using this software.
*Right click on an image to rename/delete it.
*Click "save all" to save images to the desired directory.
*To save images in a different format, right click, select export, and choose your desired format (e.g. JPEG, TIFF, or the format required for your analysis software).

Enumeration of Virus-Like Particles

- 32 With the rod all of the way in, use the fine adjustment knob to locate and resolve the viruses into focus.
- 33 Count all of the viruses in the grid. The ideal concentration is between 20-40 viruses per grid so that you have enough time to count the viruses before they fade.
- 34 After counting all of the viruses in the field, look up from the microscope, and move the stage to a different part of the slide to count viruses in another field. It is important to do this without looking so that you can randomly sample your filter and not incur a bias in your counts.
- 35 Repeat this process until you have counted enough fields to total between 200-300 viruses. Calculate the average, and determine the VLP/mL by the equation:



(Avg VLP/field)*(28352.9)*(Dil. factor)*(1000/volume filtered (uL))*((1000+glu. volume (uL))/1000)