

Jun 10, 2025

 **eRNA and eDNA co-extraction from water samples filtered through 47mm diameter filters (NucleoMag DNA/RNA Water Kit - MACHEREY NAGEL).**



DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1dr97vr2/v1>

Marine Vautier<sup>1</sup>, Cecile Chardon<sup>1</sup>, Cyrielle GALIEGUE<sup>1</sup>, Amelie Bernard<sup>1</sup>, Isabelle Domaizon<sup>1</sup>

<sup>1</sup>INRAE - Univ. Savoie Mont Blanc - UMR CARRETEL - Pole R&D Ecla

CARRETEL



Marine Vautier

INRAe

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1dr97vr2/v1>

**Protocol Citation:** Marine Vautier, Cecile Chardon, Cyrielle GALIEGUE, Amelie Bernard, Isabelle Domaizon 2025. eRNA and eDNA co-extraction from water samples filtered through 47mm diameter filters (NucleoMag DNA/RNA Water Kit - MACHEREY NAGEL).. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1dr97vr2/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** February 26, 2025

**Last Modified:** June 10, 2025

**Protocol Integer ID:** 123445

**Keywords:** water, extraction, filter, RNA, eRNA, magnetapure 32 nucleic acid purification system, nucleic acid purification system, rna extraction, dna extraction, rna water kit, reversible adsorption of nucleic acid, environmental rna, nucleic acid, nucleomag dna, rna, nucleomagdna, extraction, purification, rdnase digestion step, dna, water sample, molecular biology analysis

## Abstract

The objective of this protocol is the **environmental RNA (eRNA) and DNA (eDNA) co-extraction from water samples filtered through 47 mm diameter filters.**

RNA extraction is performed using a **MagnetaPure 32 Nucleic Acid Purification System** (Dutscher) and with the **NucleoMagDNA/RNA Water Kit** (Macherey Nagel), with a rDNase digestion step.

DNA extraction is performed using a **MagnetaPure 32 Nucleic Acid Purification System** (Dutscher) and with the **NucleoMagDNA/RNA Water Kit** (Macherey Nagel).

The procedure is based on **reversible adsorption of nucleic acids to paramagnetic beads** under appropriate buffer conditions.

The benefits of using the **MagnetaPure 32 Nucleic Acid Purification System** (Dutscher) are **increased productivity and repeatability, as well as eliminating human error and the pain of repetitive work.**

This protocol is used **prior to molecular biology analysis** (e.g. qPCR, metabarcoding, ddPCR) to specifically target **both macro- and micro-organisms eRNA and eDNA extracted from water samples.**

**This protocol is optimised for rare eDNA/RNA** and the suggested **elution volumes are therefore low** (between 50 and 65  $\mu$ L), but can be increased if targets are more abundant.

## Image Attribution

Marine Vautier

Amélie Bernard

## Guidelines

### **The main steps of the protocol are:**

- Material and rDNase preparation
- Plate preparation
- Sample lysis
- eRNA/eDNA co-extraction with the MagnetaPure 32 System (Dutscher)
- eRNA/eDNA elution

## Materials

### ▪ **Materials:**

- 1000 µL pipette
- 100 µL pipette
- 10 µL pipette
- Scissors
- Vortex + benchtop centrifuge for 5 mL tubes
- Horizontal vortex with 5 mL tube holder (15 mL tube holders)
- Centrifuge for 2 mL tubes (relative centrifugal force needed: 11,000 x g)
- MagnetaPure 32 Nucleic Acid Purification System (Dutscher)
- Specific DNA/RNA-workstation (sterile area equipped with air filtration and UV systems)

### ▪ **Consumables:**

*All tubes and tips must be sterile and open under the workstation*

- 1000 µL tips with filter
- 100 µL tips with filter
- 10 µL tips with filter
- 50 mL tubes: 4 to prepare aliquots, 3 for scissors decontamination
- 5 mL tubes: 1 per 8 samples to prepare NucleoMag B-Beads and MWA2 mix, 1 for rDNase dilution
- 2 mL tubes: 1 per sample to transfer lysate + 2 to prepare aliquots
- 0.5 mL tubes: 3 per sample to transfer eluted RNA
- 96-well plate with 2 mL deep-wells, U-Bottom (Macherey Nagel - 746032.DEEP): 1 per 16 samples RNA and 1 per 16 samples DNA
- Magnetic rod cover for MagnetaPure 32 (Macherey Nagel - 747032.TC): 1 per 8 samples
- 2 Plastics film to protect the 96-well plate
- Gloves

*For any manipulation in a rare DNA/RNA room, provide complete equipment (disposable coat, cap, mask, shoe covers & gloves).*

### ▪ **Reagents:**

- rDNase Set (Macherey Nagel - 740963)
- Reaction Buffer for rDNase (Macherey Nagel - 740834.60)
- NucleoMag DNA/RNA Water kit (Macherey Nagel)

*Note: shelf life of reagents for 24 months from production*

- To clean the scissors:
  - o Solution to remove RNA (e.g. *RNase-off*/*DNA-off*, *RNA/DNA away*)

- o H<sub>2</sub>O<sub>2</sub> 10 %
- o Ethanol 96% for molecular biology
- o Molecular biology grade water

- **Samples to be extracted:**

- Filters in 5 mL tubes and stored at -80°C until RNA/DNA extraction

## Troubleshooting

## Before start

- **Filtration and preservation of water samples through 47 mm diameter filters.**
- **As eRNA degrades very rapidly, the water samples must be filtered and preserved very quickly after sampling** (less than one hour if possible). If filtration is carried out in the field, the filters should be immediately frozen in liquid nitrogen and then stored at  $-80^{\circ}\text{C}$  until the eRNA/eDNA are extracted. If filtration is carried out in the laboratory, the filters can be stored at  $-80^{\circ}\text{C}$  immediately after filtration and until the eRNA/eDNA are extracted.

*Note: As eRNA is very sensitive to RNase, which is active at temperatures above  $-70^{\circ}\text{C}$ , it is very important to freeze the filters in liquid nitrogen or at  $-80^{\circ}\text{C}$  very quickly after filtration.*

This protocol is suitable for the RNA extraction from different types of 47 mm diameter filters:

Filter Type
Polyethersulfone (PES)
Polycarbonate (PC)
Cellulose nitrate (NC)
Cellulose mixed esters (CM)
Cellulose acetate (CA)

- The filters should be preserved in 5 mL tubes for this RNA/DNA co-extraction protocol.

- **The following precautions must be applied:**

- Wear gloves throughout the extraction process
  - - Clean the bench with a RNA/DNA-removing solution (e.g. RNase-off/DNase-off).
- Use tips with filters to avoid contaminations
- All steps have to be performed under a specific DNA/RNA-work station (sterile area equipped with air filtration and UV systems)

*For any manipulation in a rare DNA/RNA room, provide complete equipment (disposable coat, cap, mask, shoe covers & gloves).*

- **Pre and post extraction equipment decontamination:**

- Clean a specific DNA/RNA-work station and apply UV
- MagnetaPure 32 System (Dutscher): Visual check for residues to be removed and UV decontamination – see *instrument manual*

- **Scissors decontamination** (*to be done before starting the protocol, and between each filter cutting*)
  - Wash them with RNase-off/DNase-off solution and rinse with ethanol
  - Prepare :
    1. one 50 mL tube with H<sub>2</sub>O<sub>2</sub> at 10%
    2. one 50 mL tube with molecular biology grade water
    3. one 50 mL tube with ethanol
  - Successively dip, opening and closing the scissors into each tube.

*Note: The cutting of filters is not a requirement, but it does have an impact on the extraction efficiency.*



## Material and rDNase preparation

### 1 Material preparation

- *To limit contamination of the kit buffers, it is recommended to aliquote them:*
  - Into 50 mL tubes for MWA1, MWA2, MWA3 and MWA4
  - Into 2 mL tubes for NucleoMag B-Beads solution and RNase-free H<sub>2</sub>O
  
- *Tubes annotation*
  - one 2 mL tube per sample for lysate collection
  - Two 5 mL tubes per 8 samples for the NucleoMag B-Beads and MWA2 mix preparation
  - Two 1.5 mL tubes per sample for RNA and DNA collection
  
- *Scissors decontamination (to be done before starting the protocole, and between each filter cutting)*
  - Wash them with RNase-off/DNase-off solution and rince with ethanol
  - Prepare :
    1. one 50 mL tube with H<sub>2</sub>O<sub>2</sub> at 10%
    2. one 50 mL tube with molecular biology grade water
    3. one 50 mL tube with ethanol
  - Successively dip the scissors into each tube, opening and closing the scissors in each tube.

*Note: The cutting of filters is not a requirement, but it does have an impact on the extraction efficiency.*

### 2 rDNase preparation

5m

#### 1. Reconstitution of lyophilized rDNase:

- Add  4 mL of **reaction buffer for rDNase** to the rDNase vial
- Incubate for  00:02:00 -  00:03:00 at  Room temperature
- Gently shake the vials to completely dissolve the rDNase

*Note: Be careful not to mix the rDNase vigorously because the rDNase is sensitive to mechanical agitation*

#### 2. rDNase Dilution:

- Transfer the **reconstituted rDNase** to a 50 mL tube
- Add  28 mL **of reaction buffer for rDNase**
- Gently shake the tube

*Note: The resulting rDNase solution will be sufficient for 96 samples. If fewer than 96 samples are to be processed, it is advisable to aliquot the rDNase as it should not be frozen and thawed more than three times. Aliquots can be store at  -20 °C for at least 6 months.*

**Note: Put one aliquot in the fridge a few hours before to allow slow defrosting**

## eRNA Plate preparation 1/2

3

30m

*In this step, the buffers provided by the kit are distributed in a plate.*

*For RNA extraction from filtered water samples, the 12 columns of the plate are divided into 2 sections of 6 columns each, allowing up to 16 samples to be extracted per plate.*

- Annotate the plate as recommended below:

	1	2	3	4	5	6	7	8	9	10	11	12
Solution	rDNase reaction mixture	MWA3	MWA3	NucleoMag B-Beads and MWA2 mix + Lysate	MWA4	RNase-free H <sub>2</sub> O	rDNase reaction mixture	MWA3	MWA3	NucleoMag B-Beads and MWA2 mix + Lysate	MWA4	RNase-free H <sub>2</sub> O
Column name	D	3	3	L	4	E	D	3	3	L	4	E

Samples between 1 to 8
Samples between 9 to 12

**Table 1:** Recommended eRNA plaque annotation

*Note: It is useful to mark the dividing line between columns 6 and 7 with a marker pen to provide a visual cue for filling the plate.*

- Add the appropriate buffers into the appropriate wells of the plate (Table 2):

**1<sup>st</sup>column / 7<sup>th</sup>column:** Will be filled during plate preparation 2/2

**2<sup>nd</sup>column / 8<sup>th</sup>column:**  850 µL of **MWA3** (1<sup>st</sup> wash)

**3<sup>rd</sup>column / 9<sup>th</sup>column:**  850 µL of **MWA3** (2<sup>nd</sup> wash)

**4<sup>th</sup>column / 10<sup>th</sup>column:** Will be filled during plate preparation 2/2

**5<sup>th</sup>column / 11<sup>th</sup>column:**  850 µL of **MWA4** (3<sup>rd</sup> wash and bead drying)

Note: The choice of elution volume is based on the expected eRNA amount. The smaller the amount, the smaller the elution volume in order to obtain more concentrated RNA.

	1	2	3	4	5	6	7	8	9	10	11	12
Solution	Plate preparation 2/2	MWA3	MWA3	Plate preparation 2/2	MWA4	Plate preparation 2/2	Plate preparation 2/2	MWA3	MWA3	Plate preparation 2/2	MWA4	Plate preparation 2/2
Volume (µL)		850	850		850			850	850		850	
Column name		3	3		4			3	3		4	

Samples between 1 to 8
Samples between 9 to 12

**Table 2:** eRNA Plate preparation 1/2

- Film and reserve the plate at Room temperature

### MWA2 preparation for RNA and DNA plate

- Preparation of the NucleoMag B-Beads and MWA2 mix:
    - Prepare two 5 mL tube / maximum 8 samples (one for the RNA plate and one for the DNA plate)
    - Add **MWA2** only, the NucleoMag B-beads will be added during plate preparation 2/2 (*allow a margin of one sample or few percent for the mix preparation. For example: plan a mix for 9 samples if 8 samples are to be extracted, or add a margin of 10 %*)

Sample Number	Volume NucleoMag B-Beads (µL)	Volume MWA2 (µL)
1	25	475
2	50	950
3	75	1425
4	100	1900
5	125	2375
6	150	2850
7	175	3325
8	200	3800
9	225	4275

**Table 3:** Volume required to prepare NucleoMag B-beads and MWA2 mix (no margin)

## Sample Lysis for co-extraction

50m 5s

5 *During this step, a mechanical and chemical lysis of the sample is performed.*

50m 5s

▪ **For samples prefilled with buffer:**

- Collect the 5 mL tubes containing the filters

Note: If the tubes containing the filters are frozen, defrost them  00:10:00 at

 Room temperature

- Cut the filter into small pieces directly into the tube using decontaminated scissors

Note: Scissors must be decontaminated between each sample (see section 1 "material preparation")

- Place the tubes on the vertical vortex  00:05:00 at maximum speed

- Place the tubes into the benchtop centrifuge

- Pipette the lysate from the tube and transfer it into a 2 mL tube

- Centrifuge at  11000 rpm  00:30:00

- Replace the tubes into the rack and reserve them at  Room temperature until their distribution into the plate

▪ **For samples without buffer:**

- Collect the 5mL tubes containing the filters from the freezer and place them

 On ice

- Add immediately  1100  $\mu\text{L}$  of **MWA1** buffer into each tube

- Cut the filter into small pieces directly into each tube using decontaminated scissors

*Note: Scissors must be decontaminated between each sample (see section 1 "material preparation")*

- Place the tubes on the vertical vortex  00:00:05 at median speed

- Place the tubes on the horizontal vortex  00:05:00 at maximum speed

- Place the 5 mL tubes into the benchtop centrifuge

- Pipette the lysate from the tube (approximately  900  $\mu\text{L}$  ) and transfer it into a 2 mL tube

- Centrifuge at  11000 x g, 00:00:30

- Replace the tubes into the rack and reserve them at  Room temperature until their distribution into the plate

## eRNA Plate preparation 2/2

6 ▪ Preparation of the NucleoMag B-Beads and MWA2 mix

*NucleoMag B-Beads and MWA2 mix sediment quickly, vortex between each samples to ensure homogeneity*

- Vigorously vortex **NucleoMag B-Beads** tube
- For each tube containing MWA2 buffer (previously prepared), add the appropriate volume of **NucleoMag B-Beads** (see Table 4 above)
- Vortex
- NucleoMag B-Beads - MWA2 mix and lysate distribution
- Remove the film from the plate
- Add the appropriate solution into each well of the plate

**4<sup>th</sup> column / 10<sup>th</sup> column:**

 500 µL of **NucleoMag B-Beads and MWA2 mix**

 450 µL of **Lysate** (supernatant from the 2 mL tubes)

**6<sup>th</sup> column / 12<sup>th</sup> column:**  50 µL or  65 µL of **RNase-free H<sub>2</sub>O** (RNA elution)

- rDNase mix distribution
- Slowly agitate the **rDNase mix** tube (previously prepared- see step 2)
- Add the appropriate solution into each well of the plate

**1<sup>st</sup> column / 7<sup>th</sup> column:**

 300 µL of **rDNase mix**

	1	2	3	4	5	6	7	8	9	10	11	12
Solution	rDNase reaction mixture	Plate preparation 2/2	Plate preparation 2/2	NucleoMag B-Beads and MWA2 mix	Plate preparation 2/2	DNase-free H <sub>2</sub> O	rDNase reaction mixture	Plate preparation 2/2	Plate preparation 2/2	NucleoMag B-Beads and MWA2 mix	Plate preparation 2/2	DNase-free H <sub>2</sub> O
Volume (µL)	300			500		50 or 65	300			500		50 or 65
Solution			Lysate							Lysate		
Volume (µL)			450							450		
Column name	D		L			E	D			L		E

Samples between 1 to 8
Samples between 9 to 12

Table 4: eRNA Plate preparation 2/2

- Film the plate and proceed immediately to the next step.

## eRNA extraction with the MagnetaPure 32 System (Dutscher) 1/2

35m

### 7 ■ Extraction part 1

35m

- Place the plate into the MagnetaPure 32 System and insert the magnetic rod coverfor – see *instrument manual*

#### - Select the appropriate program for RNA extraction part 1

Step	Well	Name	Mix time (min)	Magnet (sec)	Wait time (min)	Volume (µL)	Mix speed	Temp (°C)
1	4	Binding	8	10	0	950	8	OFF
2	3	MWA3	2,3	5	0	850	8	OFF
3	2	MWA3	2,3	5	5	850	8	OFF
4	1	Digest DNA	2,3	0	15	300	8	37

Table 5: MagnetaPure 32 System program for NucleoMag DNA/RNA Water Kit - RNA extraction part 1

- Start the run (*The run lasts approximately*  00:35:00 )

## eDNA Plate preparation 1/2

8 In this step, the buffers provided by the kit are distributed in a plate. For DNA extraction from filtered water samples, the 12 columns of the plate are divided into 2 sections of 6 columns each, allowing up to 16 samples to be extracted per plate.

9 Annotate the plate as recommended below:

	1	2	3	4	5	6	7	8	9	10	11	12
Solution	Plate preparation 2/2	MWA3	MWA3	MWA4	*	Plate preparation 2/2	Plate preparation 2/2	MWA3	MWA3	MWA4	*	Plate preparation 2/2
Column name		3	3	4	*			3	3	4	*	

Samples between 1 to 8
Samples between 9 to 12

Table 6: Recommended eDNA plaque annotation

Note: It is useful to mark the dividing line between columns 6 and 7 with a marker pen to provide a visual cue for filling the plate.

- Add the appropriate buffers into the appropriate wells of the plate.

**1<sup>st</sup> column / 7<sup>th</sup> column:** Will be filled during plate preparation 2/2

**2<sup>nd</sup> column / 8<sup>th</sup> column:**  850 µL of **MWA3** (1<sup>st</sup> wash)

**3<sup>rd</sup> column / 9<sup>th</sup> column:**  850 µL of **MWA3** (2<sup>nd</sup> wash)

**4<sup>th</sup> column / 10<sup>th</sup> column:**  850 µL of **MWA4** (3<sup>rd</sup> wash and bead drying)

**5<sup>th</sup> column / 11<sup>th</sup> column:** not used

**6<sup>st</sup> column / 12<sup>th</sup> column:** Will be filled during plate preparation 2/2

*Note: The choice of elution volume is based on the expected eDNA amount. The smaller the amount, the smaller the elution volume in order to obtain more concentrated DNA.*

	1	2	3	4	5	6	7	8	9	10	11	12
Solution	Plate preparation 2/2	MWA3	MWA3	MWA4	*	Plate preparation 2/2	Plate preparation 2/2	MWA3	MWA3	MWA4	*	Plate preparation 2/2
Volume (µL)		850	850	850	*			850	850	850	*	
Column name		3	3	4	*			3	3	4	*	

Samples between 1 to 8
Samples between 9 to 12

Table 7: eDNA Plate preparation 1/2

- Film and reserve the plate at  Room temperature

## eRNA extraction with the MagnetaPure 32 System (Dutscher) 2/2

30m

- 10
- At the end of DNA digestion
    - Remove the plate and add the appropriate solution into each well of the plate

30m

**1<sup>st</sup> column / 7<sup>th</sup> column:**  350 µL of **MWA2**

	1	2	3	4	5	6	7	8	9	10	11	12
Solution	MWA2						MWA2					
Volume (µL)	350						350					
Column name	D						D					

Samples between 1 to 8
Samples between 9 to 12

Table 9: Addition of MWA2 buffer after DNA digestion (eRNA plate)

▪ **Extraction part 2**

- Place the plate back into the MagnetaPure 32 System (it is not necessary to change the magnetic rod coverfor)

- **Select the appropriate program for RNA extraction part 2**

Step	Well	Name	Mix time (min)	Magnet (sec)	Wait time (min)	Volume (µL)	Mix speed	Temp (°C)
1	1	Rebind	8	10	0	650	8	OFF
2	5	MWA4	2,3	5	15	850	8	OFF
3	6	Elution	5	20	0	50 or 65	8	OFF
4	5	Release	0,5	0	0	850	10	OFF

Table 10: MagnetaPure 32 System program for NucleoMag DNA/RNA Water Kit - RNA extraction part 2

- Start the run (*The run lasts approximately*  00:30:00 )

## eDNA Plate preparation 2/2

30s

11 ▪ Preparation of the NucleoMag B-Beads and MWA2 mix

30s

NucleoMag B-Beads and MWA2 mix sediment quickly, vortex between each samples to ensure homogeneity

- centrifuge the lysate (prepare in step 5) at  11000 rpm during  00:00:30

- Vigorously vortex **NucleoMag B-Beads** tube

- For each tube containing MWA2 buffer (previously prepared), add the appropriate volume of **NucleoMag B-Beads** (see Table 3 above)

- Vortex

▪ NucleoMag B-Beads - MWA2 mix and lysate distribution

- Remove the film from the plate

- Add the appropriate solution into each well of the plate

**1<sup>st</sup>column / 7<sup>th</sup>column:**

 500 µL of **NucleoMag B-Beads** and **MWA2 mix**

 450 µL of **Lysate** (supernatant from the 2 mL tubes) same as eRNA lysate (step 5)

**6<sup>th</sup>column / 12<sup>th</sup>column:**  50 µL or  65 µL of **DNase-free H<sub>2</sub>O** (DNA elution)

	1	2	3	4	5	6	7	8	9	10	11	12
Solution	NucleoMag B-Beads and MWA2 mix	Plate preparation 1/2	Plate preparation 1/2	Plate preparation 1/2	x	DNase-free H2O	NucleoMag B-Beads and MWA2 mix	Plate preparation 1/2	Plate preparation 1/2	Plate preparation 1/2	x	DNase-free H2O
Volume (µL)	500				x	50 or 65	500				x	50 or 65
Solution	Lysate				x		Lysate				x	
Volume (µL)	450				x		450				x	
Column name	L				x	E	L				x	E

Samples between 1 to 8
Samples between 9 to 12

Table 11: eDNA Plate preparation 2/2

- Film and reserve the plate at Room temperature

## Transfer of eRNA extracts

30m

- 12
- At the end of the run, remove the plate and place it On ice into the DNA/RNA-workstation
  - Remove the magnetic rod cover and start UV for decontamination – *see instrument manual*
  - In the DNA/RNA-workstation, transfer each **eRNA extract** into a 1.5 mL tube previously annotated.

*Note: It is recommended to aliquot the eRNA extracts directly to avoid freezing and thawing cycles that could degrade the eRNA.*

*Note: eRNA concentration and quality can be measured at this step (e.g. Nanodrop)*

- Store eRNA extracts at -80 °C

## eDNA Extraction step performed in the MagnetaPure 32 System

40m

- 13
- Place the plate into the MagnetaPure 32 System and insert the magnetic rod cover for – see instrument manual
  - **Select the appropriate program to the chosen elution volume and elution temperature**

40m

Step	Well	Name	Mix time (min)	Magnet (sec)	Wait time (min)	Volume (µL)	Mix speed	Temp (°C)
1	1	Binding	8	10	0	950	8	OFF
2	2	MWA3	2,3	5	0	850	8	OFF
3	3	MWA3	2,3	5	0	850	8	OFF
4	4	MWA4	2,3	5	15	850	8	OFF
5	6	Elution	5	20	0	50 or 65	8	56°C or OFF
6	4	Release	0.5	0	0	850	10	OFF

**Table 12:** MagentaPure 32 System program for NucleoMag DNA/RNA Water Kit DNA extraction

Note: Heating to 56°C during elution gives a higher yield of DNA, but there is a risk of evaporation which reduces the volume of eluted DNA recovered

- Start the run (The run lasts approximately  00:40:00 )

## Transfer of eDNA extracts

- 14
- At the end of the run, remove the plate and place it into the DNA-workstation
  - Remove the magnetic rod cover and start UV for decontamination – see instrument manual
  - In the DNA-workstation, transfer each **eDNA extract** into a 1.5 mL tube previously annotated

Note: DNA concentration and quality can be measured at this step (e.g. Nanodrop)

- Store DNA extracts at  4 °C for immediate use, or at  -20 °C . or  -80 °C for long-term preservation

## Protocol references

NucleoMag DNA/RNA Water kit (Macherey Nagel) manual : <https://www.mn-net.com/media/pdf/ce/b5/38/Instruction-NucleoMag-DNA-RNA-Water.pdf>

MagnetaPure 32 Nucleic Acid Purification System (Dutscher) manual : [https://pdf.dutscher.com/doc/255743/255743\\_MEen.pdf](https://pdf.dutscher.com/doc/255743/255743_MEen.pdf)

eDNA extraction from water samples filtered through 47 mm diameter filters (NucleoMag DNA/RNA Water Kit - MACHEREY NAGEL) : <https://www.protocols.io/view/edna-extraction-from-water-samples-filtered-through-j8nlkox7dv5r/v1>

eRNA extraction from water samples filtered through 47 mm diameter filters (NucleoMag DNA/RNA Water Kit - MACHEREY NAGEL) : <https://www.protocols.io/view/erna-extraction-from-water-samples-filtered-through-5qpvo3eq9v4o/v1>