



Dec 05, 2023

Essential surgeries for the electrophysiological recording from a behaving non-human primate brain 2 (Craniotomy)

DOI

dx.doi.org/10.17504/protocols.io.n92ldmz48l5b/v1

Daisuke Kase¹, Witold J Lipski¹, Devin R Harsch¹, Robert S Turner¹

¹Department of Neurobiology, University of Pittsburgh



Daisuke Kase

Department of Neurobiology, University of Pittsburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n92ldmz48l5b/v1>

Protocol Citation: Daisuke Kase, Witold J Lipski, Devin R Harsch, Robert S Turner 2023. Essential surgeries for the electrophysiological recording from a behaving non-human primate brain 2 (Craniotomy). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.n92ldmz48l5b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 01, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 85783

Keywords: ASAPCRN, Non-human primate, Craniotomy, Electrophysiology, essential surgeries for the electrophysiological recording, cranial implant, craniotomy, electrophysiological recording, essential surgery, implant surgery, implant, recording chamber, existing recording chamber, procedure, brain, place upon recovery

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research (MJFF)

Grant ID: ASAP-020519

NIH

Grant ID: R01 NS093164

NIH

Grant ID: R01 NS070865

NIH

Grant ID: R01 NS091853

NIH

Grant ID: R01 NS117058

Abstract

Outlined below are the steps to complete a craniotomy in an existing recording chamber placed on a cranial implant. This procedure can be done concurrently with the implant surgery (not recommended), or can take place upon recovery from the implant surgery.

Materials

	Item	Specifications	Vendor	Note
	Drill bit	40 Gauge	McMaster-Carr	https://www.mcmaster.com/3584A225/
	Pin vise		McMaster-Carr	https://www.mcmaster.com/8455A31/
	Bone rongeur	Jaw width: 1.3mm Length: 5.5"	Roboz	http://shopping.robz.com/micro-scissors-micro-forceps-groups/Rongeurs/Micro-Friedman-Rongeur-5-5-Very-Delicate-Curved
	Kerrison bone punches	Jaw width: 4mm Jaw position: 90 deg.	Aesculap	https://www.aesculapusa.com/en/healthcare-professionals/or-solutions/or-solutions-neurosurgical-instruments/kerrison-pneumatic-bone-punch.html
	Excavator spoon	#1 and #153-154	Nordent instruments	https://nordent.com/product/excavator-1-2/ https://nordent.com/product/excavator-english-pattern-153-154-2/



	Item	Specifications	Vendor	Note
	Gelfoam	12-7 mm	Pfizer	https://www.pfizermedicalinformation.com/en-us/gelfoam-absorbable-gelatin-powder
	HemaBlock		HemaBlock	https://hemablock.com/

List of key materials for the craniotomy surgery

Note: Bone punches with larger jaw width (5mm or 6mm) works well to expand the craniotomy

Troubleshooting

Preparation in the prep room on the day of surgery

- 1 Sedate the animal.
- 2 Intubate the animal, start the gas anesthesia, and bring it to the surgical suite.

Note

The intubation should be performed by a skilled veterinarian or technician.

Procedure in the surgical suite

- 3 Secure the animal's head on the stereotaxic frame.
- 4 Open the recording chamber.
- 5 Clean the surface of the implants (chamber, head fixation post, dental acrylic, etc.) with Betadine and ethanol.
- 6 Scrub in for the surgery. Wash hands (5 min per hand) and don a surgical gown and sterile gloves.
- 7 Cover the animal's whole body with a sterile drape.
- 8 Put a sterile drill bit in the pin vise, and expose the drill bit about 3-4 mm from the pin vise.

Note

You may remove some dental acrylic with dental drill if there is a thick layer of dental acrylic inside the chamber.



- 9 Slowly drill the first hole in the chamber with a pin vise until the pin vise touches the skull. The drill should not pierce the skull if the length is appropriately adjusted.

Note

Irrigate the skull and drill bit with sterile saline during this process.

- 10 Extend the drill bit a small amount, approximately 0.5 mm, then return to the hole that you drilled.

- 11 Deepen the hole.

Note

Continue to irrigate the skull and drill bit with sterile saline during this process.

- 12 Move to the next step if the skull is pierced and you can see the surface of the dura (it will look like a fascia of a piece of chicken). If the skull has not been pierced yet, return to Step 10.

- 13 Make many holes to draw a circle of holes in the recording chamber.

- 14 Remove the bone between holes with the bone rongeur.

Note

The excavator spoon is also useful for removing the skull in a narrow space.

- 15 If bones between all holes are entirely removed, there should be a circle-shaped bone in the middle of the recording chamber.

- 16 Raise the circle-shaped bone with tweezers and/or excavator spoon slowly and carefully.

Note

Do NOT push the circle-shaped bone that you have cut out down onto the dura. Be careful to ensure that no lever is created that will push onto the dura.

17 Remove the circle-shaped bone with tweezers carefully.

18 Expand the hole with the bone punches.

Note

1: Shave the edge of the hole with an excavator spoon to expand the hole if the hole is not large enough to insert the lower jaw of the bone punches.

2: Ensure there is no sharp bone piece in the hole when inserting the lower jaw of the bone punches into the space between the dura and skull.

19 Scrape the sharp edge of the skull with the excavator spoon.

Note

Bone punches cannot make a completely smooth, round shaped hole. If there is any sharp edge on the hole, such edges may push the dura and potentially cause the leak of the CSF in the future.

20 Rinse the small pieces of bone from the chamber with saline.

21 Measure the distance from the top edge of the chamber to the surface of the dura in multiple locations inside the chamber (ex. most medial, lateral, anterior, posterior corners, and center of the chamber).

22 Place the chamber's cap onto the chamber and secure it.

23 Stop the gas anesthesia, remove the intubation when the animal starts showing the reflex. Then, return the animal to the home cage.

