

Jun 13, 2019

Ethanol precipitation of RNA from small or large volumes

DOI

dx.doi.org/10.17504/protocols.io.36fgrbn

Stephen Floor¹

¹University of California, San Francisco

Stephen Floor Lab



Stephen Floor

University of California, San Francisco

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.36fgrbn>

Protocol Citation: Stephen Floor 2019. Ethanol precipitation of RNA from small or large volumes. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.36fgrbn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 13, 2019

Last Modified: June 13, 2019

Protocol Integer ID: 24487

Keywords: RNA, ethanol precipitation, ethanol precipitation of rna, rna with ethanol, procedure for large volume precipitation, lorsch rna purification, large volume precipitation, precipitation method, rna, initial precipitation, ethanol, purification, smaller tube, procedure, large tube, small tube

Abstract


How to precipitate RNA with ethanol and resuspend it! Two separate protocols - one for small and one for large volumes. The procedure for large volume precipitation is similar, but begins with an initial precipitation and transfers to a small tube for the final steps to avoid loss in the large tube. Strongly recommended to move to a smaller tube.

Recommended reading: Walker & Lorsch [RNA purification--precipitation methods](#).

Materials

MATERIALS

 GlycoBlue™ Coprecipitant **Thermo Scientific Catalog #AM9516**

 70% Ethanol

 100% Ethanol (KOPTEC) **VWR International (Avantor) Catalog #89125-186**

0.3 M NaOAc pH 5.2

3 M NaOAc pH 5.2

Note: the protocol requires 70% ethanol pre-chilled to -20 C and requires 100% ethanol pre-chilled to -20 C for large volumes.

Troubleshooting

Precipitating RNA (< 500 ul sample volume)

1 Record the sample volume. If unknown, you can estimate it using the side of the tube or measure it by pipetting it.

2 Add NaOAc pH 5.2 to sample to a final concentration of 0.3M

Note

This can be accomplished by diluting the sample by at least five volumes with 0.3M NaOAc or adding 3M NaOAc to a final concentration of 0.3M.

[CRITICAL] Do not exceed 500 ul NaOAc-diluted volume per 1.5 ml tube. For large volumes, see alternate protocol below.

3 Add 1 ul glycoblue per 250 ul sample volume

Note

volume here is the volume after adding NaOAc, not the original volume.

4 Vortex briefly

5 Add 2 volumes 100% ethanol

Note

volume here is the volume after NaOAc, not the original volume. For example, add 1 ml ethanol to 500 ul NaOAc solution.

6 Vortex briefly

7 Incubate at -80 degrees for at least 1 hour or in a dry ice ethanol bath for 30 minutes.

Note

[PAUSE] RNA is more stable for long periods when stored in ethanol at low temperature than as a pellet. Shorter incubation times are probably ok - glycobule protocol recommends 30 minutes at -20 C.

- 8 Pellet RNA at  16000 x g , in tabletop microfuge 16k g for 25 minutes at 4 °C

Note

Tip: orient the tubes so the hinge is pointed outwards and the pellet will be easier to identify along the hinge-side of the tube.

- 9 Aspirate supernatant carefully without dislodging pellet, first using a P-1000 for most volume and then a P-20 for the end. Using a vacuum aspirator is not recommended as it is easy to aspirate the pellet.

- 10 Gently add 500 ul of 70% ethanol at -20 °C

Note

the warmer the 70% ethanol is, the more RNA will re-solubilize. Take the 70% out from the -20 right before adding it and then return it to the -20.

- 11 Pellet RNA at  16000 x g , in tabletop microfuge 16k g for 5 minutes at 4 °C

- 12 Aspirate supernatant carefully without dislodging pellet, first using a P-200 for most volume and then a P-20 for the end.

- 13 Quick spin to bring down remaining ethanol

- 14 Remove as much ethanol as possible without aspirating pellet using a P-20

- 15 Dry for 5 minutes at RT and not longer - fully dried pellets are harder to resuspend



- 16 Add the desired volume of water or buffer
- 17 Scrape the tube along a tube rack to dislodge pellet
- 18 Quick spin
- 19 Measure concentration using nanodrop or qubit

Precipitation of RNA (large sample volume)

- 20 Add NaOAc to sample to a final concentration of 0.3M

Note

This can be accomplished by diluting the sample by at least five volumes with 0.3M NaOAc or adding 3M NaOAc to a final concentration of 0.3M.

- 21 Add 1 ul glycoblue per 250 ul sample volume (after adding NaOAc)
- 22 Vortex briefly
- 23 Add 2 volumes 100% ethanol (1 ml 100% ethanol to 500 ul NaOAc-diluted sample)
- 24 Vortex briefly
- 25 Incubate at -80 degrees 1 hour to overnight.

Note

[PAUSE] RNA is more stable for long periods when stored in ethanol at low temperature than as a pellet.



Pelleting and resuspension of RNA

26 Centrifuge in a conical tube at $> \sim 5000g$, for 25 minutes

27 Aspirate supernatant carefully with a serological pipet or vacuum aspirator

Note


don't worry about getting all the ethanol off - trying to get all the ethanol off will likely remove some pellet as this is a loose pellet

28 Add 750 μ l of 100% ethanol at -20°C

29 Resuspend chunks of pellet in the 100% ethanol by pipetting up and down

30 Transfer ethanol and chunks to a microfuge tube

31 Spin for 5 minutes at 16k g at 4°C

32  [go to step #12](#) to wash and resuspend pellet