

Apr 24, 2019

Euplotes crassus transfection through microinjection into the macronucleus

 Forked from [Euplotes crassus transfection through microinjection into the macronucleus](#)

 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.2argad6

Rachele Cesaroni¹

¹Universität Bern

Protist Research to Opti...



Angela Piersanti

University of Camerino



OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.2argad6

Protocol Citation: Rachele Cesaroni 2019. Euplotes crassus transfection through microinjection into the macronucleus. protocols.io <https://dx.doi.org/10.17504/protocols.io.2argad6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 24, 2019

Last Modified: April 24, 2019

Protocol Integer ID: 22577



- 1 Dilute 1:10 *Euplotes crassus* cultures of two different mating types in artificial sea water (prepare 20 ml culture for each mating type), and feed them with *E.coli* (3 ml for each mating type). For bacteria preparation see protocol 'Culturing *Euplotes crassus* to high densities using *E.coli* as the only food source'.

Note

Recipe for complete seawater (1 L):

36 g Reef Crystals
1 ml Walne's solution
1 ml of 10 µg/ml FeSO₄
0.2 ml of 2 mg/ml thiamine (light sensitive; store at 4°C)
Add distilled water up to 1 L

- 2 Grow cells in 250 ml flat-bottomed flasks at 24°C for 4 days with a 12h light/12h dark cycle, and then mix the same number of cells of both mating types in a 500 ml flat-bottomed flask at room temperature (the optimal cell density for conjugation is ~1000 cells/ml). Provide no aeration in both steps.
- 3 Isolate single *Euplotes crassus* cells with a donut shape after 2 days into artificial sea water with 2% BSA, in order to prepare drops for microinjection (ideally one cell each drop).
- 4 When drops are ready, cover them with a thin layer of Mineral Oil to not let them evaporate.
- 5 Inject into the macronucleus DNA (3 to 5 µg/µl) using Eppendorf Femtotips I injection needle.
- 6 Recover each cell individually in 500 µl of artificial sea water plus 0.25 µl of *E.coli* at 24°C (for bacteria preparation see protocol 'Culturing *Euplotes crassus* to high densities using *E.coli* as the only food source').