

Mar 30, 2017 Version 2

Euplotes crassus transfection through microinjection into the macronucleus V.2

DOI

dx.doi.org/10.17504/protocols.io.hanb2de



Rachele Cesaroni

Protist Research to Opti...



Rachele Cesaroni

Universität Bern

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.hanb2de

Protocol Citation: Rachele Cesaroni 2017. Euplotes crassus transfection through microinjection into the macronucleus. protocols.io <https://dx.doi.org/10.17504/protocols.io.hanb2de>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: March 03, 2017

Last Modified: December 06, 2017

Protocol Integer ID: 5166



- 1 Dilute 1:10 *Euplotes crassus* cultures of two different mating types in artificial sea water (20 ml total volume for each mating type) and feed them with *E. coli* (3 ml for each mating type). Before to add *E. coli* to the *Euplotes crassus* cells, pellet them and wash them once with ddH₂O (for bacteria preparation see protocol 'Culturing *Euplotes crassus* to high densities using *E. coli* as the only food source').

Note

Recipe for complete seawater (1 L):

36 g Reef Crystals
1 ml Walne's solution
1 ml of 10 µg/ml FeSO₄
0.2 ml of 2 mg/ml thiamine (light sensitive; store at 4°C)
Add distilled water up to 1 L

- 2 Grow cells at 24°C for 4 days with a 12h light/12h dark cycle even without aeration and then mix the two mating types at room temperature (the optimal cell density for conjugation is ~1000 cells/ml and the volume <2.5 cm in high).
- 3 Isolate single *Euplotes crassus* cells with a donut shape after 2 days into artificial sea water with 2% BSA in order to prepare drops for microinjection (ideally one cell each drop).

Note

Cells have a donut shape (anlagen) and macronucleus is easier to inject.
50h after formation of the pairs a round of amplification of the genome occurs.

- 4 When drops are ready cover them with a layer of Mineral Oil to not let them evaporate.
- 5 Inject into the macronucleus DNA in the concentration of 3 to 5 µg/µl using Eppendorf Femtotips I injection needle.
- 6 Recover each cell individually in 500 µl of artificial sea water plus 0.25 µl of *E. coli* at 24°C (for bacteria preparation see protocol 'Culturing *Euplotes crassus* to high densities using *E. coli* as the only food source').