



Jun 13, 2023

Ex vivo culture of SPMs

DOI

dx.doi.org/10.17504/protocols.io.8epv5jpdjl1b/v1

michela.deleidi^{1,2}, Bianca Marchetti^{3,4}, Federico Bertoli^{1,2}, Carmela Giachino³

¹Mitochondria and Inflammation in Neurodegenerative Diseases, DZNE, Tübingen-Germany;

²Hertie Institute for Clinical Brain Research, University of Tübingen;

³Neuropharmacology Laboratory, Oasi Research Institute-IRCCS, Troina, Italy;

⁴Biomedical and Biotechnological Sciences, Pharmacology Section, University of Catania-Italy



carmela giachino

institut imagine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5jpdjl1b/v1>

Protocol Citation: michela.deleidi, Bianca Marchetti, Federico Bertoli, Carmela Giachino 2023. Ex vivo culture of SPMs.
protocols.io <https://dx.doi.org/10.17504/protocols.io.8epv5jpdjl1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 31, 2023

Last Modified: June 13, 2023

Protocol Integer ID: 79805

Keywords: Ex vivo culture, spleen derived macrophages, derived macrophage, spm

Funders Acknowledgements:

ASAP-Aligning Science Across Parkinson's

Grant ID: ASAP-000420

Abstract

Ex vivo culture of spleen derived macrophages.

Attachments



[k4dpb48if.docx](#)

13KB

Materials

Materials

 RPMI 1640 **Coris bioConcept Catalog #1-41F01-I**

 DMEM High Glucose (45g/l) **Coris bioConcept Catalog #1-26F03-I**

- 10% FBS
- 2mM L-Glutamine
- Penicillin-Streptomycin Pen 10'000 IU/ml
- amphotericin B (BioConcept 4-01F00-H)
- microglia BV2 cells from Elabscience (No.: EP-CL-0493)

Troubleshooting

Ex vivo culture of spleen derived macrophages

- 1 Dissect spleens from abdominal cavity and filter it through a 40- μ m nylon strainer.
- 2 Use red cell lysis buffer to remove red cells.
- 3 Then, obtain a single splenic cell suspension. Culture cells in Roswell Park Memorial Institute (RPMI) medium 1640 RPMI 1640 (BioConcept 1-41F01-I) supplemented with 10% FBS, [IM] 2 millimolar (mM) L-Glutamine and antimicrobials (Penicillin-Streptomycin Pen 10'000 IU/ml Strep  10 mg/mL and amphotericin B ( 250 μ g/mL BioConcept).
- 4 Culture mouse microglia BV2 cells from Elabscience (No.: EP-CL-0493) in parallel for each spleen culture preparation as controls.
- 5 Briefly, maintain BV2 cells in Roswell Park Memorial Institute (RPMI) medium 1640 supplemented (BioConcept 1-41F01-I) with 10% FBS (FBS-02-0500), [IM] 2 millimolar (mM) L-Glutamine 5-10K50-H) and antimicrobials (Penicillin-Streptomycin Pen 10'000 IU/ml Strep  10 mg/mL and amphotericin B ( 250 μ g/mL) (BioConcept 4-01F00-H).
- 6 Spleen macrophages (SPMs) differentiate into the M1 phenotype after stimulation with LPS (100 ng/ml) $\pm$ IFN- γ ( 10 ng/mL).
- 7 For BV2 stimulation, replace RPMI by Dulbecco's Modified Eagle Medium (DMEM) High Glucose. (BioConcept, 1-26F03-I).