



Oct 07, 2022

Ex Vivo Electrophysiology

DOI

dx.doi.org/10.17504/protocols.io.81wgby7zovpk/v1

Harry Xenias¹, Savio Chan¹, Loukia Parisiadou¹

¹Northwestern University



divya.darwinarulseeli

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.81wgby7zovpk/v1>

Protocol Citation: Harry Xenias, Savio Chan, Loukia Parisiadou 2022. Ex Vivo Electrophysiology. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.81wgby7zovpk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 07, 2022



Last Modified: May 31, 2024

Protocol Integer ID: 71033

Keywords: ASAPCRN, ex vivo electrophysiology this protocol, ex vivo electrophysiology, preparation of brain slice, setup of the electrophysiology rig, electrophysiology rig, cell recording, brain slice, vivo, preparation, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's through the Michael J. Fox Foundation for Parkinson's Research (MJFF)

Grant ID: [ASAP-020600]

Abstract

This protocol describes the preparation of brain slices, setup of the electrophysiology rig, and solutions for collecting whole-cell recordings.

Troubleshooting

Setup

- 1 First prepare a 10x stock aCSF solution by first add about 200 mL of ddH₂O water into a clean 2L flask. Then add each of the following powders:
§NaCl: **1250 mM**
§KCl: **25 mM**
§NaHCO₃: **250 mM**
§NaH₂PO₄: **12.5 mM**
- 2 Next, fill the flask approximately three quarters with ddH₂O and using a magnetic stirrer, allow the solution to be thoroughly dissolved. The mixed solution should then be brought to a final 2L volume and stored in a 2L glass bottle and refrigerated at 4°C.

Note: the 10x stock aCSF can be used for up to a month and should be remade fresh after that time.
- 3 Prepare 1 liter of 1X aCSF (to be used as both perfusate and cutting solution):
 - Fill a 2L flask with about **450 mL** of ddH₂O.
 - Add **4 mL** of CaCl₂ solution to the bottle.
 - Add **4 mL** of MgCl₂ solution to the bottle.
 - Add **5.02 g** of glucose (final concentration of 13.93 mM)
 - Add **200 mL** of a 10X stock aCSF solution.
 - Swirl mixture around well by hand for a few seconds.
 - Fill remainder of flask to **2 L** with ddH₂O.
 - Carbogenate solution for at least 10 minutes before using.
- 4 Using continuously carbogenated 1x aCSF: §Fill a small slice holding chamber with the 1x aCSF and add **5 mM** of L-glutathione at approximately **1:1000x** and **1 mM** of Na-Pyruvate at **1:300x**.
 - Set the chamber aside near a **37 °C** heated water bath while the solution inside the holding chamber is continuously carbogenated.
 - Additionally, decant approximately **100 mL** of 1X aCSF into a small glass beaker kept cold on ice while carbogenated (to be use for perfusate)
- 5 Next, decant approximately **100 mL** of the 1X solution into a **300 mL** L plastic bucket and place the bucket in a **-20 °C** freezer for **70—80 minutes**. This solution is to be used for cutting brain tissues and should be frozen over but not frozen solid. After taking the bucket out from freezer:
 - Using a large spatula, break up ice and stir into a slurry.
 - Add approximately **100 mL** of 1X aCSF kept at room temperature.
 - Mix solution with a handheld blender until forming an easily flowing slush.
 - Keep the bucket of the 1x aCSF slush cutting solution on ice while being carbogenated.

-The slush solution should consist of approximately one fifths of liquid solution and be settled to the bottom of the bucket.

- 6 Overdose the mouse with a **1 mL** intraperitoneally injection of a **150 mg/kg ketamine and 30 mg/kg xylazine** mix. While mouse is overdosing: place the slice chamber in a heated water bath set to **37 °C** (water in the bath should come up the side of the chamber to approximately the same height as the aCSF). Next, after mouse is completely anesthetized:

- Transcardially perfuse mouse with the ice-cold 1x aCSF.
- Rapidly decapitate and extract brain in 1x aCSF slush.
- Using razor blade, cut brain down the midline.
- Glue brain medial side down on cutting block and place block quickly into cutting chamber of a Vibratome.
- Fill cutting chamber with aCSF slush solution, kept continuously carbogenated in chamber.
- Cut slices at **240 µm** thick
- Quickly transfer each cut slice into hold chamber kept in heated bath.
- Remove holding chamber from bath 30 minutes after last slice transfer is made.
- Allow holding chamber to equilibrate to room temperature (approximately 20 minutes) before transferring slices to recording chamber of electrophysiological rig.

Whole-cell Recordings

- 7 Using a pipette puller, glass pipettes should be pulled to ensure a pipette resistance of **3.2—3.8 MΩ**. For whole-cell voltage clamp recordings, a KMeSO₄ solution is used containing the following:

- KMeSO₄: **135 mM**
- KCl: **5 mM**
- CaCl₂: **0.5 mM**
- HEPES: **5 mM**
- EGTA-K: **5 mM**
- ATP-Mg: **2 mM**
- GTP-Na: **0.5 mM**
- Biocytin: 0.20% (w/v in grams)

- 8 Transferred slices in recording chamber are continuously perfused with 1x aCSF that is kept at room temperature and continuously carbogenated. Neurons for recording are identified and recorded as follows:

- Pipettes are backfilled with KMeSO₄ solution and inserted into headstage of amplifier.
- Pipettes are then pressurized to approximately **56 millibars**.
- Pipettes area offset before cell attachment.
- Neurons are clamped at **−80 mV** before whole-cell access is achieved.

- 9 Stimulus generation and data acquisition are performed using an amplifier (Molecular Devices), a digitizer (Molecular Devices), and pClamp (Molecular Devices).
- 10 For current-clamp recordings
Adjust the amplifier bridge circuit to compensate for electrode resistance and subsequently monitor it. Filter the signals at 1 kHz and digitize them at 10 kHz.
KMeSO₄ and Na₂-GTP were from ICN Biomedicals and Roche, respectively. All other reagents were obtained from Sigma-Aldrich.

Excitability of SPNs:

- 11 Examine the frequency-current (F-I) relationship of each cell with current-clamp recordings as follows.
- Apply a series of 500 ms current steps of n beginning at -150 pA and incremented at 25 pA for each consecutive sweep.
 - A **30 second** intertrial interval is used.
 - The current steps are continued until depolarization block is reached.
 - Monitor resting membrane potential was monitored for stability, and exclude cells that varied 20% from mean baseline from the analysis.

Corticoatrial responses recorded in voltage-clamp

- 12
- Perform electrical stimulation using parallel bipolar tungsten electrodes (FHC) placed in layer 5 of the cortex.
 - Adjust stimulus width and intensity via a constant current stimulator (Digitimer) to evoke a first excitatory postsynaptic current (EPSC) with an amplitude of 200–400 pA in the presence of the GABA_A receptor antagonist SR95531 (10 μM) and CGP55845 (1 μM).
 - Monitor whole-cell access was monitored with a -5 mV pulse throughout the recording. Determine off line the membrane capacitance (C_m) as $C_m = Q_t * V_{test}$, where Q_t was calculated as the integral of the transient current elicited by V_{test}, a 10 mV voltage step.
 - Take the average of the ratios of the second EPSC amplitude to the first EPSC amplitude for each recording sweep to calculate the paired-pulse ratio (PPR) for a given cell
 - Exclude data if the series resistance of the patch pipette differed by > 20% between the two recordings.