



Oct 16, 2023

Expression and purification for cryo-EM samples of FIP200NTD:ATG13(363-517)-ULK1(MIT) complexes

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwjxw7lmk/v1>

Minghao Chen¹, Xuefeng Ren¹

¹Department of Molecular and Cell Biology, University of California, Berkeley, Berkeley, CA 94720, USA



Minghao Chen

University of California, Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwjxw7lmk/v1>

Protocol Citation: Minghao Chen, Xuefeng Ren 2023. Expression and purification for cryo-EM samples of FIP200NTD:ATG13(363-517)-ULK1(MIT) complexes. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwjxw7lmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 05, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82914

Keywords: ASAPCRN, em samples of fip200ntd, purification for cryo, fip200ntd, purification, cryo

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-000350

Abstract

Expression and purification for cryo-EM samples of FIP200NTD:ATG13(363-517)-ULK1(MIT) complexes



Expression

2d 0h 40m

- 1 Transfect HEK GNTi cells at concentration of 2×10^6 cells/ml
- 2 Dilute PEI with Warm Hybridoma-SFM(1X)
- 3 In a separate tube, dilute DNA with Hybridoma-SFM(1X)
- 4 Add PEI to DNA dilution. Incubate mixture for 00:30:00 at 37 °C
- 5 Add mixture to cells. Let cells grow for 48:00:00
- 6 Harvest Cells 500 rpm , 4 °C , 00:10:00
- 7 Wash pellet with cold PBS. Store pellet at -80 °C until purification.

30m

2d

10m

Purification

2h 50m

- 8 Cell pellets were lysed at room temperature for 20 min with lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 1 mM TCEP, 10% Glycerol) with 5 mM EDTA, 1% Triton X-100 and protease inhibitor cocktail (Thermo Scientific) for 00:15:00
- 9 Clarify lysate for 17000 rpm for 00:35:00 at 4 °C .
- 10 Wash strep-tactin resin (IBA Lifesciences, Germany) into lysis buffer (without Triton). Load clarified lysate onto resin

15m

35m



- 11 Rock supernatant with equilibrated Strep-Tactin Sepharose resin for 01:00:00 at 4 °C 1h
- 12 Wash with 5CV lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 1 mM TCEP, 5 mM EDTA, 10% Glycerol)
- 13 Elute with lysis buffer plus 4 mM desthiobiotin for STREP resin
- 14 His6-TEV cleavage at 4 °C Overnight , 1h
- 15 Concentrate elution and inject onto pre-equilibrated Superose 6 Increase 10/300 GL column (Cytiva)
(25 mM HEPES pH 7.5, 150 mM NaCl, 1 mM MgCl₂, 1 mM TCEP)
- 16 Pool peak fractions, concentrate, snap freeze, and store at -80 °C