

Oct 16, 2023

Expression and purification of PI3KC3-C1 complex

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvorjezv4o/v1>

Minghao Chen¹, Xuefeng Ren¹

¹Department of Molecular and Cell Biology, University of California, Berkeley, Berkeley, CA 94720, USA



Minghao Chen

University of California, Berkeley

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvorjezv4o/v1>

Protocol Citation: Minghao Chen, Xuefeng Ren 2023. Expression and purification of PI3KC3-C1 complex. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvorjezv4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 05, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 82916

Keywords: ASAPCRN, c1 complex expression, complex expression, purification, c1

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-000350

Abstract

Expression and purification of PI3KC3-C1 complex



Expression

2d 0h 40m

- 1 Transfect HEK GNTi cells at concentration of 2×10^6 cells/ml
- 2 Dilute PEI with Warm Hybridoma-SFM(1X)
- 3 In a separate tube, dilute DNA with Hybridoma-SFM(1X)
- 4 Add PEI to DNA dilution. Incubate mixture for  00:30:00 at  37 °C
- 5 Add mixture to cells. Let cells grow for  48:00:00
- 6 Harvest Cells  500 rpm ,  4 °C ,  00:10:00
- 7 Wash pellet with cold PBS. Store pellet at  -80 °C until purification.

30m

2d

10m

Purification

2h 50m

- 8 Pellets were homogenized 20 times by Pyrex douncer (Corning) in lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 10% Glycerol) with 25mM TCEP/proteinase inhibitors (Thermo Scientific), and then add 10% Triton X-100 stock to final 1% concentration.
- 9 Rocking at  4 °C for  01:00:00
- 10 Clarify lysate for  17000 rpm for  00:40:00 at  4 °C .

1h

40m



- 11 Wash strep-tactin resin (IBA Lifesciences, Germany) into lysis buffer (without Triton). Load clarified lysate onto resin
- 12 Rock supernatant with equilibrated Strep-Tactin Sepharose resin for  Overnight at  4 °C
- 13 Wash with 5CV lysis buffer (25 mM HEPES pH 7.5, 200 mM NaCl, 2 mM MgCl₂, 25 mM TCEP, 10% Glycerol)
- 14 Elute with lysis buffer plus 4 mM desthiobiotin for STREP resin
- 15 Concentrate elution and inject onto pre-equilibrated Superose 6 Increase 10/300 GL column (Cytiva)
(25 mM HEPES pH 7.5, 300 mM NaCl, 1 mM MgCl₂, 1 mM TCEP)
- 16 Pool peak fractions, concentrate, snap freeze, and store at  -80 °C

1h