

Jan 27, 2023

Expression and Purification of recombinant Bst DNA polymerase (Bst)

DOI

dx.doi.org/10.17504/protocols.io.dm6gpb2e8lzp/v1

Diana A Tapia-Sidas¹, Brenda Vargas-Hernández¹, José Abrahán Ramírez-Pool¹, Leandro A Nuñez-Muñoz¹, Berenice Calderón-Pérez¹, Rogelio González-González¹, Luis Gabriel Briebe², Rosalía Lira-Carmona³, Eduardo Ferat-Osorio⁴, Constantino López-Macías⁴, Roberto Ruiz-Medrano¹, Beatriz Xoconostle-Cázares¹

¹Departamento de Biotecnología y Bioingeniería, Centro de Investigación y de Estudios Avanzados, Mexico City, Mexico;

²Laboratorio Nacional de Genómica para la Biodiversidad, Centro de Investigación y de Estudios Avanzados, Irapuato, Guanajuato, Mexico;

³Unidad de Investigación Médica en Enfermedades Infecciosas y Parasitarias, Hospital de Pediatría, Centro Médico Nacional Siglo XXI, Instituto Mexicano del Seguro Social, Mexico City, Mexico;

⁴División de Investigación en Salud, UMAE Hospital de Especialidades "Dr. Bernardo Sepúlveda Gutiérrez", Centro Médico Nacional Siglo XXI, Instituto Mexicano del Seguro Social, Mexico City, Mexico



Diana A Tapia-Sidas

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpb2e8lzp/v1>

External link: <https://doi.org/10.1371/journal.pone.0279681>

Protocol Citation: Diana A Tapia-Sidas, Brenda Vargas-Hernández, José Abrahán Ramírez-Pool, Leandro A Nuñez-Muñoz, Berenice Calderón-Pérez, Rogelio González-González, Luis Gabriel Briebe, Rosalía Lira-Carmona, Eduardo Ferat-Osorio, Constantino López-Macías, Roberto Ruiz-Medrano, Beatriz Xoconostle-Cázares 2023. Expression and Purification of recombinant Bst DNA polymerase (Bst). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpb2e8lzp/v1>

Manuscript citation:

D. Tapia-Sidas, B. Vargas-Hernández, J. Ramírez-Pool, L. Nuñez-Muñoz, B. Calderón-Pérez, R. González-González, L. Briebe, R. Lira-Carmona, E. Ferat-Osorio, R. Ruiz-Medrano, B. Xoconostle-Cázares. (2022). Starting from scratch: step-by-step development of a diagnostic test for SARS-CoV-2 detection by RT-LAMP. PLoS One.

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: May 27, 2022

Last Modified: January 27, 2023

Protocol Integer ID: 63361

Keywords: Bacillus stearothermophilus, polymerase, Bst, protein expression, protein purification, purification of recombinant bst dna polymerase, recombinant bst dna polymerase, bst polymerase, bst enzyme, isothermal nucleic acid amplification, gold standard in isothermal nucleic acid amplification, dna polymerase, mediated isothermal amplification, thermostable reverse transcriptase, bst, geobacillus stearothermophilus, isothermal amplification, bacillus steothermophilus, useful for rna virus, rna template, purification procedure by fplc, implementation in diagnostic technique, diagnostic technique, purification procedure, mass production of diagnostic test, nucleic acid, rna virus, diagnostic test, transcript detection, rna

Funders Acknowledgements:

Asociación Mexicana de Cooperación para el Desarrollo (AMEXCID) - Secretaría de Relaciones Exteriores (SRE), Mexico

Grant ID: AMEXCID/2020-1

Consejo Nacional de Ciencia y Tecnología (CONACyT)

Grant ID: D.A.T-S (754736) and J.A.R.-P. (483659)

Abstract

Bst is a type I DNA polymerase with strong strand displacement activity isolated from *Geobacillus stearothermophilus* (previously *Bacillus steothermophilus*). Bst enzyme is a gold standard in isothermal nucleic acid amplification diagnostic techniques, especially in Loop-mediated isothermal amplification (LAMP). LAMP is a low-cost technique, provides a visual detection (when combined with pH indicators) and does not require the use of thermal cyclers. Also, Bst combined with thermostable reverse transcriptase can amplify RNA templates, in a technique known as RT-LAMP. RT-LAMP is useful for RNA virus and transcript detection and can be employed in circumstances that require mass production of diagnostic tests or limited availability of resources. This protocol shows the expression and purification procedure by FPLC of the Bst polymerase for its implementation in diagnostic techniques such as end-point colorimetric and real time fluorometric LAMP and RT-LAMP.



Guidelines

During the process of protein purification maintain all samples that contain the protein of interest in a cold environment to avoid protein degradation.



Materials

Reagents:

-  Chaperone Plasmid Set **Takara Bio Inc. Catalog #3340**
-  Chemically Competent E. coli One Shot™ BL21(DE3) **Invitrogen - Thermo Fisher Catalog #C600003**
- [M] 100 ng/μL pColdI-Bst plasmid Step 3.1
- Stock [M] 30 mg/mL chloramphenicol
- Stock [M] 100 mg/mL carbenicillin
- Stock [M] 1 Molarity (M) IPTG
- Stock [M] 10 % (v/v) Triton X-100
-  250 μL of SOC medium Step 1.4
-  8 mL of TFB I solution Step 2.4
-  2.5 mL of TFB II solution Step 2.5
- Cryotubes with  500 μL of 30% v/v glycerol (sterilized)
- LB agar plates
- Tubes with  3 mL Luria-Bertani (LB) medium
- Tubes with  5 mL Luria-Bertani (LB) medium
- Tubes with  5 mL Terrific Broth (TB) medium
- Flasks with  20 mL LB medium
- Flasks with  50 mL LB medium
- Flasks with  50 mL TB medium
- Flask with  100 mL LB medium
- Flasks with  1 L LB medium
-  250 mL lysis buffer A (LB-A) Step 5.5
-  250 mL Saline Buffer (SB) Step 9.2
-  500 mL Washing Buffer-AI (WB-AI) Step 10.2
-  500 mL Elution Buffer-AI (EB-AI) Step 10.5
-  1 L Desalting Buffer-A (DB-A) Step 11.2
-  500 mL Washing Buffer-AII (WB-AII) Step 12.2
-  500 mL Elution Buffer-AII (EB-AII) Step 12.5
-  500 mL Storage Buffer-A (SB-A) Step 13.1
-  HisTrap HP 5mL **Cytiva Catalog #17524801**

-  HiPrep 26/10 Desalting Column **Cytiva Catalog #17508701**
-  HiTrap Heparin HP affinity column **Cytiva Catalog #17040701**
-  Quick Start™ Bradford 1x Dye Reagent **Bio-Rad Laboratories Catalog #5000205**
-  Precision Plus Protein™ Unstained Protein Standards **Bio-Rad Laboratories Catalog #1610363**
- Tricine-SDS-PAGE electrophoresis solutions (Step 6.2)
- 8% polyacrylamide gels for Tricine-SDS-PAGE (Step 6.2)

Equipments:

- Thermomixer

Equipment

Thermomixer® R

NAME

Dry block heating and cooling shaker, 120 V, 60 Hz, 1/cs

TYPE

Eppendorf

BRAND

T3317

SKU

- Orbital shaker

Equipment

MaxQ™ HP Incubated Tabletop Orbital Shaker

NAME

MaxQ™ HP, 120 V 60 Hz, 6,5 A o 230 V 50/60 Hz, 3,2 A

TYPE

Thermo Scientific

BRAND

SHKE420HP

SKU

<https://www.thermofisher.com/order/catalog/product/SHKE420HP>

LINK

- Floor model orbital shaker

Equipment

MaxQ™ HP Incubated and Refrigerated Console Shakers NAME

MaxQ™ 481 HP, 230 V, 50 Hz TYPE

Thermo Scientific BRAND

SHKE481HP SKU

<https://www.thermofisher.com/order/catalog/product/SHKE481HP> LINK

- Centrifuge

Equipment

Sorvall™ Legend™ XT/XF Centrifuge Series NAME

Thermo Scientific BRAND

75004541 SKU

<https://www.thermofisher.com/order/catalog/product/75004541> LINK

- Ultrasonic Processor 130W

Equipment

Ultrasonic Processor

NAME

130-Watt Ultrasonic Processor

TYPE

Cole-Parmer

BRAND

ML-04714-52

SKU

<https://www.coleparmer.com/p/cole-parmer-130-watt-ultrasonic-processors/44347>^{LINK}

- Ultrasonic Processor 750W

Equipment

750-Watt Ultrasonic Processor

NAME

CPX750

TYPE

Cole-Parmer

BRAND

ML-04711-60

SKU

<https://www.coleparmer.com/p/cole-parmer-500-and-750-watt-ultrasonic-processors/16401>^{LINK}

-Nanodrop

Equipment

NanoDrop™ One UV-Vis Spectrophotometer

NAME

spectrophotometer

TYPE

Thermo Scientific

BRAND

ND-ONE-W

SKU

<https://www.thermofisher.com/order/catalog/product/ND-ONE-W>

LINK

Sample Volume (Metric): Minimum 1µL; Spectral Bandwidth: ≤1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm

SPECIFICATION

- FPLC system

Equipment

ÄKTA pure

NAME

Protein purification system

TYPE

Cytiva

BRAND

29046665

SKU

<https://www.cytivalifesciences.com/en/us/support/products/akta-pure-150-l-29046665>

LINK

- Spectrophotometer UV/Vis
- Incubator (37°C)
- Water bath (60°C)
- Ultra Low-Temperature Freezer (-80°C)
- Freezer -20°C
- Refrigerator (4°C)
- Analytical balance

Other:

- Ice bath
- Microcentrifuge tubes
- Sterile 0.45 µm membrane filter
- 150 mL Superloop (Cytiva)
- Dialysis membrane
- Ultrafiltration tube (Amicon Ultra-15)

Equipment	
Amicon Ultra-15	NAME
PLTK Ultracel-PL membrane, 15 ML - 30 kDa cutoff	TYPE
Millipore	BRAND
UFC903024	SKU
https://www.merckmillipore.com/MX/es/product/Amicon-Ultra-15-Centrifugal-Filter-Unit,MM_NF-UFC903024?ReferrerURL=https%3A%2F%2Fwww.bing.com%2F&bd=1	LINK

- Image Lab 6.1 Software (Bio-Rad)



Protocol materials

⊗ Chaperone Plasmid Set **Takara Bio Inc. Catalog #3340**

⊗ HiPrep 26/10 Desalting Column **Cytiva Catalog #17508701**

⊗ HisTrap HP 5mL **Cytiva Catalog #17524801**

⊗ Chemically Competent E. coli One Shot™ BL21(DE3) **Invitrogen - Thermo Fisher Catalog #C600003**

⊗ Quick Start™ Bradford 1x Dye Reagent **Bio-Rad Laboratories Catalog #5000205**

⊗ HiTrap Heparin HP affinity column **Cytiva Catalog #17040701**

⊗ Precision Plus Protein™ Unstained Protein Standards **Bio-Rad Laboratories Catalog #1610363**

⊗ Chemically Competent E. coli One Shot™ BL21(DE3) **Invitrogen - Thermo Fisher Catalog #C600003**

⊗ Chaperone Plasmid Set **Takara Bio Inc. Catalog #3340**

⊗ HisTrap HP 5mL **Cytiva Catalog #17524801**

⊗ HiPrep 26/10 Desalting Column **Cytiva Catalog #17508701**

⊗ Quick Start™ Bradford 1x Dye Reagent **Bio-Rad Laboratories Catalog #5000205**

⊗ HiTrap Heparin HP affinity column **Cytiva Catalog #17040701**

⊗ Precision Plus Protein™ Unstained Protein Standards **Bio-Rad Laboratories Catalog #1610363**

Troubleshooting

Before start

Ensure to have all the necessary materials and reagents already cleaned, sterilized and filter (in case of the purification solutions).



Preparation of Bst expression cells

20h 45m

1 Transformation of chemically competent BL21 (DE3) cells with pKJE7 plasmid.

1.1 Add 1 μ L of plasmidic DNA consisting of the pKJE7 plasmid from the

20m

Chaperone Plasmid Set Takara Bio Inc. Catalog #3340 to

50 μ L of competent cells

Chemically Competent E. coli One Shot™ BL21(DE3) Invitrogen - Thermo
Fisher Catalog #C600003

. Mix the cells gently and incubate On ice for 00:20:00 .

1.2 Transfer the cells to a heat block at 42 °C and incubate for 00:00:53 .

53s

1.3 Transfer the cells immediately to an ice bath and incubate On ice for

5m

00:05:00 .

1.4 Add 250 μ L of SOC medium at room temperature to the transformed cells and
incubate at 225 rpm, 37°C, 01:00:00 .

1h

**Note****SOC medium composition**

	A	B
	Tryptone	2%
	Yeast extract	0.5%
	NaCl	10 mM
	KCl	2.5 mM
	MgCl ₂	10 mM
	MgSO ₄	10 mM
	Glucose	20 mM

Adjust to pH 7 and sterilize by filtration.

- 1.5 Plate 25 µL of transformed cell culture onto LB agar with the corresponding selective agent. Incubate the plates Overnight at 37 °C .

18h

Note

The pKJE7 plasmid requires 30 µg/mL chloramphenicol as selective agent.

- 1.6 Select a single colony of transformed cells and inoculate in 3 mL Luria-Bertani (LB) medium supplemented with the selective antibiotic. Incubate Overnight at 180 rpm, 37°C .

18h

- 1.7 Centrifuge the cell culture at 10000 x g, 4°C, 00:05:00 . Remove the supernatant and resuspend the cell pellet in 500 µL LB medium .

5m



1.8 Add 500 μ L of 30% v/v glycerol , mix by pipetting up and down and store at -80 $^{\circ}$ C .

5m

2 Preparation of chemically competent BL21 (DE3) cells harboring pKJE7 plasmid.

2.1 Take BL21(DE3) cells harboring pKJE7 plasmid from a frozen glycerol stock using a bacterial inoculating loop and inoculate 3 mL LB liquid medium with 30 μ g/mL of chloramphenicol . Incubate Overnight at 180 rpm, 37 $^{\circ}$ C

18h

2.2 Inoculate 1 mL overnight culture in 100 mL LB medium with 30 μ g/mL chloramphenicol and incubate at 180 rpm, 37 $^{\circ}$ C, 03:00:00 .

3h

Note

Monitor the cell growth by measuring the optical density (**OD**) at **600 nm** and remove the cells from incubation when the OD reaches **0.3 to 0.4**.

2.3 Chill the cell culture On ice for 00:10:00 and centrifuge the cells at 4000 x g, 4 $^{\circ}$ C, 00:10:00 .

20m

2.4 Gently resuspend the cell pellet in 8 mL of TFBI solution pre-chilled and incubate On ice for 00:45:00 . Centrifuge the cells at 4000 x g, 4 $^{\circ}$ C, 00:05:00 .

50m

**Note****TFBI Solution Composition**

	A	B
	Potassium acetate	30 mM
	Rubidium chloride	100 mM
	Calcium chloride	10 mM
	Manganese chloride	50 mM
	Glycerol	15% v/v

Adjust to pH 5.8 with 1M acetic acid and sterilize by filtration.

- 2.5 Gently resuspend the cell pellet in  2.5 mL of TFBII solution pre-chilled and incubate  On ice for  00:05:00 .

5m

Note**TFBII Solution Composition**

	A	B
	MOPS	10 mM
	Rubidium chloride	10 mM
	Calcium chloride	75 mM
	Glycerol	15% v/v

Adjust to pH 6.5 with 1M sodium hydroxide and sterilize by filtration.



- 2.6 Prepare aliquots of  50 μ L of competent cells using microcentrifuge tubes previously chilled on an ice bath. Place the aliquots on dry ice until frozen and store at  -80°C .

30m

3 Transformation of chemically competent BL21 (DE3)/pKJE7 cells with the pColdI-Bst plasmid, the expression vector for the large fragment of the *Bacillus stearotherophilus* DNA polymerase (Bst).

- 3.1 Add  1 μ L of plasmidic DNA of  100 ng/ μ L pColdI-Bst expression vector to  50 μ L of competent cells BL21 (DE3)/pKJE7. For the transformation procedure  [go to step #1](#).

Note

The pColdI-Bst plasmid require  100 μ g/mL carbenicillin as selective agent and the pKJE7 plasmid require  30 μ g/mL chloramphenicol. Use LB medium supplemented with both antibiotics as selective media.

Small-scale screening cultures

4 Preparation of bacterial cultures for Bst expression.

- 4.1 Inoculate  5 μ L glycerol stock of **BL21(DE3)/pKJE7/pColdI-Bst** or **BL21(DE3)/pKJE7** cells in  5 mL culture medium (LB or TB) supplemented with selection agent(s). Incubate  Overnight at  200 rpm, 37°C .
- 4.2 Inoculate  500 μ L overnight culture in  50 mL culture medium with antibiotic(s). Use LB or TB according to the medium used for the overnight culture. Incubate  200 rpm, 37°C for approximately  03:00:00.

18h

3h

Note

Monitor the cell growth by measuring the optical density **(OD) at 600 nm** and remove the cells from incubation when the OD is between **0.6**.

- 4.3 Once the culture reaches an OD600 of 0.6, incubate the cell cultures  On ice for  00:20:00 before adding the inducer (IPTG).

20m

5 Small-scale Bst expression under different induction conditions.

- 5.1 Induce the expression of the Bst under different conditions. Each treatment should be evaluated in triplicate. For example:

30m

	Strain	[IPTG]	Temperature	Medium
	BL21(DE3)/pKJE7*	0.5 mM	16°C	LB
	BL21(DE3)/pKJE7*	0.5 mM	37°C	LB
	BL21(DE3)/pKJE7*	0.5 mM	16°C	TB
	BL21(DE3)/pKJE7/pColdl-Bst	0 mM	16°C	LB
	BL21(DE3)/pKJE7/pColdl-Bst	0.1 mM	16°C	LB
	BL21(DE3)/pKJE7/pColdl-Bst	0.5 mM	16°C	LB
	BL21(DE3)/pKJE7/pColdl-Bst	1.0 mM	16°C	LB
	BL21(DE3)/pKJE7/pColdl-Bst	0 mM	37°C	LB
	BL21(DE3)/pKJE7/pColdl-Bst	0.5 mM	37°C	LB
	BL21(DE3)/pKJE7/pColdl-Bst	0.5 mM	16°C	TB

* BL21(DE3)pKJE7 strain is used as negative expression control.



5.2 Incubate at 16 °C or 37 °C according to each treatment at 180 rpm for 16:00:00 . 16h

5.3 Centrifugate the cell cultures at 6000 x g, 4°C, 00:10:00 . Discard the supernatant, remove all the liquid and leave the cell pellet as dry as posible. 10m

5.4 Weigh the centrifugation tube with the cell pellet (total weight). 10m

Note

Weigh the empty tube prior centrifugation and subtract it to the total weight to calculate the weight of the cell pellet and hence the biomass produced.

5.5 Resuspend the cell pellet in 5 mL lysis buffer A (LB-A) (pre-chilled). 5m

Note

Lysis buffer A (LB-A) composition

	A	B
	Tris-HCl pH 7.5	50 mM
	EDTA	0.5 mM
	2-mercaptoethanol	10 mM
	Tergitol NP-40	0.1% v/v
	Tween-20	0.1% v/v
	PMSF	3 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

5.6 Disrupt cells by ultrasonication at an amplitude of 40%. Apply five cycles of 00:00:15 on and 00:00:30 off. 3m 45s



Note

Place the tubes  On ice while processing.

Equipment

Ultrasonic Processor

NAME

130-Watt Ultrasonic Processor

TYPE

Cole-Parmer

BRAND

ML-04714-52

SKU

<https://www.coleparmer.com/p/cole-parmer-130-watt-ultrasonic-processors/44347>

LINK

5.7 Centrifugate at  6000 x g, 4°C, 00:15:00 . Recover the supernatant (soluble protein fraction) and discard the pellet.

15m

5.8 For protein clarification, incubate the supernatant at  60 °C for  00:20:00 in a water bath and centrifugate at  14500 x g, 4°C, 00:15:00 . Recover the clarified supernatant and place it  On ice .

35m

6 Analysis of Bst expression.

6.1 Measure total protein concentration by measuring absorbance at 280 nm in a NanoDrop spectrophotometer.

3m



Equipment

NanoDrop™ One UV-Vis Spectrophotometer

NAME

spectrophotometer

TYPE

Thermo Scientific

BRAND

ND-ONE-W

SKU

<https://www.thermofisher.com/order/catalog/product/ND-ONE-W>

LINK

Sample Volume (Metric): Minimum 1µL; Spectral Bandwidth: ≤1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm

SPECIFICATIONS

- 6.2 Analyze all clarified supernatant samples by Tricine-SDS-PAGE electrophoresis through a 8% polyacrylamide gel. Load  100 µg protein sample per well.

4h

Citation

Hermann Schägger (2006). Tricine-SDS-PAGE. Nature Protocols.

[10.1038/nprot.2006.4](https://doi.org/10.1038/nprot.2006.4)

LINK

- 6.3 Select the best conditions for protein expression according to the results analysis (biomass, total protein production and electrophoretic profile).

Large-scale production of Bst

2h 6m 20s

7 Expression of recombinant Bst.



7.1 Inoculate  10 μ L glycerol stock of **BL21(DE3)/pKJE7/pColdI-Bst** expression cells in  20 mL LB medium supplemented with  100 μ g/mL carbenicillin and  30 μ g/mL chloramphenicol. Incubate  Overnight at  180 rpm, 37°C.

18h

7.2 Inoculate  10 mL overnight culture in  1 L LB medium with  100 μ g/mL carbenicillin and  30 μ g/mL chloramphenicol. Incubate at  200-220 rpm, 37°C, 03:00:00.

3h

Note

Monitor the cell growth by measuring the optical density (**OD**) at **600 nm** and remove the cells from incubation when the OD reaches **0.6**.

7.3 Place the inoculum  On ice for  00:30:00 and then add  0.5 millimolar (mM) IPTG for induction.

30m

Note

Do not add any additional inducers. For expression of the chaperones contained in pKJE7 plasmid, the basal expression is enough to promote correct Bst enzyme folding.

7.4 Incubate at  180 rpm, 16°C, 16:00:00 for recombinant protein expression.

16h

8 Soluble protein fraction recovery and clarification.

8.1 Centrifugate at  6000 x g, 4°C, 00:12:00 to harvest cells. Discard the supernatant ensuring to remove all the liquid and leave the cell pellet as dry as possible.

12m

8.2 Weigh the centrifugation tube with the cell pellet (total weight).

3m

**Note**

Weigh the empty tube prior centrifugation and subtract it to the total weight to calculate the weight of the cell pellet and hence the biomass produced.

8.3 Store the cell pellet at -80 °C until use (just in case that the purification step is not performed immediately after expression).

8.4 Resuspend the cell pellet in 50 mL lysis buffer A (LB-A) (pre-cooled). If necessary, defroze the cell pellet in an ice bath before adding the lysis buffer.

5m

Note

Lysis buffer A composition (LB-A)

	A	B
	Tris-HCl pH 7.5	50 mM
	EDTA	0.5 mM
	2-mercaptoethanol	10 mM
	Tergitol NP-40	0.1% v/v
	Tween-20	0.1% v/v
	PMSF	3 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

8.5 Disrupt cells by ultrasonication with an ultrasonic processor at an amplitude of 40% applying pulses of 00:00:10 of ultrasonication and 00:00:10 of pause during 00:04:00 .

4m 20s



Note

Place the sample  On ice and keep it cold while processing.

Equipment

750-Watt Ultrasonic Processor

NAME

CPX750

TYPE

Cole-Parmer

BRAND

ML-04711-60

SKU

<https://www.coleparmer.com/p/cole-parmer-500-and-750-watt-ultrasonic-processors/16401>

LINK

8.6 Centrifugate at  11000 x g, 4°C, 00:30:00 . Recover the supernatant (soluble protein fraction) and discard the pellet.

30m

8.7 Incubate the supernatant in a water bath at  60 °C for  00:20:00 for protein clarification.

20m

8.8 Centrifugate at  11000 x g, 4°C, 00:30:00 . Recover the clarified supernatant and discard the pellet.

30m

Note

Place the supernatant  On ice or store at 4°C until use.

Purification of recombinant Bst by FPLC

9 Sample preparation.

Note

Keep all protein samples  On ice during the purification process to avoid degradation.

9.1 Filter the clarified supernatant through a 0.45 µm membrane.

9.2 Dilute the clarified supernatant with saline buffer at a ratio of 1:1.

Note

Saline buffer composition (SB).

	A	B
	Tris-HCl pH 7.5	50 mM
	NaCl	100 mM
	PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

9.3 Load the diluted fraction onto a 150 mL Superloop (Cytiva). Store at  4 °C until use.

10 Immobilized metal affinity chromatography (Ni²⁺-IMAC).

10.1 Connect a  HisTrap HP 5mL **Cytiva Catalog #17524801** to a FPLC system.

Equipment

ÄKTA pure	NAME
Protein purification system	TYPE
Cytiva	BRAND
29046665	SKU
https://www.cytivalifesciences.com/en/us/support/products/akta-pure-150-l-29046665	LINK

- 10.2 Equilibrate the column with 8 column volumes (CV) of washing buffer-AI (WB-AI) at a flow of 2.5 mL/min.

Note

Washing buffer-AI composition (WB-AI).

	A	B
	Tris-HCl pH 7.5	50 mM
	NaCl	100 mM
	Imidazole	10 mM
	PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

- 10.3 Connect the Superloop charged with the protein fraction and load the sample onto the column at a flow of 2.5 mL/min.
- 10.4 Wash the column with 10 CV of WB-AI at a flow of 2.5 mL/min.

- 10.5 Wash the column with 10 CV of 2% elution buffer-AI (EB-AI) at a flow of 2.5 mL/min.

Note

Elution buffer-AI composition (EB-AI).

	A	B
	Tris-HCl pH 7.5	50 mM
	NaCl	100 mM
	Imidazole	500 mM
	PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

- 10.6 Elute the proteins by passing 5 CV of 100% EB-AI through the column using a flow of 2.5 mL/min.

- 10.7 Analyze all recolected fractions by a 8% Tricine-SDS-PAGE gel electrophoresis. Load  10 µL protein sample per well.

- 10.8 Pool all elution fractions carrying the recombinant Bst protein. Store at  4 °C until use.

11 Desalting step.

- 11.1 Connect a  HiPrep 26/10 Desalting Column **Cytiva Catalog #17508701** to the FPLC system.
- 11.2 Wash the column with 2.5 CV of Mili-Q water. Then, equilibrate the column with 2 CV of desalting buffer-A (DB-A). For both steps use a flow of 10 mL/min.

Note

Desalting buffer-A composition (DB-A).

	A	B
	Tris-HCl pH 7.5	50 mM
	KCl	20 mM
	EDTA	1 mM
	DTT	1 mM
	PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

11.3 Load the sample onto the column at a flow of 5 mL/min.

11.4 Wash the column with 2 CV of DB-A for protein elution at a flow of 10 mL/min.

11.5

Analyze all collected fractions by qualitative Bradford assay using the

 Quick Start™ Bradford 1x Dye Reagent **Bio-Rad Laboratories Catalog #5000205** .

Pool the fractions with higher protein concentration.

11.6 Load the pool of desalted fractions onto a 150 mL Superloop (Cytiva). Store at  4 °C until use.

12 **Heparin affinity chromatography.**

12.1 Connect a  HiTrap Heparin HP affinity column **Cytiva Catalog #17040701** to the FPLC system.



12.2 Equilibrate the column with 10 CV of washing buffer-All (WB-All) at a flow of 2 mL/min.

Note

Washing buffer-All composition (WB-All).

	A	B
	Tris-HCl pH 7.5	50 mM
	EDTA	1 mM
	DTT	1 mM
	PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

12.3 Connect the Superloop charged with the protein fraction and load the sample onto the column at a flow of 2 mL/min.

12.4 Wash the column with 5 CV of WB-All at a flow of 2 mL/min.

12.5 Elute proteins by washing the column with a linear gradient of 10 CV of elution buffer-All (EB-All). Use a flow of 2 mL/min.

Note

Elution buffer-All composition (EB-All).

	A	B
	Tris-HCl pH 7.5	50 mM
	KCl	1 M
	EDTA	1 mM
	DTT	1 mM
	PMSF	1 mM

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

- 12.6 Analyze all collected fractions by Tricine-SDS-PAGE electrophoresis through a 8% polyacrylamide gel. Load  10 µL protein sample per well.
- 12.7 Pool all elution fractions carrying the recombinant Bst protein. Store at  4 °C until use.
- 13 **Purified Bst enzyme concentration and formulation.**
- 13.1 Load the purified Bst enzyme pool onto a dialysis membrane (pre-hydrated). Place the membrane into a beaker with precooled storage buffer-A (SB-A) at a ratio 1:50 (v/v).

**Note**

Storage buffer-A composition (SB-A).

	A	B
	Tris-HCl pH 7.5	10 mM
	KCl	50 mM
	EDTA	0.1 mM
	DTT	2 mM
	Glycerol	50% v/v

Prepare the buffer with Milli-Q water and adjust to pH 7.5. Store at 4°C.

13.2 Dialyze  Overnight at  4 °C with slow agitation.

13.3 Recover the dialized protein, load it onto an **Amicon Ultra-15ML - 30 kDa cutoff centrifugal filter**. Concentrate until a concentration equal or higher than  1 mg/mL .



Equipment

Amicon Ultra-15

NAME

PLTK Ultracel-PL membrane, 15 ML - 30 kDa cutoff

TYPE

Millipore

BRAND

UFC903024

SKU

https://www.merckmillipore.com/MX/es/product/Amicon-Ultra-15-Centrifugal-Filter-Unit,MM_NF-UFC903024?ReferrerURL=https%3A%2F%2Fwww.bing.com%2F&bd=1

LINK

Note

Monitor protein concentration measuring absorbance at 280 nm using a NanoDrop spectrophotometer.

Equipment

NanoDrop™ One UV-Vis Spectrophotometer

NAME

spectrophotometer

TYPE

Thermo Scientific

BRAND

ND-ONE-W

SKU

<https://www.thermofisher.com/order/catalog/product/ND-ONE-W>

LINK

Sample Volume (Metric): Minimum 1 µL; Spectral Bandwidth: ≤1.8 nm (FWHM at Hg 254 nm); System Requirements: Windows™ 8.1 and 10, 64 bit; Voltage: 12 V (DC); Wavelength Range: 190–850 nm

SPECIFICATIONS



- 13.4 Prepare aliquots of 50 μL of concentrated Bst enzyme .
- 13.5 Add 0.1 % (v/v) triton X-100 to the enzyme aliquots and store at -20 $^{\circ}\text{C}$.
- 13.6 Determine final protein concentration by measuring absorbance at 280 nm in a NanoDrop spectrophotometer.
- 13.7 Analyze the final Bst enzyme formulation by Tricine-SDS-PAGE electrophoresis through a 8 gel. Load 3 μL protein sample per well. Load 3 μL protein ladder
- Precision Plus Protein™ Unstained Protein Standards **Bio-Rad**
Laboratories Catalog #1610363
- .
- 13.8 Analyze the electrophoresis gel by densitometry using the **Image Lab 6.1 Software (Bio-Rad)**. Determine protein concentration for each Bst enzyme aliquot analyzed using the protein ladder as weight standard.

Note

The protein ladder



Precision Plus Protein™ Unstained Protein Standards **Bio-Rad**
Laboratories Catalog #1610363

includes three reference bands: the 50 KDa with 750 ng, the 20 KDa and 100 KDa bands with 150 ng each per 10 μL of the protein ladder mix.

Citations

Step 6.2

Hermann Schägger. Tricine-SDS-PAGE
[10.1038/nprot.2006.4](https://doi.org/10.1038/nprot.2006.4)