

Jan 22, 2023

## Expression and purification of recombinant UvsX recombinase

Forked from a private protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.5qpvobemdl4o/v1>

lucero.merino.c<sup>1</sup>, lucero.mascaro.r<sup>1</sup>

<sup>1</sup>Universidad Peruana Cayetano Heredia (UPCH, Peru)



**Lucero Merino**

Universidad Peruana Cayetano Heredia

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.5qpvobemdl4o/v1>

**Protocol Citation:** lucero.merino.c , lucero.mascaro.r 2023. Expression and purification of recombinant UvsX recombinase. protocols.io <https://dx.doi.org/10.17504/protocols.io.5qpvobemdl4o/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working



**Created:** August 11, 2022

**Last Modified:** January 22, 2023

**Protocol Integer ID:** 68488

**Keywords:** FPLC, UvsX, RPA, purification of recombinant uvsx recombinase, recombinant uvsx recombinase, recombinant uvsx, part for an isothermal dna amplification, isothermal dna amplification, recombination process, uvsx, purification, enzyme, isothermal amplification technique, uvsy, rpa reaction, other protein, protocols for the production, dna

**Funders Acknowledgements:**

ProCiencia, Peru

Grant ID: 170-2020

## Abstract

The UvsX is a enzyme that is part for an isothermal DNA amplification based on the recombination process, the RPA reaction.

RPA uses 4 enzymes: UvsX, UvsY, Bsu and Gp32. It's an isothermal amplification technique that can run at 37°C.

In this protocol we are producing a recombinant UvsX that has a 6xHIS-tag using a *E. coli* expression system.

The protocols for the production of the other proteins are also available in protocols.io.



## Materials

### Binding buffer, pH 7.2

[M] 50 millimolar (mM) Tris-HCl, pH 7.2

[M] 20 millimolar (mM) Imidazole, pH 7.2

[M] 500 millimolar (mM) KCl

[M] 5 % (v/v) Glycerol

[M] 0.1 millimolar (mM) PMSF

[M] 0.03 % (v/v) 2-Mercaptoethanol (BME)

### Elution buffer, pH 7.2

[M] 50 millimolar (mM) Tris-HCl, pH 7.2

[M] 500 millimolar (mM) Imidazole, pH 7.2

[M] 500 millimolar (mM) KCl

[M] 5 % (v/v) Glycerol

[M] 0.03 % (v/v) BME

### Storage buffer, pH 8

[M] 20 millimolar (mM) Tris-HCl, pH 8

[M] 400-500 millimolar (mM) KCl

[M] 20 % (v/v) Glycerol

[M] 1 millimolar (mM) DTT

### Ladder:

⊗ Pageruler Prestained Protein Ladder **Thermo Fisher Scientific Catalog #26616**

### Equipment:

**Sonicator** OMNI Sonic Ruptor 400

**Protein purification system** FPLC AKTA START

## Protocol materials

⊗ Pageruler Prestained Protein Ladder **Thermo Fisher Scientific Catalog #26616**

⊗ Amicon Ultra-15 Centrifugal Filter Unit **Merck Millipore (EMD Millipore) Catalog #UFC910024**

⊗ Phusion High-Fidelity DNA Polymerase - 500 units **New England Biolabs Catalog #M0530L**



## DAY1: Transformation of competent cells

1d

- 1 Quantify the plasmid containing the UvsX recombinase gene and determine the volume that contains  100 ng of the plasmid.
- 2 Defrost the aliquot of BL21(DE3) chemically competent cells  On ice . Softly pipette  100 ng of the plasmid in the aliquot and let the tube rest  On ice for  00:30:00 .
- 3 Incubate the tube at  42 °C for  00:00:30 .
- 4 Quickly return the tube  On ice and incubate for  00:05:00 .
- 5 Add the mixture to a microcentrifuge tube with  800 µL SOC medium and incubate at  37 °C for  00:45:00 .
- 6 Centrifuge the tube  4500 rpm, Room temperature, 00:08:00 .
- 7 Discard  800 µL of the supernatant and gently resuspend the pellet with the rest.
- 8 Add the resuspension to a LB agar plate with [M] 0.05 mg/mL Kanamycin and spread the recently transformed cells. Incubate plate  Overnight at  37 °C .

## DAY2: Preparation of pre-inoculum

1d

- 9 For verification that the colonies in the plate contain the desired plasmid with the protein sequence, perform a PCR colony using universal T7 primers and the PCR protocol for Phusion DNA Polymerase

 Phusion High-Fidelity DNA Polymerase - 500 units **New England Biolabs Catalog #M0530L**

. Use the following thermocycling procedures for the UvsX recombinase plasmid:



| Step                 | Temperature (°C) | Time   | Cycles |
|----------------------|------------------|--------|--------|
| Initial Denaturation | 98               | 3 min  | 1      |
| Denaturation         | 98               | 30 sec | 25     |
| Annealing            | 60               | 30 sec |        |
| Extension            | 72               | 42 sec |        |
| Final Extension      | 72               | 5 min  | 1      |
| Hold                 | 4                | ∞      |        |

Run the PCR product in a 1% agarose gel and verify if there is a band of the desired weight (UvsX insert = 1387 bp).

- 10 Select an isolated bacterial colony from the plate and inoculate a test tube with  10 mL LB medium and  0.05 mg/mL Kanamycin . Incubate the tube  Overnight at  220 rpm, 37°C .

## DAY 3-A:Protein expression in small scale

2d

- 11 Inoculate  50  $\mu$ L from the pre-inoculum to an Erlenmeyer flasks with  50 mL LB medium and  0.05 mg/mL Kanamycin . Incubate at  220 rpm, 37°C until OD<sub>600</sub> = 0.5 - 0.6 (3-4 hours).
- 12 Add IPTG to a final concentration of  0.5 millimolar (mM) and incubate at  220 rpm, 18°C, 16:00:00 .
- 13 Centrifuge the cell culture  8000 rpm, 4°C, 00:05:00 . Discard the supernatant. At this point, you may store the cells pellet at -20°C until you are ready to run the purification.

5m

## DAY 4-A:Protein purification in resin

1d

- 14 Resuspend the cell pellet in  5 mL Binding buffer . Then add lysozyme to a final concentration of  0.1  $\mu$ g/ $\mu$ L .
- 15 Incubate the cells at  220 rpm, Room temperature , 00:20:00 and add 10% SDS to a final concentration of 0.02%.



- 16 Add ~  100  $\mu\text{L}$  of glass beads and shake vigorously in a vortex for  00:20:00 at room temperature. You can do this by fixing a 15 mL tube to the vortex rubber platform with tape.
- 17 Centrifuge at  13500 rpm, 4°C, 00:07:00 . Collect the supernatant and label it as a Soluble fraction. The pellet is the Insoluble fraction. Collect small fractions of each one to run an acrylamide gel afterwards.
- 18 **Prepare the resin.** Homogenize resin with its storage buffer by shaking the bottle and transfer it to a new tube. You will use  330  $\mu\text{L}$  of resin for each  1 mL of soluble fraction. Let the slurry sediment or spin it down. Remove the storage buffer and wash the resin in **Binding buffer**. Wash the resin with the same volume as the obtained soluble fraction. Repeat this wash step 3 times.
- 19 Add the soluble fraction to the resin. Homogenize the mixture gently in an orbital shaker for 20 min (~60 RPM) at room temperature.
- 20 Let the resin sediment for 10 minutes. Collect a small fraction of the supernatant to run an acrylamide gel afterwards, and discard the remainder. Resuspend resin with 1 mL of **Binding buffer**. Homogenize the tube gently with finger taps. Don't flip the tube (1st washing step).
- 21 Spin down for a few seconds and discard supernatant. Resuspend resin with 1 mL of **Binding buffer**. Homogenize the tube gently with finger taps. Don't flip the tube (2nd washing step).
- 22 Spin down for a few seconds and discard supernatant. Resuspend resin with 1 mL of **Elution buffer (50mM Imidazole)**. Homogenize the tube gently with finger taps. Don't flip the tube. Incubate for  00:10:00 . 10m
- 23 Spin down for a few seconds and collect the supernatant. Resuspend resin with 1 mL of **Elution buffer (500mM Imidazole)**. Homogenize the tube gently with finger taps. Incubate for  00:10:00 . Collect small fractions of elutions to run an acrylamide gel afterwards. 10m

Run a 12% acrylamide gel at 200 V to evaluate all the samples you just generated: Lysis sample, Soluble fraction, Insoluble fraction, Flowthrough, 1st washing step, 2nd washing step and Eluted fraction.

## DAY 3-B: Protein expression in medium scale

2d



- 24 Inoculate  2.5 mL from the pre-inoculum to an Erlenmeyer flask with  250 mL LB medium and  0.05 mg/mL Kanamycin , use 4 flasks to obtain 1L of cell culture. Incubate at  220 rpm, 37°C until OD<sub>600</sub> = 0.5 - 0.6 (3-4 hours).
- 25 Add IPTG to a final concentration of  0.5 millimolar (mM) to each flask and incubate at  220 rpm, 18°C, 16:00:00 .
- 26 Centrifuge the cell culture  4000 rpm, 4°C, 00:20:00 . Discard the supernatant. At this point, you may store the 1-2 grams of cell pellet at -20°C until you are ready to run the purification. 20m

## DAY 4-B: Cells Lysis

- 27 Resuspend all the cell pellets in  100 mL Binding buffer . Add PMSF to a final concentration of  0.1 millimolar (mM) . Add lysozyme to a final concentration of  0.1 µg/µL .
- 28 Incubate the cells at  220 rpm, Room temperature , 00:20:00 .
- 29 Sonicate on ice until the lysate turns translucent. Use 5 cycles of  00:15:00 power ON, pulse 10 . Then  00:15:00 power OFF , with the tube on ice.
- 30 Centrifuge  6000 rpm, 4°C, 00:20:00 to separate the insoluble fraction (pellet) from the soluble fraction. Transfer the soluble fraction to a new and clean tube on ice. Collect small fractions of each one to run an acrylamide gel afterwards.

## DAY 4-B: Protein Purification with FPLC

1d

- 31 Prepare the 5 mL HisTrap column in the FPLC system. Wash the tubes, pumps system and the column with 7 column volumes (c.v.) of distilled and filtrated water. Then equilibrate the column with 7 c.v. of **Binding buffer**.
- 32 Load the soluble fraction to the FPLC system at a flow of 1 mL/min. Collect a small fraction of each step and signal change to run an acrylamide gel afterwards. Wash the column with 5 c.v. of **Binding buffer**, until the UV and conductivity signal stabilizes.
- 33 **Washing:** Load the column with 17% of pump B (**Elution Buffer**), which is equivalent to ~100 mM Imidazole, until the signal stabilizes.

**Elution:** Load the column with 48% of pump B (**Elution Buffer**), which is equivalent to ~250 mM Imidazole, until the signal stabilizes. Collect the elution in 8 mL tube fractions and in case the UV signal changes. Then load the column with 3 c.v. of 100% of pump B (**Elution Buffer**), which is equivalent to 500 mM Imidazole, until the signal stabilizes.

34 **Wash the column for storage.** Wash the FPLC system with distilled and filtrated water. Load the column with 7 c.v. of distilled and filtrated water. To storage the column, load it with 5 c.v. of **ethanol 20%** and storage it at 4°C. Finally, remove the rest of the water from the system with **ethanol 20%** and keep the system with it until next use.

35 Determine the fractions with the UvsX recombinase by running a SDS-PAGE in a 10% acrylamide gel. The UvsX recombinase weights ~44 kDa.

36 Concentrate the eluted fractions with the protein with an



Amicon Ultra-15 Centrifugal Filter Unit **Merck Millipore (EMD Millipore) Catalog #UFC910024**

**10kDa.** Reconstitute the concentrate so it is stored with the components detailed in **Storage Buffer** to decrease the Imidazol to 20 mM or less. Add glycerol to a 20%, homogenize, make aliquots of 400 µL of the protein and storage them at -80°C.

