

May 25, 2022

Expression and purification protocol of GST-mCh-FYVE

 Forked from [Expression and purification protocol of GST-NDP52](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov1n39pgr2/v1>

Liv Jensen¹, Chunmei Chang²

¹James Hurley Lab, University of California Berkeley; ²James Hurley Lab, UC, Berkeley



Liv Jensen

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov1n39pgr2/v1>

Protocol Citation: Liv Jensen, Chunmei Chang 2022. Expression and purification protocol of GST-mCh-FYVE. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1n39pgr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: May 02, 2022

Last Modified: May 31, 2024

Protocol Integer ID: 61821

Keywords: ASAPCRN, purification protocol of gst, purification of gst, purification protocol, purification, gst, fyve this protocol, mcherry, protocol detail

Abstract

This protocol details the expression and purification of GST-mCherry-FYVE.

Attachments



[246-484.docx](#)

1.4MB



Materials

General information:

	A	B
	Expression system	E.Coli BL21DE3
	Medium	Luria Bertani
	Plasmid origin	
	Backbone	pGST2
	Resistance	Amp
	Insert	
	Tags & cleavage sites	N-term GST
	Ext coeff	82320 M-1cm-1, MW 77.9 kDa

Lysis Buffer:

	A	B
	Hepes pH=7.5	50 mM
	NaCl	300 mM
	TCEP	1 mM
	Protease Inhibitors (Roche)	

Wash Buffer:

	A	B
	Hepes pH=7.5	50 mM
	NaCl	150 mM
	TCEP	1 mM

Elution Buffer:

	A	B
	Hepes pH=7.5	50 mM



	A	B
	NaCl	150 mM
	TCEP	1 mM
	Glutathione	50 mM

SEC Buffer:

	A	B
	Hepes pH=8	20 mM
	NaCl	150 mM
	TCEP	1 mM

Columns/Resin:

- Glutathione Sepharose 4B
- S6_10/300 Increase



Protein expression

20m

- 1 Transform the E.Coli BL21DE3 cells with plasmid encoding for GST-mCh-FYVE and plate them on Amp plate.
- 2 Carry out protein expression in  1 L medium, induce with [M] 200 micromolar (μM) IPTG (isopropyl- β -d-thiogalactopyranoside) to an OD_{600} of 0.8 and grow at  18 °C  Overnight . 
- 3 Harvest the cells by spinning at  4500 x g for  00:20:00 at  4 °C and stock at  -80 °C until purification. 

20m

Protein purification

5h

- 4 Follow the GST batch purification by Size Exclusion Chromatography.
- 5 Resuspend the pellets in Lysis Buffer, sonicate for cell lysis and clear at  16000 rpm at  4 °C for  01:00:00 
- 6 Incubate the supernatant with Glutathione Sepharose 4B (GE Healthcare) at  4 °C with gentle shaking for  00:30:00 , apply to a gravity column, and wash extensively with Wash Buffer. 
- 7 Elute the protein of interest with Elution Buffer and then apply onto a Superdex 6 column (10/300 Increase) pre-equilibrated in SEC Buffer at  4 °C .
- 8 Pool the peak fractions containing pure protein, snap-frozen in liquid nitrogen, and store at  -80 °C .

30m