Jun 10, 2025

Extraction of RNA/DNA from Dried Blood on Filter Papers after Long-Term Storage

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Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS

DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpbbbr5lzp/v1>

Protocol Citation: Gabriela Ulloa Urizar 2025. Extraction of RNA/DNA from Dried Blood on Filter Papers after Long-Term Storage. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gpbbbr5lzp/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: February 01, 2022

Last Modified: June 10, 2025

Protocol Integer ID: 57661

Keywords: dna from dried blood, extraction of rna, start protocol allprep dna, dried blood, rna mini kit, dna, extraction, filter papers after long, rna, filter paper, qiagen quick, blood, term storage modified protocol from the manufacturer, qiagen

Abstract

Modified protocol from the manufacturer's protocol (Qiagen Quick-Start Protocol AllPrep DNA/RNA Mini Kit)

Cutting the sample

- 1 Using a sterile punch, cut 10 to 12 3 mm circles of blood from the filter paper (the aim is to have at least the equivalent of 50ul of blood for each collection).

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Dilute the sample for extraction

- 2 Add  600 μL of RLT Plus Buffer to the 1.5 mL microtube with the blood sample on filter paper cut into small pieces (in 3 mm circles, equivalent to 50ul of whole blood), spin for  00:00:15 at 10000 rpm. 15s

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- 3 Put the 1.5 mL microtube in a ThermoMixer at 55°C for  00:30:00 at 2000 rpm. 30m

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- 4 Add  65 μL of proteinase K, vortex and spin for  00:00:15 at 10000 rpm. 15s

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- 5 Put the 1.5 mL microtube in a ThermoMixer at 37°C for at least  02:00:00 (120 minutes) at 2000 rpm. 2h

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- 6 Centrifuge the 1.5 mL tube with the sample for  00:03:00 at 13500 rpm. 3m

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Sample separation for DNA/RNA extraction

- 7 Pipette the supernatant carefully to avoid touching the pellet with the tip (approx.  650 μL). Transfer the supernatant into the DNA column (Allprep DNA) with collecting tube. Close carefully and centrifuge for  00:00:30 at 10000 rpm.

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- 8 Put the DNA column in a new collecting tube and reserve it. Recover the centrifuged sample (flow through) from the first collection tube and add the same volume of 70° ethanol (approx.  650 μL), mix well with the pipette.

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RNA extraction

- 9 Transfer  650 μL of the sample homogenate with ethanol to the RNA column with collecting tube, centrifuge for  00:00:30 at 10000 rpm. Discard the flow through and add the remaining 650ul of the sample homogenate, centrifuge for 30 seconds at 10000 rpm.
NOTE: Keep the collecting tube for the next 4 steps.

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- 10 Add  700 μL of RW1 buffer to the RNA column with collecting tube, centrifuge for  00:00:30 at 12000 rpm. Discard the flow through.

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- 11 Add  500 μL RPE buffer to the RNA column with collecting tube, centrifuge for  00:00:30 at 12000 rpm. Discard the flow through.

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- 12 Add  500 μL RPE buffer to the RNA column with collecting tube, centrifuge for  00:02:00 at 12000 rpm. Discard the flow through.

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- 13 Centrifuge the RNA column again for  00:01:00 at 13500 rpm for better membrane drying.

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- 14 Put the RNA column into a new 1.5mL tube. Add  30 μL of RNase-free water directly onto the column membrane without touching it. Wait 1 minute to hydrate the column. Centrifuge for  00:01:00 at 12500 rpm.

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- 15 Add another  30 μL of RNase-free water directly onto the column membrane without touching it. Wait 1 minute to hydrate the column. Centrifuge for  00:01:00 at 12500 rpm.

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- 16 Discard the RNA column, correctly label the 1.5 mL microtube and store the extracted RNA at -20°C for early use.

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DNA extraction

7m 30s

- 17 To the DNA column reserved in step 7. Add  500 μL of Buffer AW1. Centrifuge for  00:00:30 at 12500 rpm. Discard the flow through. 30s
NOTE: Keep the collecting tube for the next 2 steps.

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