

Mar 25, 2020

## Extraction of total RNA from *E. coli* cells

DOI

[dx.doi.org/10.17504/protocols.io.gtnbwme](https://dx.doi.org/10.17504/protocols.io.gtnbwme)

Alice Pawlowski<sup>1</sup>

<sup>1</sup>Institute of Synthetic Microbiology, CEPLAS, Heinrich Heine University Duesseldorf

Axmann Lab



Alice Pawlowski

Heinrich-Heine Universität Düsseldorf

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.gtnbwme>

**Protocol Citation:** Alice Pawlowski 2020. Extraction of total RNA from *E. coli* cells. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.gtnbwme>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 19, 2016

**Last Modified:** March 25, 2020

**Protocol Integer ID:** 4686

**Keywords:** extraction of total rna, total rna, rna, small regulatory rna, expressing small regulatory rna, trrna, extraction

## Abstract

The protocol is used for the extraction of total RNA from *E. coli* cells. It is based on the method described by Chomczynski and Sacchi, 1987 ([https://doi.org/10.1016/0003-2697\(87\)90021-2](https://doi.org/10.1016/0003-2697(87)90021-2)). Total RNA is isolated from cells expressing small regulatory RNAs (trRNAs) either constitutively (*E. coli* W3110) or after induction with anhydrotetracycline (*E. coli* W1310 Z1).

## Materials


### MATERIALS

 Roti-Aqua-P/C/I **Carl Roth Catalog #X985.2**


 Roti-C/I **Carl Roth Catalog #X984.2**

 NucleoZOL **Macherey-Nagal Catalog #REF 740404.200**

 Roti-Aqua-Phenol for RNA extraction **Carl Roth Catalog #A980.1**

 LB-Medium (Lennox) vegetal **Carl Roth Catalog #0155.1**

 IPTG **Carl Roth Catalog #CN08.1**

 Anhydrotetracycline **Catalog #2-0401-001**

additional solutions/buffers:

- Antibiotic stock solutions:
  - Kanamycin (Km): 25 mg/ml
  - Spectinomycin (Spec): 100 mg/ml
- Inducer stock solutions: 20 µg/ml Anhydrotetracycline (aTc) in EthOH, 10 mM Isopropyl-β-D-thiogalactopyranosid (IPTG ) in EthOH and a mixture of 20 µg/ml aTc and 10 mM IPTG in EthOH
- Ethanol
- Isopropanol
- 3 M Na-Acetat, pH 5.2
- RNase free molecular grade water (DEPC-treated)

## Troubleshooting

## Safety warnings



Phenol is toxic! Work under the hood, wear protective gloves (Nitril) and change gloves immediately after contamination. Collect solid and liquid waste in special waste containers.

## Before start

RNA is sensitive to degradation! Wear gloves, keep samples on ice when possible, use filter-tips and RNase free reagents. Pre-cool centrifuges and store isolated RNA-samples immediately at -20 or -80°C.

## incubation of cells

- 1
  - freshly transform *E. coli* cells with plasmids encoding trRNA sequences or empty vector
  - inoculate 3 ml LB-vegetal medium + antibiotic (25 µg/ml Km for W1310, 25 µg/ml + 100 µg/ml Spec for W1310Z1) in culture tubes with a single colony
  - incubate o/n at 37°C and 230 rpm

for constitutive RNA-synthesis (*E. coli* W3110):

- dilute 1:100 in fresh medium with 25 µg/ml Km and incubate until OD600 reaches 0.7 - 1.0 (in culture tubes or 6 well plates), proceed to step 2

for induced RNA-synthesis (*E. coli* W3110Z1):

- dilute 1:50 in fresh medium ( 200 µl o/n culture in 10 ml medium + antibiotic in 100 ml Erlenmeyer flask) and incubate until OD600 reaches 0.4 - 0.6
- if you have several different cultures that do not grow at same speed: store samples on ice until the last one reaches the defined OD (keep that one on ice as well for 10 min)
- meanwhile prepare 12 well plates with inducer: pipet in each well of a row 10 µl of either
  - 1) EthOH,
  - 2) aTC (= 200 ng/ml final concentration)
  - 3) IPTG ( = 100 µM final concentration)
  - 4) aTc + IPTG ( 200 ng/ml, 100 µM final concentration)
- add 1 ml of culture and incubate at 37°C and 230 rpm for 1 h
- proceed with step 2

## stop of RNA synthesis

- 2
  - work under the hood
  - mix 1 ml of cells with 200 µl 'stopmix'- solution (5 % phenol in ethanol) in a 2 ml safe lock tube tube → stops RNA production in the cells
  - centrifuge for 5 min at 4°C and 14000 x g
  - discard the supernatant and resuspend the pellet in 1 ml NucleoZOL (Macherey and Nagel), place on ice (better for rapid freezing: dry ice or liquid nitrogen)

→ store cells at -20 or -80 °C (only for a short time, maximum 2 weeks) or proceed to next step

## RNA-isolation

- 3
  - incubate the sample at 65 °C and 250 rpm (Thermomixer) for 10 min



- mix with 400  $\mu$ l Phenol-Chloroform/Isoamylalcohol (Roti®-)Aqua-P/C/I) by inverting for 10 s
- centrifuge at 4°C for 10 min at 14000 x g
- transfer aqueous (upper) phase to a fresh 1,5 ml safe lock reaction tube, work on ice
- mix with 450  $\mu$ l Chloroform/Isoamylalcohol (Roti®-C/I)
- centrifuge at 4°C for 10 min at 14000 x g
- transfer aqueous (upper) phase to a new reaction tube and add 1 Vol. icecold Isopropanol + 1/10 vol (e.g. 20  $\mu$ l for 200  $\mu$ l Isopropanol) 3 M Na-Acetate (pH 5.2), mix and store at least 30 min at -20 °C or -80 °C
- centrifuge at 4°C for 30 min at 14000 x g
- remove the supernatant (take care of the RNA-pellet) and add 350  $\mu$ l of icecold 75% Ethanol
- centrifuge for 5 min at 4°C and 15000 rpm
- add again 350  $\mu$ l of icecold 75% ethanol and centrifuge for 5 min at 4°C and 15000 rpm
- remove the supernatant and dry the pellet at room temperature for ca.15 min
- resuspend the pellet in 30  $\mu$ l Molecular Biology Grade Water and store at -80°C

4