



Nov 07, 2024

Extraction Protocol for untargeted LC-MS/MS - Animal tissues

 [PLOS Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr8x8blmk/v1>

Mauricio Caraballo¹, Marwa Tawfik¹, Janet Tseng¹, Monica Wu¹

¹University of California San Diego

Mauricio Caraballo: MSCollaboratory Director

Marwa Tawfik: MSCollaboratory visitor

Janet Tseng: undergraduate student

Monica Wu: Undergraduate student



Mauricio Caraballo

University of California San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

[Create free account](#)

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr8x8blmk/v1>

External link: <https://doi.org/10.1371/journal.pbio.3003705>

Protocol Citation: Mauricio Caraballo, Marwa Tawfik, Janet Tseng, Monica Wu 2024. Extraction Protocol for untargeted LC-MS/MS - Animal tissues. protocols.io <https://dx.doi.org/10.17504/protocols.io.6qpvr8x8blmk/v1>

Manuscript citation:

Carrier TJ, Rufino-Navarro A, Knoop T, Repnik U, Caraballo-Rodríguez AM, Needham DM, Bang C, Franzenburg S, Bramkamp M, Rath W, Biastoch A, Hernández JC, Hentschel U (2026) Sea urchin eggs contain a plastid-derived structure that contributes to their development. PLOS Biology 24(4). doi: [10.1371/journal.pbio.3003705](https://doi.org/10.1371/journal.pbio.3003705)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: March 19, 2024

Last Modified: November 07, 2024

Protocol Integer ID: 96958

Keywords: Metabolomics, Untargeted LC-MS/MS, untargeted metabolomics analysis, mass spectrometry, tandem mass spectrometry, extraction protocol for untargeted lc, preparation of animal tissue sample, animal tissue sample, liquid chromatography, using liquid chromatography, ms analysis, preparation for lc, tissue sample, animal tissues this protocol, m, extraction protocol, extraction, untargeted lc, metabolomic

Funders Acknowledgements:

Gordon and Betty Moore Foundation

Grant ID: GBMF12120

Abstract

This protocol details the preparation of animal tissue samples for untargeted metabolomics analysis using liquid chromatography-tandem mass spectrometry (LC-MS/MS). It covers sample collection, homogenization, metabolite extraction, and preparation for LC-MS/MS analysis.

Protocol materials

⊗ Water, Optima™ LC/MS Grade, Fisher Chemical™ **Fisher Scientific Catalog #W64**

⊗ Methanol, Optima™ LC/MS Grade, Fisher Chemical **Fisher Scientific Catalog #A456-4**



Description

- 1 This protocol details the preparation of animal tissue samples for untargeted metabolomics analysis using liquid chromatography-tandem mass spectrometry (LC-MS/MS). It covers sample collection, homogenization, metabolite extraction, and preparation for LC-MS/MS analysis. This protocol can be adapted to other sample types. This is not a final document, not peer-reviewed and it can be updated.

Solvents

- 2

	Water, Optima™ Scientific Catalog #W64	LC/MS	Grade,	Fisher	Chemical™ Fisher
	Methanol, Optima™ Scientific Catalog #A456-4	LC/MS	Grade,	Fisher	Chemical Fisher

Note

Additionally, prepare the following solutions:

- Extraction solvent: 50% methanol:water → 50:50 v/v of LC grade methanol/ LC grade water.
- Resuspension solvent (50:50 v/v methanol:water) containing 1uM sulfadimethoxine (CAS 122-11-2) as an internal standard.

Equipment

- 3

Equipment

TissueLyser II

NAME

Bead Mill

TYPE

QIAGEN

BRAND

85300

SKU

<https://www.qiagen.com/us/products/human-id-and-forensics/automation/tissuelyser-ii/#orderinginformation>

LINK

Equipment

Legend RT

Refrigerated Benchtop Centrifuge

Sorvall

sorvall legend rt

https://knowledge1.thermofisher.com/Lab_Equipment/Centrifuges/Table_Top_Centrifug

Equipment

Centrifuge 5418

NAME

Centrifuge

TYPE

Eppendorf US

BRAND

Catalog No. 5401000137

SKU

<https://www.eppendorf.com/us-en/eShop-Products/Centrifugation/Microcentrifuges/Centrifuge-5418R-p-5401000137>

LINK

Equipment

Sonicator

NAME

Branson

BRAND

B2510

SKU

Consumables

- 4 2 mL Qiagen sample tubes (catalog no. 990381)
- Stainless Steel Beads (Qiagen, Inc., 5 mm diameter, Cat. No. NC9257481)
- Shallow polypropylene 96-well plate (Cat. No. 951040048)
- Deep polypropylene 96-well plate (Cat. No. 951033405)
- Seal plates (Cat. No. ZAFPE50) **Note:** for LC-MS/MS acquisition
- AlumaSeal CS (Sealing Foil, Excel Scientific) (Cat. No. FC100) **Note:** for long time storage
- Metal tweezers
- Bead dispenser



Note

Wash the stainless steel beads with LC-MS grade water followed by LC-MS methanol in a glass Petri dish then allow to dry overnight.

Sample description

- 5 Liver, stomach and intestine preserved in Ethanol 95% from fish (threespine stickleback, *Gasterosteus aculeatus*).

Note

Worms (fish parasite) can be extracted following this protocol

Note

In case samples were put in tubes without ethanol: Shake (vortex can be used) the tube without touching/manipulation of the samples and transfer to a new empty tube.

Note

Transferring beads from the Petri dish or big beaker to the beads dispenser using aluminum foil to cover either of them and to create a small mouth to either the big beaker or the Petri dish to avoid falling of beads.

Sample preparation

- 6 Samples are received in plastic centrifugation tubes containing 1mL of 95% Ethanol.

Note

If ink labels are inside the tube, carefully remove them using tweezers disinfected with methanol every time. Make sure the tube is labeled from the outside via sticker or marker.

- 7 Dry samples using Centrivap.



- 8 Label empty sample tube.
- 9 Add 1 bead with the beads dispenser into each empty sample tube.
- 10 Add 1 mL of extraction solvent (50% methanol:water) per tube.
- 11 Homogenize at 30 Hz for 8 min in a Qiagen TissueLyzer II.
- 12 Allow to extract at 4 °C for 1 h.

Note

This step helps with protein precipitation

- 13 Centrifuge at  14000 rpm, 4°C, 00:15:00

15m

- 14 Transfer at total of 800 uL of supernatant split in 2 labeled deep polypropylene 96-well plates (400 uL for each plate).

Note

- Label plates before transferring any sample.
- One of the plates will be used as a backup plate (storage under -80 °C).

- 15 Dry plates using Centrivap.
- 16 Seal plates and store under -80 °C until LC-MS/MS analysis.



17 Prior to LC-MS/MS analysis, reconstitute one of the plates with 200uL of resuspension solvent (50% methanol/water solution containing 1uM sulfadimethoxine as an internal standard).

18 Seal plate and sonicate  00:15:00

15m

19 Hand-shake or vortex

Note

This guarantees complete solubilization of the sample

20 Transfer to a shallow polypropylene 96-well plate

Note

If particles or turbid solution are observed, centrifuge plate

 2000 rpm, 4°C, 00:15:00 and transfer to a new shallow polypropylene 96-well plate.

21 This shallow polypropylene 96-well plate will be used for LC-MS/MS acquisition.

References

- 22 <https://animalmicrobiome.biomedcentral.com/articles/10.1186/s42523-022-00182-z>
<https://protocolexchange.researchsquare.com/article/pex-2055/v1>

Protocol references

<https://animalmicrobiome.biomedcentral.com/articles/10.1186/s42523-022-00182-z>
<https://protocolexchange.researchsquare.com/article/pex-2055/v1>

Acknowledgements

This work was supported by the Gordon and Betty Moore Foundation, GBMF12120 and <https://doi.org/10.37807/GBMF12120>.