



Jun 10, 2025

Extraction Protocol for untargeted LC-MS/MS - Microbial Cultures

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov19417lr2/v1>

Mauricio Caraballo¹

¹University of California San Diego



Mauricio Caraballo

University of California San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov19417lr2/v1>

Protocol Citation: Mauricio Caraballo 2025. Extraction Protocol for untargeted LC-MS/MS - Microbial Cultures. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov19417lr2/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

This is a working protocol and can be adapted to your purposes.

Created: April 24, 2024

Last Modified: June 10, 2025

Protocol Integer ID: 98893

Keywords: Microbial culture, Metabolomics, LC-MS/MS, extraction of bacterial small molecule, microbial cultures this protocol, extraction protocol for untargeted lc, microbial culture, bacterial small molecule, extraction protocol, extraction, m, untargeted lc

Funders Acknowledgements:

Gordon and Betty Moore Foundation

Grant ID: GBMF12120

Abstract

This protocol details the extraction of bacterial small molecules for untargeted LC-MS/MS.

Guidelines

Definitions (optional)

Metadata:

Sample information.

Untargeted LC-MS/MS:

Liquid chromatography tandem mass spectrometry, not targeting specific molecules.

Materials

 Water Optima™ LC/MS Grade Fisher Chemical™ **Fisher Scientific Catalog #W6-4**

 Methanol Optima™ LC/MS Grade Fisher Chemical **Fisher Scientific Catalog #A456-4**

Extraction solvent: 50% methanol → 50:50 v/v of LC grade methanol/ LC grade water



Metadata - sample information

- 1 Prepare metadata table according to [ReDU2-MS2 Sample Information Template](#) [for more information about the ReDu, see the Preprint from our lab: [Repository-scale Co- and Re-analysis of Tandem Mass Spectrometry Data as well as the ReDu platform](#)].
- 2 Fill the metadata table with all the applicable information. Information from LC-MS/MS can be included later on after data has been acquired.
- 3 Specifically, from each sample we will need at least the NCBITaxonomy (e.g. [1313|Streptococcus pneumoniae](#)).
- 4 Additional columns can be created to include information related to the specificity of the experiments (e.g. microbial interactions: co-cultures vs monocultures).

Sample preparation and microbial extraction

1h 45m

5

Note

Exometabolome: extraction procedure performed on supernatant from microbial cultures [See [Pinu & Villas-Boas, 2017](#)].

Briefly, for **liquid cultures**:

- 5.1 Separate microbial culture from cells either by **centrifugation or filtration**.
 - 5.2 **Add internal standard** to monitor extraction procedure (although it can be optional, the addition of a known compound will be useful to guarantee the efficiency of the extraction procedure).
 - 5.3 **Quenching** (deactivate enzymatic activity) by freeze-drying or any other drying process (such as vacuum-drying): at this step, samples can be shipped to us as they are dry;
 - 5.4 If not quenching but direct extraction is performed, it can be done via liquid-liquid extraction or **SPE (solid-phase extraction)** [the extracted metabolites will depend on the sorbent material used in SPE].
- 6 For **semi-solid cultures** (eg. containing agar), the following protocol can be used:

**Note**

If the culture is performed using 96-deep well plates, then the extraction can be performed using a 96-well SPE plate.

6.1 After the incubation period **freeze/thaw** the plate (10-15mins, 3 times) to **lyse** the cells.

1h 15m

- After the incubation period, **freeze/thaw** the plate ( 00:10:00 -  00:15:00) to **lyse** the cells (1/3).
- After the incubation period, **freeze/thaw** the plate ( 00:10:00 -  00:15:00) to **lyse** the cells (2/3).
- After the incubation period, **freeze/thaw** the plate ( 00:10:00 -  00:15:00) to **lyse** the cells (3/3).

7 Add  1 mL **methanol** and **sonicate** ( 00:15:00).

15m

8 **Centrifuge** ( 865 x g, 00:15:00) to separate cell pellets from supernatant

15m



9 **Transfer the supernatant** to a new 96-well plate and store at  -80 °C if the extraction will not be performed immediately. Also, if the cell pellets will be extracted separately, they can be stored at  -80 °C until extraction.



10 Prepare the 96-well **SPE** plate before performing the extraction. See the SPE steps described below.

11 The following are the steps for **SPE**:

Note

Specific protocols can vary from providers, so follow recommendations from the providers of the purchased product: Phenomenex, Sigma-Aldrich or Waters.

Briefly, the steps are:

11.1 **Condition** the cartridge with methanol ( 2 mL).



- 11.2 **Equilibrate** with a similar solvent as the sample. In this case, water ( 2 mL) as the supernatant is aqueous.
- 11.3 **Load the sample** ( 1 mL of supernatant is fine for a 100mg SPE cartridge).
- 11.4 **Wash** the sample ( 1 mL of water).
- 11.5 **Elute** ( 1 mL of methanol).
- 11.6 Optional: **recondition** cartridges for another SPE.
- 12 Now the extracts ( 1 mL of methanol containing the extracted metabolites) are in a 96-deep well plate.
- 13 You can **dry them** under vacuum. Dried samples will be suitable for shipping to the Dorrestein Lab.

Note

Endometabolome: extraction procedure performed on filtered microbial cells (pellets) [For optional protocols, such as for *B. subtilis*, see [Meyer et al., 2013](#). This review [Pinu et al., 2017](#) might be also helpful];

- 14 Include, together with the samples, **blanks of the culture media** used in each experiment (extracted under exact conditions as the exometabolome).
- 15 At least **biological replicates (3-6)** from each sample. We run technical replicates for each sample.
Ship **dry extracts** if possible.
- 16 If shipping standards together with the microbial extracts is possible, please include them as well as their information in the metadata table.
- 17 If possible, use Deepwell Plates, 96-well, Polypropylene, 2000µL [e.g. [E951033600](#)]. Include the map



of the samples' name/label.

Protocol references

References and additional information:

<https://www.embl.org/groups/metabolomics/protocols-used-for-lc-ms-analysis/>

<https://redu.gnps2.org/selection/>

<https://onlinelibrary.wiley.com/doi/full/10.1002/rcm.8725>

<https://www.mdpi.com/2218-1989/7/3/43>

<https://microbialcellfactories.biomedcentral.com/articles/10.1186/1475-2859-12-69>

<https://www.mdpi.com/2218-1989/7/4/53>

Acknowledgements

This work was supported by the Gordon and Betty Moore Foundation, GBMF12120 and <https://doi.org/10.37807/GBMF12120>.