Nov 19, 2025

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Annie Archambault

CIEL

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

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<https://dx.doi.org/10.17504/protocols.io.36wggpk8ovk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: November 10, 2025

Last Modified: November 19, 2025

Protocol Integer ID: 232042

Keywords: fruit juice from kiwiberry, fruit par centrifugation et, extracting fruit juice, cranberry for brix determination, ml of fruit juice, préparation simple mais efficace de fruits charnus afin, robust preparation of fleshy fruit, fruit paste with centrifugation, fruit juice, simple mais efficace de fruits charnus afin, kiwiberry, fruit en environ, extraire le jus de fruit, fruit, juice, fleshy fruit, fruit paste, frozen juice, kiwi et, cranberry, initialement développé avec des kiwis et de, brix determination, des kiwis et de, détermination du brix, par le nombre de puits du rotor, votre centrifugeus, centrifuge, centrifugation, la pulpe

Funders Acknowledgements:

Programme Innovation bioalimentaire MAPAQ

Grant ID: T-09943-F7H9

Abstract

This protocol describes the simple yet robust preparation of fleshy fruits to enable the collection of sufficient volume of juice from fresh or frozen juice by sedimenting the fruit paste with centrifugation and collecting the the supernatant.

The protocol was initially developed with kiwiberry and cranberries by CIEL <https://www.ciel-cvp.ca/> in L'Assomption, Quebec.

The protocol is run by batches (only limited by the number of wells in your rotor's centrifuge) and should yield between 200 ul and 1 ml of fruit juice in approximately 2 hours.

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Materials

 Distilled water **Fisher Scientific**

Delicate Task Wipers

Pipette, single channel, 200 ul

Pipette tips 200 ul

50 ml centrifuge tubes

Centrifuge with rotor compatible for 50 ml tubes and with capacity of up to 2000 RPM

Equipment

HI96811	NAME
Refractometer	TYPE
Hanaa	BRAND
HI96811	SKU
https://hannacan.com/product/digital-refractometer-for-sugar-brix-analysis-in-wine-must-and-juice-hi96811/	LINK

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