

Jul 10, 2024

Ezrin plasmids

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwjww7lmk/v1>

Shiyi Wang¹

¹Duke University



Shiyi Wang

Duke University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwjww7lmk/v1>

Document Citation: Shiyi Wang 2024. Ezrin plasmids. protocols.io
<https://dx.doi.org/10.17504/protocols.io.e6nvwjww7lmk/v1>

License: This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Created: May 24, 2023

Last Modified: July 10, 2024



Document Integer ID: 82338

Keywords: ASAPCRN, ezrin plasmid

Funders Acknowledgements:

Aligning Science Across Parkinson's (ASAP) initiative

Grant ID: ASAP-020607

Abstract

How Ezrin plasmids were made.

1. pZac2.1-GfaABC1D-Lck-GCaMP6f was a gift from Dr. Baljit Khakh (Addgene plasmid #52924).
2. pZac2.1-GfaABC1D-BioID2-HA was generated by PCR of BioID2 from pAAV-hSyn-BioID2-Linker-Synapsin1a-HA (primers Fw: 5'-ctagcctcgagaattcaccatgttcaaaaatcttatttg-3', Rv: 5'-ccgggtcgactctagatgcgtaatccggtacatcg-3') and insertion into the EcoRI and XbaI restriction sites of pZac2.1-GfaABC1D-Lck-GCaMP6f using In-Fusion cloning (TaKaRa).
3. pHJ421(pEGFP-Ezrin WT) and pHJ423 (pEGFP-Ezrin T567D) were a gift from Stephen Shaw (Addgene plasmid # 20680 and # 20681).
4. pZac2.1-GfaABC1D-Ezrin WT-BioID2-HA and pZac2.1-GfaABC1D-Ezrin T567D-BioID2-HA were generated by PCR of Ezrin from pHJ421 or pHJ423 (primers Fw: 5'-ctagcctcgagaattcaccatgccgaaaccaatca-3', Rv: 5'-tgaacatggtgaattccgacagggcctcgaactcg-3') and insertion into the EcoRI restriction sites of pZac2.1-GfaABC1D-BioID2, respectively.
5. To generate pZac2.1-GfaABC1D-Ezrin T567A-BioID2-HA plasmid, Q5® Site-Directed Mutagenesis Kit (NEB) was used with mutagenesis primers (Fwd: CAAGTACAAGGCGCTGCGGCAGA; Rev: TCCCGGCCTTGCCTCATG)