

Jun 12, 2019

## f/2 medium with soil extract

DOI

[dx.doi.org/10.17504/protocols.io.s29egh6](https://dx.doi.org/10.17504/protocols.io.s29egh6)



Roscoff Culture Collection<sup>1</sup>

<sup>1</sup>Station Biologique, CNRS-Sorbonne Université, Roscoff France

Roscoff Culture Collection



**Daniel Vaultot**

Station Biologique, Roscoff, France , Nanyang Technologica...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.s29egh6>

**Protocol Citation:** Roscoff Culture Collection 2019. f/2 medium with soil extract. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.s29egh6>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 28, 2018

**Last Modified:** June 12, 2019

**Protocol Integer ID:** 15169

**Keywords:** dinoflagellate, soil extract, phytoplankton species, soil

## Abstract

Medium to grow phytoplankton species in particular dinoflagellates which often require soil extract.

## Troubleshooting

## Before start

Please refer to our general recommendations to grow cultures :

<https://www.protocols.io/private/A48906DC1374AD6281495CB86A8F092F>

## Prepare f/2 medium

1

### Protocol



NAME

f/2 medium

CREATED BY

Daniel Vaultot

[Preview](#)

- 1.1
- Filter 1L of old seawater of at least two months on prefilter and 0.2  $\mu\text{m}$  filter
  - Heat seawater during 10/20 min at 100°C (pasteurisation)
  - Let the water cool down
  - Under laminar flow hood, , add 25 mL of Guillard's (f/2) Marine Water Enrichment Solution (Sigma-Aldrich G9903)
  - Filter the medium on 0,2  $\mu\text{m}$
- Guillard's (F/2) Marine Water Enrichment Solution **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G0154-500ML**
- 1.2 See: [https://ncma.bigelow.org/media/wysiwyg/Algal\\_recipes/NCMA\\_algal\\_medium\\_f\\_2\\_1.pdf](https://ncma.bigelow.org/media/wysiwyg/Algal_recipes/NCMA_algal_medium_f_2_1.pdf)

## Add soil extract

2

### Protocol



NAME

Soil extract for algal media

CREATED BY

Daniel Vaultot

[Preview](#)

- 2.1 Sample 10 g of dry soil from a place that does not contain any pesticide nor pollutant and where you can be sure to be able to come back in order to always use the same soil. The



soil must be dry because the boiling step will be longer if the soil is not dry.

2.2 Add the 10 g soil to 400ml of milliQ water.

2.3 Boil during 1 hour.

2.4 Filter on 0.2  $\mu\text{m}$  or first through 0.8/0.4  $\mu\text{m}$  and then through 0.2 $\mu\text{m}$  in order not to clog the filter.

2.5 Aliquot and freeze.

2.6 Add 30 à 50ml of soil extract per 1 L of medium.