



Sep 04, 2024

FACS Nuclei Isolation

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2qbwql1y/v1>

Naeem Nadaf¹

¹Broad Institute of MIT and Harvard



Emily Finn

Broad Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2qbwql1y/v1>

Protocol Citation: Naeem Nadaf 2024. FACS Nuclei Isolation . [protocols.io](https://dx.doi.org/10.17504/protocols.io.4r3l2qbwql1y/v1)

<https://dx.doi.org/10.17504/protocols.io.4r3l2qbwql1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 14, 2024



Last Modified: September 04, 2024

Protocol Integer ID: 105343

Keywords: facs nuclei isolation, isolation of nuclei, such as 10x genomic, minimizing rna degradation, 10x genomic, integrity of nuclei, rna, rna degradation, sequencing, nuclei, frozen tissue, droplet

Disclaimer

This protocol provides a meticulous approach to extracting and processing nuclei from frozen tissues, ensuring high-quality samples suitable for droplet-based single-nuclei sequencing.

Abstract

This protocol details the isolation of nuclei from frozen tissue for subsequent single-nuclei sequencing using a droplet-based platform, such as 10X Genomics. The procedure emphasizes maintaining the integrity of nuclei and minimizing RNA degradation, using a series of specialized buffers and careful handling techniques.

Nuclei Extraction protocol, optimized for small tissue pieces

1 Prepare Reagents and Buffers

Note

ALL buffers to be cooled to 4 °C before use

1.1 Dissociation Buffer (DB) Preparation: Dissolve these components in ultrapure, nuclease-free water to a final volume of 500 mL

- 5.83 g Na₂SO₄
- 2.615 g K₂SO₄
- 0.905 g Glucose
- 1.2 g HEPES
- 2.5 mL MgCl₂

1.2 Extraction Buffer (ExB):

- 15 mL DB (Dissociation Buffer)
- 150 mg Kollidon VA64
- 150 µL TX-100 (final concentration 1%)
- 15 µL 10% BSA (final concentration 0.01%)
- RNase inhibitor: 1 tube (Lucigen, 10,000 units)

1.3 Wash Buffer (WB): Prepare 30 mL per sample

- 30 mL DB (Dissociation Buffer)
- 30 µL 10% BSA
- 50 µL of RNase inhibitor (Lucigen, 10,000 units)

1.4 FACS Capture Buffer (CDB):

- 20 mL Dissociation Buffer (DB)
- 20 µL 10% BSA
- RNase inhibitor: 1 tube (Lucigen, 10,000 units)

1.5 Prepare: 5% BSA-DB for FACS (Prepare 200 µL per sample)

2

Procedure**2.1 Pre-Cooling:**

- Ensure all buffers are cooled to 4 °C before use
- Pre-cool centrifuge, tubes, well plates, 26-gauge needle, 40 µm cell strainer, and syringe at 4 °C for at least 20 minutes.

2.2 Tissue Dissociation:

- Transfer the frozen tissue sample to the prepared 6-well plate using 150 µL of ExB buffer.
- Triturate by pipetting the 1 mL volume slowly up and down using a 1 mL Rainin pipette tip (#30389212) without creating froth or bubbles. Perform 25-30 strokes, pause for 2 minutes, and repeat 4-5 times. Monitor progress; add a fifth pass only if dissociation appears incomplete.
- Pass the entire volume through a 26-gauge needle twice into the same well.

2.3 Sample Processing:

- Transfer the extracted sample (~1 mL) into a 50 mL Falcon tube. Add 30 mL of wash buffer to dissociate the nuclei, then split the volume between two pre-coated Falcon tubes (15+ mL each)
- Centrifuge at 600xg for 10 minutes at 4°C.

2.4 Supernatant Removal:

- Aspirate the supernatant, leaving 500 µL in each tube. Pool the two samples to obtain 1 mL of nuclei suspension.
- Filter the sample through a pre-cooled 40 µm cell strainer using gravity only (no pressure) into a new, clean, pre-cooled Falcon tube.

2.5 Staining:

- Transfer and measure the volume of the filtrate into a pre-cooled 1.5 mL Eppendorf tube.
- Stain the nuclei by adding DAPI (Thermo, #62248) at a 1:1000 dilution.

2.6 FACS Preparation:

- Coat a 0.2 mL PCR tube with 5% BSA. Use a chilled 96-well FACS plate (Sony M800 FACS sorter) to capture nuclei in the PCR tube. Pre-fill the PCR tube with 20 µL of CDB to cushion the FACS stream/nuclei.

2.7 FACS Sorting:

- Perform FACS enrichment for singlets using the singlet DAPI peak. Sort at speed 6 with a forward scatter gain of 1% on the DAPI gate. Use the "purity" mode. Avoid spinning down after FACS unless using the gentle spin protocol.



Nuclei FACS enrichment for singlets on singlet DAPI peak

3



Note

We FACS on a Sony SH800 sorter, using a 70um chip, with these settings:

FACS Settings:

- **Instrument:** Sony SH800/SH900 cell sorter
- **Chip:** 70 μ m chip
- **Gating:** Focus on the singlet DAPI peak to enrich for singlet nuclei.
- DAPI peak: y axis = events, x axis = brilliant violet 421-A
- DAPI gate on **purity mode**
- Sample pressure of **3-6**
- Cell size on **regular**
- Sort rate at **4eps**

Sample Group - 1 - Measurement Settings

Parameter Settings

	Marker	Fluorochrome	Acquisition Select		
			Area	Height	Width
FSC:	---	---	☒	☒	☒
BSC:	---	---	☒	☐	☐
FL1:		Brilliant Violet 421 v	☒	☐	☐
FL2:		FITC v	☒	☐	☐
FL3:		PE v	☒	☐	☐
FL4:		APC v	☒	☐	☐
FL5:		PerCP-Cy5.5 v	☒	☐	☐
FL6:		PE-Cy7 v	☒	☐	☐

If this Sample Group has a recorded tube, you cannot change Acquisition Select Settings.

Instrument Settings

Laser:

405nm: On

488nm: On

561nm: On

638nm: On

Threshold:

Channel: FSC

Value: 1.00%

Sensor Gain:

FSC: 4

BSC: 27.0%

FL1: 40.0%

FL2: 40.0%

FL3: 40.0%

FL4: 40.0%

FL5: 40.0%

FL6: 40.0%

Sample Pressure: 6


AD Advanced Settings:

Forward Window Extension: 50

Back Window Extension: 50

Note

We gate on the singlet DAPI peak:



4 Collection of FACS-Sorted Samples:

To ensure the integrity and minimal collection volume of FACS-sorted nuclei, samples are collected in PCR tubes positioned within a chilled (-20°C) 96-well cold block. This setup is essential to maintain nuclei quality and prevent degradation. A suitable cold block can be found here: Eppendorf PCR Coolers.

Determination of Nuclei Concentration Post-FACS:

5 1. Volume Measurement:

- Use a pipette to accurately determine the volume of the collected nuclei suspension in each PCR tube post-FACS.

2. Dilution Preparation:

- Prepare a 1:10 dilution by combining 18 μL of chilled Dissociation Buffer (DB) with 2 μL of the nuclei suspension in a PCR tube. Mix thoroughly.

3. Counting Nuclei:

- Load the diluted sample onto a Fuchs-Rosenthal hemocytometer (16 chambers).
- Visualize using a fluorescent microscope, capturing images under both brightfield and DAPI excitation to ensure debris removal and accurate nuclei identification.
- Count the nuclei across all 16 large squares for precise concentration determination.



- Calculate the average nuclei count per square, multiply by 10 (to account for the 1:10 dilution), and then multiply by 5 (the hemocytometer factor) to determine the final nuclei concentration in nuclei/ μL .

Proceeding to the 10X Genomics Protocol:

- 6
 - For the 10X Genomics v3 workflow, the input volume is 46.6 μL . To maximize data output, the optimal nuclei concentration should be 364 nuclei/ μL , translating to an input of approximately 17,000 nuclei per 46.6 μL next gem kit.

This protocol provides a meticulous approach to extracting and processing nuclei from frozen tissues, ensuring high-quality samples suitable for droplet-based single-nuclei sequencing.