

Apr 09, 2025

FAMES analysis

 [Nature Cancer](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8e5b7l2w/v1>

Laura Antonucci¹, David A. Scott², Michael Karin¹

¹University of California, San Diego;

²Cancer Metabolism Core, Sanford Burnham Prebys Medical Discovery Institute, La Jolla CA 92037, USA



Laura Antonucci

University of California, San Diego

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8e5b7l2w/v1>

External link: <https://doi.org/10.1038/s43018-025-01003-3>

Protocol Citation: Laura Antonucci, David A. Scott, Michael Karin 2025. FAMES analysis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8e5b7l2w/v1>

**Manuscript citation:**

Antonucci, L., Li, N., Duran, A. *et al.* Self-amplifying NRF2–EZH2 epigenetic loop converts KRAS-initiated progenitors to invasive pancreatic cancer. *Nat Cancer* **6**, 1263–1282 (2025). <https://doi.org/10.1038/s43018-025-01003-3>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 09, 2025

Last Modified: April 09, 2025

Protocol Integer ID: 126404

Keywords: assessment of de novo fatty acid synthesis, fames analysis fatty acid methyl ester, carbon incorporation into fatty acid, de novo fatty acid synthesis, fatty acid, pec organoid, organoid, glucose, derived carbon incorporation

Abstract

Fatty acid methyl ester (FAMES) analysis was used to measure glucose-derived carbon incorporation into fatty acids in *KrasG12D/PEC* organoids. Organoids were labeled with [U-¹³C₆] glucose, extracted, derivatized, and analyzed by GC-MS. Quantification was done using internal standards and MetaQuant software, allowing assessment of de novo fatty acid synthesis.

FAMES analysis

- 1 Cell labeling with [U-¹³C₆] glucose was performed as previously described¹ with some modification for the organoid cultures. A total of 3×10^5 *Kras*^{G12D/PEC} cells for each experimental point were plated in 96 96-well plate, in a dome 1:1 of Matrigel (Corning, cat #356231) and DMEM no glucose (Gibco, 11966-025). Organoids were treated and cultured as described above for 3 days. At the start of the labeling period, (t) = 0 h, 4 days after plating in Matrigel, the media was replaced with a growth media containing DMEM no glucose, 1/3 Advanced DMEM/F-12 supplemented with 1 g/L [U-¹³C₆] glucose (Sigma-Aldrich). After 24 hours, organoids were harvested with Cell Recovery Solution (Corning, 354253), as described previously^{2, 3}. Cells or organoids were extracted with 0.45 ml 1:1 water: methanol plus 0.225 ml chloroform containing 20 μM D31-palmitic acid (Cayman). Samples were dried using a Speedvac centrifugal evaporator. Blank tubes were included to control for contaminant fatty acids in tubes.
For quantification, FAMES-37 (Sigma) standards mixture (50 μl) was mixed with 225 μl 200 μM D31-palmitic acid. Varying amounts of this mixture were dispensed and dried. HCl 0.5N in methanol (50 μl) was added to dried samples or standards, and these were heated for 30 min at 60°C, before drying and resuspension (with heating to 60°C) in 50 μl pyridine and transfer to autosampler vials with inserts. GC-MS: Samples and standards were run on a Shimadzu QP2010 Plus GC-MS. Injection volume was 1-2 μl with 1/10 split, injection temperature of 250°C; GC oven temperature was initially 80°C for 3.5 min, rising to 280°C at 10°C/min with a final hold at this temperature for 2 min. GC flow rate with helium carrier gas was 50 cm/s. The GC column used was a 15 m x 0.25 mm x 0.25 μm RXI-5ms (Restek). GC-MS interface temperature was 300°C and (electron impact) ion source temperature was 200°C, with 70 eV/ 70 μA ionization voltage/current. The mass spectrometer was set to scan m/z range 100-500, with varied detector sensitivity. Quantification was done against standard curves constructed in Metaquant⁴, and amounts were adjusted for recovery of the internal standard. Correction for isotopic natural abundance and calculation of acetyl unit (acetyl-CoA) fractional labeling and fatty acid synthesis is done in MS Excel, using raw intensity m/z data ranging between unlabeled and fully labeled ions for palmitic acid and other fatty acids¹.

Protocol references

1. Scott, D.A. *et al.* Comparative metabolic flux profiling of melanoma cell lines: beyond the Warburg effect. *J Biol Chem* **286**, 42626-42634 (2011).
2. Boj, S.F. *et al.* Organoid models of human and mouse ductal pancreatic cancer. *Cell* **160**, 324-338 (2015).
3. Huch, M. *et al.* Unlimited in vitro expansion of adult bi-potent pancreas progenitors through the Lgr5/R-spondin axis. *EMBO J* **32**, 2708-2721 (2013).
4. Bunk, B. *et al.* MetaQuant: a tool for the automatic quantification of GC/MS-based metabolome data. *Bioinformatics* **22**, 2962-2965 (2006).