

Oct 28, 2025

FASP (for fChEP samples)

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5rry6g1b/v1>



Victoria Munro¹

¹University of Edinburgh



Victoria Munro

University of Edinburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5rry6g1b/v1>

Protocol Citation: Victoria Munro 2025. FASP (for fChEP samples). protocols.io

<https://dx.doi.org/10.17504/protocols.io.8epv5rry6g1b/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 06, 2024

Last Modified: October 28, 2025

Protocol Integer ID: 101329

Keywords: ChEP, FASP, fchep sample, fasp, sample preparation, fchep, sample, preparation, containing sd, filter

Abstract

Filter aided sample preparation (FASP) clean-up protocol optimised for fChEP samples containing 1% SDS.

Guidelines

Day 1 spin time = 3 hours (+ 40 min for reduction and alkylation)

Day 2 spin time = 40 minutes

ABC per sample = 500ul

Urea/ABC per sample = 1ml

Materials

Vivacon 500, Hydrosart, 30,000 MWCO (APPLETON WOODS LTD, VN01H22)

Urea (Sigma-Aldrich, U5378-500G)

Ammonium bicarbonate (Sigma-Aldrich, A6141-500G)

Dithiothreitol (DTT; Promega, V3151)

Iodoacetamide (IAA; ?)

Trypsin (Thermo, 90058)

Safety warnings

 Note: this protocol has 2 extra urea/ABC washes compared to the original and is for samples which have not yet been reduced and alkylated. The number of washes will be sample and buffer dependent and may need optimising.

Sample preparation

- 1 Start with cell lysate with sheared DNA and of known protein concentration.
- 2 Heat at 95°C for 5 min.
- 3 Spin at 16,000 x g for 5 min to clarify the lysate

Load to filter, wash with urea, reduction and alkylation

- 4 Combine pre-determined amount of sample with 8M urea, 50 mM ammonium bicarbonate in the 30kDA MWCO filter unit so that total volume is 200ul. Spin at 14,000 rcf for 20 minutes.
- 5 Add 200 µl 8M urea, 50 mM ammonium bicarbonate to filter, and spin at 14,000 rcf for 20 minutes. Discard flow through
- 6 Add 100 µl 10 mM DTT in 8M urea, 50 mM ammonium bicarbonate. Incubate at room temperature for 20 minutes
- 7 Then add 100 ul 55mM IAA (27.5 mM final) in 8M urea, 50mM ABC. Incubate in dark at room temperature for 20 minutes.
- 8 Spin 14,000 rcf for 20 minutes
- 9 Add 100 µl 8M urea, 50 mM ABC to filter. Spin at 14,000 rcf for 20 minutes
- 10 100 µl 8M urea, 50 mM ammonium bicarbonate, and spin again
- 11 100 µl 8M urea, 50 mM ammonium bicarbonate, and spin again
- 12 100 µl 8M urea, 50 mM ammonium bicarbonate, and spin again



Deplete urea

13 Add 100 μ l of 50 mM ammonium bicarbonate. Spin at 14,000 rcf for 20 minutes

14 Repeat 100 μ l 50 mM ammonium bicarbonate, and spin again

Digestion

15 Transfer the filter unit to a new collection tube

16 Add 1/20 trypsin in 100ul 50mM ABC to each sample. Parafilm tube lids shut. Digest 4 - 18 hrs (usually 16 h) in 37°C incubation hood with wet paper towel.

Peptide elution

17 Collect filtrate by spinning at 14,000 rcf for 20 minutes. If dried out, add 100ul 50mM ABC first. (Instead of 100 ul add 50 ul of 50 mM ABC so that the final volume isn't that high).

18 Add 100uL of ABC and spin again for 20min @ 14,000 rcf. Collect again.

19 Proceed to C18