

Oct 24, 2025

 FASP-IP

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7156kgwz/v1>



Cristina CARDENAL PERALTA¹

¹University of Edinburgh



Cristina CARDENAL PERALTA

University of Edinburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7156kgwz/v1>

Protocol Citation: Cristina CARDENAL PERALTA 2025. FASP-IP. protocols.io

<https://dx.doi.org/10.17504/protocols.io.q26g7156kgwz/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 24, 2024

Last Modified: October 24, 2025

Protocol Integer ID: 102307

Keywords: FASP, CO-IP, Tryptic Digestion, efficient protein digestion, membrane filter, using membrane filter, fasp, mass spectrometry, aided sample preparation, filter, ip the filter, contaminant, eliminating detergent

Abstract

The Filter-Aided Sample Preparation (FASP) protocol enables efficient protein digestion and cleanup using membrane filters, eliminating detergents and contaminants prior to mass spectrometry.

Guidelines

Ensure all liquid has passed through the filter before proceeding to the next step. The membrane should remain visibly moist — it must have a wet sheen but never be completely dry. Drying can cause protein loss, membrane damage, or reduced recovery during subsequent washes and digestion steps.

Materials

Solutions:

DTT solution: 25 mM DTT in 50 mM ABC.

Precipitation buffer: 8M urea in 50 mM ABC

IAA solution: 55 mM IAA in 8M urea in 50 mM ABC.

ABC solution: 50 mM ABC

Materials:

FASP tubes: VN01H22 Sartorius

Safety warnings

- ! ABC must be made fresh the same day of the experiment.
Urea will precipitate if exposed to high temperatures and prevent carbamylation



Reduction/Alkylation

1h 5m

- 1 Add DTT solution to the sample for a final concentration of 25 mM. Incubate at 80°C for 15 min. 15m
- 2 Wait for the sample to cool down. Add urea to the sample to have a final concentration of 8M.
This can be either directly adding the urea powder to the solution, or diluting 5x with 8M urea in 50 mM ABC. 15m
- 3 Spin sample through a 30k spin column until all the liquid goes through the filter, 15000r cf for 15. 15m
- 4 Add 100 µL of 55 mM IAA to the column. Incubate at  Room temperature for 20 min in the dark. 20m
- 5 Spin through column at 15000 rcf for 15min. 15m

Wash

45m

- 6 Wash with 100 µL of 8 M urea in ABC. Spin through column at 15000 rcf for 15-20 min. 15m
- 7 Wash twice with 100 µL of ABC. Spin through column at 15000 rcf for 15min. 30m

Digestion

- 8 Make trypsin solution by adding 20 µL of 0.1 % TFA to a trypsin vial. Add 980 µL of ABC.
- 9 Place column into a new tube and add 100 µL of freshly prepared trypsin solution. Wrap the top of the tubes with parafilm and incubate at 37 °C O/N. Use a wet paper towel underneath the Eppendorf rack to minimise evaporation. 