

Oct 24, 2025

FASP-No IP

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk8d26l5r/v1>



Cristina CARDENAL PERALTA¹

¹University of Edinburgh



Cristina CARDENAL PERALTA

University of Edinburgh

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk8d26l5r/v1>

Protocol Citation: Cristina CARDENAL PERALTA 2025. FASP-No IP. protocols.io

<https://dx.doi.org/10.17504/protocols.io.j8nlk8d26l5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: June 24, 2024

Last Modified: October 24, 2025

Protocol Integer ID: 102301

Keywords: FASP, Mass Spectrometry, efficient protein digestion, fasp, membrane filter, using membrane filter, mass spectrometry, aided sample preparation, filter, no ip the filter, ip the filter, contaminant, eliminating detergent

Abstract

The Filter-Aided Sample Preparation (FASP) protocol enables efficient protein digestion and clean-up using membrane filters, eliminating detergents and contaminants prior to mass spectrometry.

Guidelines

Ensure all liquid has passed through the filter before proceeding to the next step. The membrane should remain visibly moist — it must have a *wet sheen* but never be completely dry. Drying can cause protein loss, membrane damage, or reduced recovery during subsequent washes and digestion steps.

Materials

Buffers:

Precipitation buffer: 8M urea in 50 mM ABC

DTT solution: 10 mM DTT in 8M urea in 50 mM ABC

IAA solution: 55 mM IAA in 8M urea in 50 mM ABC

ABC: 50 mM of ammonium bicarbonate

Materials:

FASP tubes: VN01H22 Sartorius

Safety warnings

! Always avoid urea > 40 °C (prevents carbamylation).



Sample reception

15m

- 1 Combine 30 μg worth of protein with 150 μL of precipitation buffer in the filter unit.
- 2 Spin at 15000 rcf for 15 min
- 3 Add 200 μL of precipitation buffer to the filter unit. Spin at 15000 rcf for 15 min.

Reduction and Alkylation

15m

- 4 If the sample was lysed without DTT, add 100 μL of DTT solution. Incubate at  Room temperature for 20min.
- 5 Add 100 μL of IAA solution. Incubate in the dark for 20 min.
- 6 Spin at 15000 rcf for 15 min.

Washes

1h

- 7 Add 100 μL of precipitation buffer to filter, spin at 15000 rcf for 15 min.
- 8 Add 100 μL of 50 mM ABC to filter, spin at 15000 rcf for 15 min.
These washes serve the purpose of eliminating urea from the sample, as trypsin is not active in high urea concentrations.
- 9 Repeat previous step 2 times. (If the sample was lysed with SDS, increase the number of washes by two or three).
- 10 Transfer the filter unit to a new collection tube. Label both tube and filter.

Digestion



- 11 Add 20 μL of 0.1 % TFA to the trypsin vial, then add 980 μL of ABC.
- 12 Add 100 μL of the newly made trypsin solution to each filter unit. Incubate O/N at 37 °C. Wrap each tube with parafilm and put a wet paper towel underneath the rack to ensure the incubation environment remains moist.
- 13 Collect filtrate by spinning at 15000rcf for 15 minutes.
- 14 Add 100ul of ABC and spin again for 15 min at 15000 rcf. Collect again.



Sample acidification

- 15 Add 20 μL of 10 % TFA to acidify the samples. Proceed to StageTips