



Aug 02, 2022

FAST Dil Injection Protocol from Neural Tracing in Porcine SAN

DOI

<https://dx.doi.org/10.17504/protocols.io.4r3l2okqqv1y/v1>

Donald Hoover¹, Peter Hanna², Zixi (Jack) Cheng²

¹East Tennessee State University; ²University of California, Los Angeles

SPARC

Tech. support email: info@neuinfo.org



Chloe Garbe

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2okqqv1y/v1>

Protocol Citation: Donald Hoover, Peter Hanna, Zixi (Jack) Cheng 2022. FAST Dil Injection Protocol from Neural Tracing in Porcine SAN . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2okqqv1y/v1>

**Manuscript citation:**

Hanna P, Dacey MJ, Brennan J, Moss A, Robbins S, Achanta S, Biscola NP, Swid MA, Rajendran PS, Mori S, Hadaya JE, Smith EH, Peirce SG, Chen J, Havton LA, Cheng ZJ, Vadigepalli R, Schwaber J, Lux RL, Efimov I, Tompkins JD, Hoover DB, Ardell JL, Shivkumar K. Innervation and Neuronal Control of the Mammalian Sinoatrial Node a Comprehensive Atlas. Circ Res. 2021 Apr 30;128(9):1279-1296. doi: 10.1161/CIRCRESAHA.120.318458. Epub 2021 Feb 25. PMID: 33629877; PMCID: PMC8284939.

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 18, 2022

Last Modified: August 02, 2022

Protocol Integer ID: 66991

Keywords: fast dii injection protocol from neural tracing, fast dii into the sinoatrial node, fast dii injection protocol, porcine test subject, injecting fast dii, neural tracing, sinoatrial node, porcine san, tracing purpose

Funders Acknowledgements:

Kalyanam Shivkumar

Grant ID: OT2 OD023848

Abstract

Protocol for injecting FAST Dil into the sinoatrial node of porcine test subjects for neural tracing purposes.

Materials

FAST Dil from Invitrogen, catalog D7756

CitiFluor AF1 Mounting Medium Solution (Electron Microscopy Sciences, Hatfield, PA, USA, Cat. #17970-100

Injection of FAST DiI

- 1 Sedated, anesthetized, and intubated pig (induction: ketamine 10 mg/kg IM/midazolam 1 mg/ kg IM, maintenance: isoflurane 1-2% inhalation).
- 2 Performed right unilateral thoracotomy by dividing the pectoral muscle.
- 3 Made a small incision in the pericardium and exposed the right atrial (RA)-superior vena cava junction.
- 4 Used 27-gauge needle to inject 5 mg of an analog of 1,1'-dioctadecyl-3,3,3',3'-tetramethylindocarbocyanine perchlorate (FAST DiI, Invitrogen, D7756) in 250 μ L of 100% methanol (MeOH) into the sinoatrial nodal (SAN) region.
- 5 Placed chest tube and closed incision.
- 6 Aspirated chest tube immediately prior to its removal.

Tissue Dissection and Preparation

- 7 Harvested tissue in a terminal procedure at least 3 weeks after injection procedure.
- 8 Immersed tissue in cold 4% PFA for 24 hours.
- 9 Washed tissue in PBS 4 times.
- 10 Transferred tissue to cold 0.01M PBS + 0.02% Na azide + 20% sucrose solution for shipment to East Tennessee State University.
- 11 Dissected tissue as needed to isolate specific regions.
- 12 Froze tissues on dry ice and stored at -80°C until ready to be sectioned.



Tissue Sectioning Using the Cryostat

- 13 Removed tissues from -80°C freezer and mounted onto specimen plate using Tissue-Tek® O.C.T. Compound (Sakura Finetek USA, Cat. No. 4583).
- 14 Cut tissues into 30µm sections at -20 to -25°C using a Leica CM3050S cryostat (Leica Microsystems Inc., Bannockburn, IL, USA).
- 15 Collected sections on charged slides in a sequence that yields at least 4 sets of tissue sections per region of the heart, with each set spanning the entire thickness of the specimen.
- 16 Stored sets of tissues in slide boxes wrapped in aluminum foil at -20°C until further processing.

Immunohistochemistry

- 17 Removed every second slide from a set to perform immunolabel for PGP9.5. (Described in separate protocol, found here: dx.doi.org/10.17504/protocols.io.8rmhv46)
- 18 Applied cover glass to adjacent tissue slides using a drop of CitiFluor AF1 Mounting Medium Solution (Electron Microscopy Sciences, Hatfield, PA, USA, Cat. #17970-100).

Imaging

- 19 Imaged tissue sections using fluorescence illumination with an Olympus BX41 microscope equipped with an Olympus DP74 digital camera and cellSens software (Olympus America Inc., Center Valley, PA; RRID:SCR_016238).
- 20 Loaded slide onto stage and collected representative images of ganglia.

Quantification

- 21 Counted Dil-labeled neurons and total number of neurons (i.e., PGP9.5-positive).