

Dec 16, 2024

Fatty Acid Induction to Hippocampal Neurons

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj93jxlk5/v1>

Mukesh Kumar¹, Timothy A. Ryan^{1,2}

¹Department of Biochemistry, Weill Cornell Medicine, New York, NY 10065;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, Maryland 20815, USA



Mukesh Kumar

Weill Cornell Medical College

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj93jxlk5/v1>

Protocol Citation: Mukesh Kumar, Timothy A. Ryan 2024. Fatty Acid Induction to Hippocampal Neurons. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj93jxlk5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: November 27, 2024

Last Modified: December 16, 2024

Protocol Integer ID: 114387

Keywords: ASAPCRN, fatty acid induction to hippocampal neuron, acid induction to primary hippocampal neuron, fatty acid induction, bsa complexes for lipid droplet induction, fatty acid, lipid droplet induction, hippocampal neurons this protocol, preparing homogeneous fatty acid, hippocampal neuron, primary cultured hippocampal neuron, primary hippocampal neuron, homogeneous fatty acid, acid induction, neuron, bsa complex

Funders Acknowledgements:

NIH

Grant ID: NS036942

NIH

Grant ID: NS11739

ASAP

Grant ID: ASAP-024404

Abstract

This protocol provides a stepwise procedure for preparing homogeneous fatty acid–BSA complexes for lipid droplet induction in primary cultured hippocampal neurons. The protocol also describes the fatty acid induction to primary hippocampal neurons for electrical stimulation experiments.

Guidelines

The protocol needs prior approval by the users' Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee.



Materials

Materials:

-  Oleic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #O1383**
-  Palmitic acid **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P0500**
-  BSA **Merck MilliporeSigma (Sigma-Aldrich) Catalog ##A8806**
- Tris Buffer: 0.1 M Tris (pH 7.4 using NaOH)
- NaOH Solution: 0.1 M NaOH (200 mg NaOH in 50 ml ultrapure water)
- NaCl Solution: 0.9% NaCl (450 mg NaCl in 50 ml ultrapure water)
-  BODIPY[®] 493/503 (4,4-Difluoro-1,3,5,7,8-Pentamethyl-4-Bora-3a,4a-Diaza-s-Indacene) **ThermoFisher Catalog #D3922**

Equipments:

- Heating block (temperature controlled to 40°C and 70°C as required)
- Syringe filter: 0.22 µm pore size
- Culture hood
- Aluminum foil
- Ultrapure water
- Vortex heater



Procedure to prepare Oleic Acid–BSA (OA-BSA) Complex

1d 3h 10m

1 Prepare 0.1 M Tris Buffer (pH 7.4):

Dissolve Trizma in ultrapure water and adjust pH to 7.4 using NaOH.

2 Dissolve BSA:

2.1 Pour  4 mL of [M] 0.1 Mass Percent Tris into a 50 ml tube. Overlay  697.5 mg of FA-free BSA on top.

2.2 Place the tube on a heating block maintained at  42 °C and stir gently until the BSA is completely dissolved. Adjust the final volume to  4.98 mL with [M] 0.1 Mass Percent Tris buffer.

3 Add Oleic Acid:

Thaw OA at  37 °C . Add  19.83 µL of OA to the BSA solution.

4 Incubate:

Wrap the tube with aluminum foil to protect from light. Incubate at  37 °C with gentle shaking for  02:00:00 -  03:00:00 . Ensure that the solution becomes homogeneously clear (turbidity should disappear).

5 Filter the Complex:

Filter the solution through a 0.22 µm syringe filter in a sterile culture hood. Aliquot  250 µL portions and store at  4 °C .

6 Induction of Lipid droplets:

Add OA-BSA complex to 14-21 DIV culture primary hippocampal neurons/fibroblasts to the final concentration of  200 µL -  400 µL for  24:00:00 .



7 Detection of lipid droplets:

10m

Add  1 μL BODIPY to the culture media and return the dishes to incubator for  00:10:00 . Wash the cells with PBS and mount the glass coverslip to image on confocal microscope.



Procedure to prepare Palmitic Acid–BSA (PA–BSA) Complex

4h 20m

8 Prepare Solutions:

8.1 Prepare  0.1 Mass Percent NaOH by dissolving  200 mg NaOH in  50 mL ultrapure water.



8.2 Prepare 0.9% NaCl Solution by dissolving  450 mg NaCl in  50 mL ultrapure water.



9 Dissolve BSA:

9.1 Pour  3.5 mL of 0.9% NaCl solution into a 50 ml tube. Overlay  555.55 mg of FA-free BSA on top.

9.2 Place the tube on a heating block maintained at  40 °C and stir gently until the BSA is completely dissolved.

9.3 Adjust the final volume to  4 mL with 0.9% NaCl solution. Aliquot  800 μL portions into 5 Eppendorf tubes and maintain at  42 °C on a vortex heater.

10 Dissolve Palmitic Acid:

Weigh  15.38 mg of PA and dissolve in  1 mL of  0.1 Mass Percent NaOH. Heat the mixture to  70 °C on a heating block, carefully monitoring the temperature to prevent overheating.

11 Complexing PA with BSA:





Immediately add  200 μL of the PA solution to each  800 μL aliquot of BSA maintained at  40 $^{\circ}\text{C}$ while vortexing.

Note

PA may solidify if the temperature drops.

12 Incubate:

20m

Leave PA-BSA solution under agitation for  00:15:00 -  00:20:00 . The solution may appear turbid/whitish at this stage.



13 Filter the Complex:

Filter the solution through a 0.22 μm syringe filter in a sterile culture hood. Prepare the complex fresh before use and maintain the temperature at  37 $^{\circ}\text{C}$ until used.

14 Induction of neurons:

4h

Add PA-BSA complex to the media of primary hippocampal neurons to the final concentration of  150 micromolar (μM) for  04:00:00 .



Note

- Protect fatty acid solutions from light during incubation to prevent photooxidation.
- Ensure precise temperature control during the preparation of both complexes to avoid solidification or denaturation.
- Always use sterile techniques and equipment to prevent contamination.