

Apr 24, 2024

Fecal Output Protocol

DOI

dx.doi.org/10.17504/protocols.io.rm7vzj3j5lx1/v1

Adam Hamilton¹, Ian N Krout¹, Tim Sampson¹

¹Emory University



Ian N Krout

Emory University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.rm7vzj3j5lx1/v1>

Protocol Citation: Adam Hamilton, Ian N Krout, Tim Sampson 2024. Fecal Output Protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.rm7vzj3j5lx1/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 23, 2024



Last Modified: May 31, 2024

Protocol Integer ID: 98674

Keywords: ASAPCRN, fecal output protocol this assay, fecal output protocol, number of fecal pellet, fecal pellet, measure of colonic motility, assay, colonic motility, mice, protocol, individual clear plastic beaker

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP020527

Abstract

This assay is used to quantify the number of fecal pellets produced over a short period of time, which serves as a measure of colonic motility. Mice are placed into individual clear plastic beakers and the number of fecal pellets in the container are counted every 5 minutes over a 30 minute period.

Materials

1L Translucent Beakers

Troubleshooting




Prior to Assay

- 1 Prepare clean 1 liter (~12cm x 25cm) translucent beakers- sterilize if necessary, with aluminum foil covers. These should be autoclaved before use if collection of fecal pellets for microbiome analysis is required.

Day of Set-up





1h

- 2 Bring mice to testing room for at least  01:00:00 prior to assay, to acclimate (if testing mice in the same room they are housed no adjustment period is necessary).

1h

Day of Assay

50m



- 3 Separate 1L beakers and cut/rip foil into lids to cover each beaker individually.
- 4 Place each mouse individually into a beaker, cover with foil, and start timer. Generally, 8-10 animals can be handled at once before start times begin to overlap.
- 5 ***Do not leave animals unattended while in beakers- they will attempt to jump out.*** 
- 6 In  00:05:00 intervals, gently lift each cylinder, and count the number of fecal pellets present on the beaker floor and walls. Record the number of pellets for a minimum of  00:15:00 , up to  00:30:00 .
- 7 Because mice are coprophagic, some pellets may disappear (be eaten) over this experiment. To be consistent, record only the cumulative number of pellets observed. Pellets produced cannot "decrease" so the value should either stay level or increase in each 5 minute bin over the course of this assay.
- 8 Return animals to home cage, collect fecal pellets if needed.

50m

Analysis

35m



- 9 Comparisons should be made between genotypes and treatment groups as to the amount of pellets produced in  00:30:00 . There should also be comparisons of the amount of pellets produced at each  00:05:00 interval.

35m