

Nov 20, 2024

FFPE Tissue Microarray protocol

DOI

dx.doi.org/10.17504/protocols.io.eq2ly6pppgx9/v1


Mackenzie Bolen¹, Qing-Shan Xue¹, Malu G Tansey²

¹University of Florida; ²Indiana University School of Medicine



Mackenzie Bolen

University of Florida



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly6pppgx9/v1>

Protocol Citation: Mackenzie Bolen, Qing-Shan Xue, Malu G Tansey 2024. FFPE Tissue Microarray protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.eq2ly6pppgx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 20, 2024

Last Modified: November 20, 2024

Protocol Integer ID: 112395

Keywords: ASAPCRN, TMA , NanoString CosMx SMI, ffpe tissue microarray protocol protocol, tissue microarray, nanostring cosmx smi, preparation for nanostring cosmx smi, tma, ffpe block

Funders Acknowledgements:

Parkinson's Foundation

Grant ID: PF-VSA 1293721

Aligning Science Across Parkinson's

Grant ID: ASAP-020527

Abstract

Protocol generated to as a step by step walk through of how to go from FFPE blocks to a tissue microarray (TMA) in preparation for NanoString CosMx SMI.



Materials

	Item	Cat #	Vendor
	No charge slides *do not add a sticker label either hand write or print directly on*	9035	Tissue-Tek [®] SmartWrite [®] Frosted Slides
	T-sue Microarray 1mm Mold	69132-07	T-Sue Microarray Mold Kit, 1mm, 170 Cores : Amazon
	Biopsy punch *need 1 per sample* *Tansey lab optimized with 1mm and were able to fit 40 samples per slide for CosMx*	MT3331AAP25	MILTEX Sterile Disposable Biopsy Punch with Plunger, 1mm diameter, 25/box. MFID: 33-31AA-P/25: HOSPEQ
	Cassette *if possible use smaller cassettes, generally green as opposed to the white, to decrease the depth of the paraffin*	18000-004	VWR [®] Premium Tissue Cassettes: Standard Cassettes
	Liquid paraffin/embedder	3801360	Paraplast X-Tra or Blue-ribbon paraffin from Leica Biosystems

Troubleshooting

Block generation

- 1 Put the T-sue mold in oven for 30 minutes at 70 °C
- 2 Once warm, bring tissue mold to embedder and dispense liquid paraffin into mold until the top of the core is fully submerged. Remove any bubbles with heated forceps
- 3 Place a cassette flat side down on top of mold, pick up mold with gloved hands (caution will be hot) and dispense liquid paraffin directly over the cassette until paraffin fills ½ of the depth of the cassette.
- 4 Put mold with cassette on top on cooling plate for 1 hour
- 5 After mold has cooled, carefully peel off the T-sue green mold from the cassette
- 6 Trim any excess wax around the edge of the cassette

Extracting and placing tissue into the block:

- 7 Retrieve sample blocks and place on horizontal surface. Hover the punch needle perpendicular to the preferred sampling area, then slowly insert the punch needle into the block
- 8 Retrieve the recipient block that was made, slowly push the core punched from the sample block into desired hole. Repeat steps 7-8 for all samples.
- 9 Once all samples have been loaded, make sure the top surface area is completely flat, if not carefully trim off excess tissue
- 10 Place the recipient block containing samples face down on a glass slide and incubate at 37 °C for at least 24 and up to 48 hours.

Preparing the block for sectioning: **if block is for CosMx do NOT cut until two weeks or less before experiment per NanoString manual**



- 11 Remove block from incubator and remove slide by placing on a warm plate glass slide down. Be very careful to remove block from glass slide as fast as possible. Allow block to cool.
- 12 TMA is now ready for sectioning. For 5um sections, allow block to rehydrate by placing face down on paper towel on top of ice with 1/3 volume water added.
- 12.1 Make sure to orient filled tissue area to centered approved imaging area if preparing for NanoString CosMx.
- 13 After 4 hours of rehydration, fix block to microtome and section the block