


Oct 23, 2019 Version 1

# FFPE Tissue Sectioning for Staining V.1

DOI

[dx.doi.org/10.17504/protocols.io.8kvhuw6](https://dx.doi.org/10.17504/protocols.io.8kvhuw6)



Marda Jorgensen<sup>1</sup>

<sup>1</sup>University of Florida

Human BioMolecular Atlas Program (HuBMAP) Method Development Community  
Tech. support email: [Jeff.spraggins@vanderbilt.edu](mailto:Jeff.spraggins@vanderbilt.edu)



Franchesca Farris

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.8kvhuw6](https://dx.doi.org/10.17504/protocols.io.8kvhuw6)

**Protocol Citation:** Marda Jorgensen 2019. FFPE Tissue Sectioning for Staining. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8kvhuw6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 22, 2019

**Last Modified:** December 04, 2019

**Protocol Integer ID:** 29045

## Abstract

Method for sectioning paraffin blocks for staining



## Guidelines

FFPE tissues sectioned onto poly-lysine-coated coverslips can be stored at 4° C for up to six months prior to staining. It is critical not to exceed a thickness of 10 µm because it can disrupt the autofocusing capabilities of the microscope. For best results, the tissue should be devoid of folds and tears. To ensure the integrity of the tissue slices, it is critical that they are not stacked on top of one another after they have been placed on the coverslips.

- 1 Prepare a water bath at 42°C and place it next to the Microtome.
- 2 Label positive charged microscopes slides.
- 3 Insert a new blade for sectioning each block
- 4 Section the tissue between 4-5  $\mu\text{m}$  thick.
- 5 Place the sectioned tissue in the water bath for a few seconds and observe it expanding.
- 6 When the observed expansion is enough to remove folds and wrinkles from the tissue, quickly place positive charged microscope slide in the water bath and gently move it upward towards the tissue.
- 7 Doing so, the tissue will lay on the slide as this is moved out from the water bath.
- 8 Put the slides upright in a rack to let air dry overnight or dry in 50°C oven for 1 hour.