

Sep 25, 2019 Version 1

Field sampling of root-associated microbes for DNA/RNA extraction V.1

DOI

dx.doi.org/10.17504/protocols.io.qprdvm6

Roey Angel¹

¹Soil and Water Research Infrastructure

Anaerobic and Molecular Microbiology Lab, Biology Centre CAS
Tech. support email: eva.petrova@bc.cas.cz



Roey Angel

Soil and Water Research Infrastructure

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.qprdvm6

Protocol Citation: Roey Angel 2019. Field sampling of root-associated microbes for DNA/RNA extraction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.qprdvm6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: June 04, 2018

Last Modified: September 25, 2019

Protocol Integer ID: 12753

Keywords: soil, sampling, root, rhizosphere, RNA, DNA, LifeGuard, preservation solution, field work, plant

Abstract

This protocol describes a procedure for sampling plant roots in the field for future DNA and RNA extraction for microbiome analysis. The protocol is deliberately designed to be simple and requires no electronic equipment. Root samples are preserved in LifeGuard Soil Preservation Solution for protecting against nucleic acid degradation.



Materials

MATERIALS

- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ LifeGuard Soil Preservation Solution **Qiagen Catalog #12868-100**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Technical-grade ethanol (70%) **Carl Roth Catalog #T913.1**
- ✕ Paper towels **Carl Roth Catalog #Y03.1**
- ✕ Microcentrifuge tubes 2 ml **Carl Roth Catalog #CK06.1**
- ✕ Garden trowel **Amazon**
- ✕ Disposable pasteur pipettes **Carl Roth Catalog #EA61.1**
- ✕ Tweezers set **Carl Roth Catalog #PX40.1**
- ✕ Cooling box **Carl Roth Catalog #AA46.1**
- ✕ Cooling packs **Carl Roth Catalog #E447.1**

STEP MATERIALS

- ✕ Garden trowel **Amazon**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Paper towels **Carl Roth Catalog #Y03.1**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Paper towels **Carl Roth Catalog #Y03.1**
- ✕ Tweezers set **Carl Roth Catalog #PX40.1**
- ✕ Microcentrifuge tubes 2 ml **Carl Roth Catalog #CK06.1**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ LifeGuard Soil Preservation Solution **Qiagen Catalog #12868-100**
- ✕ Disposable pasteur pipettes **Carl Roth Catalog #EA61.1**
- ✕ Cooling box **Carl Roth Catalog #AA46.1**
- ✕ Cooling packs **Carl Roth Catalog #E447.1**
- ✕ Garden trowel **Amazon**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Paper towels **Carl Roth Catalog #Y03.1**



- ✂ Scissors **Carl Roth Catalog #HCT7.1**
- ✂ Paper towels **Carl Roth Catalog #Y03.1**
- ✂ Tweezers set **Carl Roth Catalog #PX40.1**
- ✂ Microcentrifuge tubes 2 ml **Carl Roth Catalog #CK06.1**
- ✂ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✂ LifeGuard Soil Preservation Solution **Qiagen Catalog #12868-100**
- ✂ Disposable pasteur pipettes **Carl Roth Catalog #EA61.1**
- ✂ Cooling box **Carl Roth Catalog #AA46.1**
- ✂ Cooling packs **Carl Roth Catalog #E447.1**



Protocol materials

- ✕ Tweezers set **Carl Roth Catalog #PX40.1**
- ✕ Disposable pasteur pipettes **Carl Roth Catalog #EA61.1**
- ✕ Disposable pasteur pipettes **Carl Roth Catalog #EA61.1**
- ✕ Paper towels **Carl Roth Catalog #Y03.1**
- ✕ Microcentrifuge tubes 2 ml **Carl Roth Catalog #CK06.1**
- ✕ Tweezers set **Carl Roth Catalog #PX40.1**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Microcentrifuge tubes 2 ml **Carl Roth Catalog #CK06.1**
- ✕ Garden trowel **Amazon**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ Cooling box **Carl Roth Catalog #AA46.1**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ Cooling packs **Carl Roth Catalog #E447.1**
- ✕ Garden trowel **Amazon**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ Scissors **Carl Roth Catalog #HCT7.1**
- ✕ LifeGuard Soil Preservation Solution **Qiagen Catalog #12868-100**
- ✕ LifeGuard Soil Preservation Solution **Qiagen Catalog #12868-100**
- ✕ Garden trowel **Amazon**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ Paper towels **Carl Roth Catalog #Y03.1**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ Technical-grade ethanol (70%) **Carl Roth Catalog #T913.1**
- ✕ LifeGuard Soil Preservation Solution **Qiagen Catalog #12868-100**
- ✕ Cooling box **Carl Roth Catalog #AA46.1**
- ✕ Cooling packs **Carl Roth Catalog #E447.1**
- ✕ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ✕ Paper towels **Carl Roth Catalog #Y03.1**



- ☒ Cooling packs **Carl Roth Catalog #E447.1**
- ☒ Paper towels **Carl Roth Catalog #Y03.1**
- ☒ Cooling box **Carl Roth Catalog #AA46.1**
- ☒ Microcentrifuge tubes 2 ml **Carl Roth Catalog #CK06.1**
- ☒ Paper towels **Carl Roth Catalog #Y03.1**
- ☒ Tweezers set **Carl Roth Catalog #PX40.1**
- ☒ Disposable pasteur pipettes **Carl Roth Catalog #EA61.1**
- ☒ Garden trowel **Amazon**
- ☒ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ☒ Scissors **Carl Roth Catalog #HCT7.1**
- ☒ Paper towels **Carl Roth Catalog #Y03.1**
- ☒ Scissors **Carl Roth Catalog #HCT7.1**
- ☒ Paper towels **Carl Roth Catalog #Y03.1**
- ☒ Tweezers set **Carl Roth Catalog #PX40.1**
- ☒ Microcentrifuge tubes 2 ml **Carl Roth Catalog #CK06.1**
- ☒ Micro-spatula set **Carl Roth Catalog #AT16.1**
- ☒ LifeGuard Soil Preservation Solution **Qiagen Catalog #12868-100**
- ☒ Disposable pasteur pipettes **Carl Roth Catalog #EA61.1**
- ☒ Cooling box **Carl Roth Catalog #AA46.1**
- ☒ Cooling packs **Carl Roth Catalog #E447.1**

Before start

Clean spatulas using 70% ethanol

G

M

T

<input type="checkbox"/>	Detect language	English					
--------------------------	-----------------	---------	--	--	--	--	--

Text-to-speech function is limited to 200 characters

<input checked="" type="checkbox"/>	Options : History : Feedback : Donate	Close
-------------------------------------	---	-------



- 1 **Sample a triplicate of plant individuals spaced out a few metres apart from each other**Make sure you are reall sampling individual plants and not offshoots of the same plant
- 2 **Using a garden trowel, carefully dig out the plants while keeping the root system intact** (as much as possible, of course)

G
M
T

<input type="checkbox"/>	Detect language	English					
--------------------------	-----------------	---------	--	--	--	--	--

Text-to-speech function is limited to 200 characters

<input checked="" type="checkbox"/>	Options : History : Feedback : Donate	Close
-------------------------------------	---	-------

 Garden trowel **Amazon**

- 3 **While holding the plant by the shoot, shake the root system hard enough so that all loose soil is removed from it.**Take care to damage the plant as little as possibleYou can use a spatula to remove large soil aggregates that are attached to the roots

 Micro-spatula set **Carl Roth Catalog #AT16.1**

- 4 **From the remaining root system (plus soil particles plus attached to the roots), trim a 'representative sample' of roots using scissors or scalpell**It is usually best to trim the



roots onto a piece of paper towel

 Scissors Carl Roth Catalog #HCT7.1

 Paper towels Carl Roth Catalog #Y03.1


5 Cut the trimmed out roots a little so that they fit into a 2.0 ml tube

 Scissors Carl Roth Catalog #HCT7.1

 Paper towels Carl Roth Catalog #Y03.1

6 Place about 2-3 g of that cut out sample into a 2.0 ml tube

 Tweezers set Carl Roth Catalog #PX40.1

 Microcentrifuge tubes 2 ml Carl Roth Catalog #CK06.1

 2 g

Note

The root tissue should make up at least half or more of the mass, while the remaining attached soil should make up the rest


7 Press the sample a little into the bottom of the tube to decrease its volume

 Micro-spatula set Carl Roth Catalog #AT16.1


Note

Make sure the roots do not take up more than 1/2-2/3 of the volume
It is of course possible to split each sample into several separate tubes, depending on the specific type of roots, and submerge each with LifeGuard solution

8 Add as much LifeGuard solution so that the sample is submerged in about twice of its volume (about 1.0 – 1.5 ml) Best is to use a disposable Pasteur pipette for dispensing the solution

 LifeGuard Soil Preservation Solution Qiagen Catalog #12868-100

 Disposable pasteur pipettes Carl Roth Catalog #EA61.1

 1.5 mL



- 9 **Place the tubes in cooling (around 4 °C) and keep them cooled until you reach the lab. The solution will protect nucleic acids even at room temperature for several days, but cooling is preferred**

 Cooling box Carl Roth Catalog #AA46.1

 Cooling packs Carl Roth Catalog #E447.1

 4 °C

- 10 **In the lab, store the samples in a freezer (-20 – -80 °C)**

 -20 °C or -80 °C