

Jun 11, 2024

Filter trap assay for the detection of alpha-synuclein aggregation

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1>

Cole S Sitron¹, Victoria A Trinkaus¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Cole Sitron

Max Planck Institute of Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1>

Protocol Citation: Cole S Sitron, Victoria A Trinkaus, F Ulrich Hartl 2024. Filter trap assay for the detection of alpha-synuclein aggregation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.x54v92re4l3e/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 17, 2024



Last Modified: June 11, 2024

Protocol Integer ID: 101394

Keywords: synuclein aggregation, synuclein in cell lysate, synuclein aggregation this protocol, synuclein, detection of alpha, assay

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol describes an assay that detects aggregated alpha-synuclein in cell lysates.

Materials

Cell culture

Cells with alpha-synuclein aggregates, induced as in [[dx.doi.org/10.17504/protocols.io.eq2lyjbbplx9/v1](https://doi.org/10.17504/protocols.io.eq2lyjbbplx9/v1)], growing in 12-well plates.

Reagents

- TrypLE™ Express Thermo Fisher Scientific Catalog #12605010
- Pierce™ Rapid Gold BCA Protein Assay Kit Thermo Fisher Catalog #A53225

Buffers

- 10% FBS (Fetal Bovine Serum Gibco - Thermo Fisher Scientific Catalog #10270106) in 1X PBS (diluted from PBS (10X), pH 7.4 Thermo Fisher Scientific Catalog #70011051)
- 1X PBS
- Lysis Buffer:

	A	B
	Tris-HCl pH 7.5	50 mM
	NaCl	150 mM
	Triton-x100	1% (v/v)
	cOmplete protease inhibitor tablet	1 mini tablet per 15 ml buffer (Sigma-Aldrich cat. no. 4693159001)
	Phos-STOP phosphatase inhibitor tablet	1 tablet per 15 ml buffer (Sigma-Aldrich cat. no. 4906845001)
	Benzonase ((Max Planck Institute of Biochemistry Core Facility)	7.5 U/ml

- Wash Buffer:

	A	B
	Tris-HCl pH 7.5	50 mM
	NaCl	150 mM
	Triton-x100	1% (v/v)

Consumables

- 0.2 micrometer pore size cellulose acetate membrane (e.g. GE cat. no. 10404131)
-  PARAFILM® M Merck MilliporeSigma (Sigma-Aldrich) Catalog #P7793

Equipment

Equipment	
Bioruptor® Plus sonication device	NAME
Sonication device	TYPE
Bioruptor®	BRAND
B01020001	SKU
https://www.diagenode.com/en/p/bioruptor-plus-sonication-device	LINK

- Slot blot vacuum manifold (e.g. Hoefer PR648)
- Fume hood



Harvesting cells

- 1 Prepare lysis buffer and place On ice .
- 2 After 24:00:00 of exposure to PFFs/PBS, wash wells 1X with 1 mL warm PBS. 1d
- 3 Withdraw PBS and treat wells with 300 μ L of TrypLE.
- 4 Quench dissociation reagent with 300 μ L of 10% FBS and transfer the contents of the wells to 1.5 ml tubes.
- 5 Centrifuge 1500 x g for 00:03:00 . 3m
- 6 Withdraw the supernatant and wash cells with 1 mL PBS.
- 7 Centrifuge 1500 x g for 00:03:00 . 3m
- 8 Withdraw supernatant and resuspend in 100 μ L of lysis buffer and incubate On ice for 00:05:00 . 5m
- 9 Sonicate the cells in the Bioruptor for 3 cycles of 00:00:30 on and 00:00:30 off at high power. 1m
- 10 Centrifuge the lysate at 500 x g for 00:05:00 and place the clarified lysate in a new 1.5 ml tube. 5m

**Note**

This step is crucial to ensure the removal of high-diameter material (e.g. unlysed cells) that may clog the filter paper and non-specifically trap proteins on it.

Note

If desired, this is a convenient stopping point, wherein you can snap-freeze your lysate in liquid nitrogen and proceed with the rest of the assay at a later time.

- 11 Determine the concentration of protein in each sample by Pierce Rapid Gold BCA Protein Assay.

Filter Trap

10m

- 12 Soak the cellulose acetate membrane in wash buffer for >  00:10:00 .

10m

- 13 For each well that you run, prepare  1020 μL of the sample containing 10.2 - 51 ug of lysate, with the remaining volume consisting of lysis buffer. 

- 14 Assemble the slot blot vacuum manifold apparatus with the cellulose acetate membrane.

- 15 Wash each well with  300 μL wash buffer, inspecting each well to ensure that the buffer can drain through it.  

Note

Avoid using any non-functioning wells.

- 16 Load  1000 μL (10 - 50 ug) of sample to each well. 



- 17 Once the sample has drained through the wells, wash each well 3 times with  300 μ L of wash buffer.
- 18 When all of the washes have drained through the wells, trim the membrane and proceed with standard immunodetection methods.

